

labCA Tutorial

Till Straumann



Overview

- What is labCA
- General Concepts
 - PV s, Timestamps, Timeout, Error Handling, Help
- Basic Commands
 - IcaGet, IcaPut
- Monitors
- Examples
- Summary



What is labCA

- EPICS / Channel-Access interface for matlab and scilab (open-source, more features than octave but language not matlab-compatible)
- Builds on top of modified ezca library which has been made thread-safe. Use multi-threaded EPICS-3.14.
- read, write and monitor PVs
- relation to MCA: different history. I developed labCA for scilab, as it evolved matlab support was added



labCA 3.x Features

- handle multiple PVs in a single call
- retrieve data and time-stamp
- support for string data
- error handling
- asynchronous operation
- abort labCA with <Ctrl-C>
- on-line help and detailed documentation
- tested on linux, solaris, win, matlab, scilab, 32 and 64-bit.



General Concepts: PVs

PVs are passed as a column vector of strings (PV names) (cell-array on matlab):

```
my_pvs = {'one_pv'; 'second_pv'; 'third_pv'};
lcaGet( my_pvs )
```

(no need for curly braces when dealing with single PV)

```
lcaGet( 'another_pv' )
```



General Concepts: Timestamps

- Timestamps have a 'quirky' format:
 - POSIX-timespec: number of nanoseconds since 0:0:0 1/1/1970 UTC is given as a *complex number*:

full_seconds + j nanoseconds_remainder

- Rationale:
 - I wanted to avoid more complex objects and stick with numbers
 - easy separation of seconds and sub-second parts (LCLS, SSRL embed 'pulse-ID' in nanosecond part)
 - standard POSIX format; utilities for conversion to local time available.



General Concepts: Timeout

- Synchronous Ica routines block until request completes or times out
 - IcaGet, IcaPut, IcaNewMonitorWait, ...
- Total timeout is product of 'timeout' * 'retry_count'. Inspect and change with

```
lcaGetTimeout(), lcaSetTimeout( new_timeout )
lcaGetRetryCount(), lcaSetRetryCount( newcnt )
```

Matlab is polled for <Ctrl-C> activity every 'timeout' (recommended: timeout=0.5s..1s, retry_count = total_timeout / timeout)



General Concepts: Errors

- Errors are reported by 'throwing' exceptions; use try/catch/lasterror.
- If an EPICS PV read with lcaGet() has a non-OK error status then a warning message is printed.
- If an EPICS PV has severity INVALID then IcaGet() returns NAN.
- Both of these features may be disabled and/or tuned; consult lcaSetSeverityWarnLevel().



General Concepts: Help

- On-line documentation for all lcaxxx() calls available
- RT(F)M; available in PDF and HTML.
- latest versions always @

www/~strauman/labca



Basic Commands: lcaGet/lcaPut

Read a PV:

```
[ value, timestamp ] = lcaGet( <pv_name> )
```

Write a PV

```
lcaPut( <pv_name>, <value> )
```

Additional, optional arguments available (RTM)



Asynchronous Operation

- IcaPut () is synchronous, i.e., waits for completion of the write operation on the server (IOC). In some cases (e.g., driving stepper motor to target position) this may take a very long time.
- lcaPutNoWait() is asynchronous. It does not wait for completion (and you won't be notified if the operation fails). Merely sends the 'write' command and returns (AKA, 'posted-write').



Asynchronous Op.: Monitors

- lcaSetMonitor(<pvs>) registers a subscription
 (monitor) for a PV. Every time the value of the PV
 changes labCA is notified and a background task reads
 the new value into a internal buffer
- lcaGet() of a monitored PV reads from internal buffer
- Check if the background buffer has new data (since last lcaGet()):
 - lcaNewMonitorValue(<pvs>) returns flag vector
 (nonzero value flags new data)
 - lcaNewMonitorWait(<pvs>) blocks until all PVs have
 new data available.

Feb. 14, 2008



Monitors (cont.)

■ Subscription remains active until connection is released (lcaClear (<pv>))



Monitors Example

Trigger instrument and read new data (w/o monitor: how do we know if we have old or fresh data?)

```
my_pvs = {'det_1'; 'det_2' }; // helper

lcaSetMonitor( my_pvs ); // establish monitor

lcaNewMonitorWait( my_pvs ); // initial buffer fill

lcaGet( my_pvs ); // discard initial value

trigger_data_acq(); // trigger;

lcaNewMonitorWait( my_pvs ); // wait for new data

data = lcaGet( my_pvs ); // read new data

lcaClear( my_pvs ); // clear channel/monitor
```



Connection Management

- Background:
 - EPICS-Channel Access maintains one TCP connection for each IOC <-> client connection
 - multiple PVs on same IOC share one connection
 - multiple clients (matlab instances) to same IOC have separate connections
 - ■labCA establishes connections transparently and ondemand but never automatically releases them (except if matlab terminates). [Connection 'aging' and garbagecollection could be feature-addition.]



Connection Management (cont)

- labCA can release connections explicitly:
- lcaClear() releases all connections and monitors
- IcaClear(<pvs>) releases monitors and channels on named PVs – connection is released when last channel (=PV) using it is cleared.
- labCA cannot release monitor w/o clearing channel
- next time a cleared PV is used a new channel is created transparently (but takes a little bit longer).
- avoid frequent clearing/re-creation of channels; if done automatically at high frequency broadcast storms may result. Occasional clearing/re-connecting is OK.

