
Bergoz Instrumentation
Espace Allondon Ouest
01630 Saint Genis Pouilly, France
Tel. : +33-450.426.642
Fax : +33-450.426.643



Instrumentation

BCM-RF-E

Beam Charge Monitor - RF receiver

User's Manual

Rev. 1.1

U.S.A. :

GMW Associates
www.gmw.com
955 Industrial Rd.
San Carlos, CA 94070,
U.S.A.
Tel.: (650) 802-8292
sales@gmw.com

Japan :

REPIC Corp.
www.repic.co.jp
28-3 Kita Otsuka 1-Chome
Toshima-ku, Tokyo 170-0004,
Japan
Tel.: 03-3918-5326
sales@repic.co.jp

India :

GEEBEE International
www.geebinternational.com
12 - A, Lane W - 16,
Sainik Farms, New Delhi,
India
+91 98107 54666
sahil@geebinternational.com

Record of updates

Version:	Release date:	Change:
0.0	Dec. 2012	First review
0.1	Jan. 2013	Connectors pins allocation : corrected & updated – Summary : cleaning – Sensitivity : small change
0.2	Mar. 2013	DB9 Remote control section removed – Architecture : DB9 remote control removed - Pin allocations table updated – Summary updated
0.3	Jul. 2013	USB update.
1.0	Feb. 2014	SENSITIVITY : Qcal and Ucal added, log10() instead of exp(). GUI section added, FIRMWARE section added. CALIBRATION section added. ARCHITECTURE : DB9,6ADC.TRG.OUT added.
1.1	Dec. 2014	General improvements and corrections

SUMMARY

INITIAL INSPECTION, WARRANTY, ASSISTANCE	2
YOU JUST RECEIVED A BEAM CHARGE MONITOR	3
Determine your BCM-RF operating mode (charge or current measurement)	3
ACCESS to SWITCHES	4
OPERATING MODE SELECTION.....	5
QUICK CHECK	7
Front panel	8
Charge measurement, select Sample & Hold mode (S&H), to measure bunch charge	9
Waveforms	10
ADC trigger waveforms	11
Sample-to-sample noise	12
Current measurement, select Track-Continuous mode (T-C), to measure CW beam or long macropulses	
average current	13
Waveforms and reaction time.....	14
Noise	15
Wide-band output power spectrum	18
I/O AND SWITCHES.....	19
GENERAL DESCRIPTION	20
SPECIFICATIONS.....	22
Charge measurement:	22
Current measurement:	22
Input / output interfaces and others . . .	23
Power supply and fuses	23
SENSITIVITY of the TURBO-ICT and BCM-RF.....	24
GRAPHICAL USER INTERFACE	25
BCM-RF FIRMWARE.....	27
CONNECTORS PINS ALLOCATION	29
BCM CABLE LAYOUT INSTALLATION	30
ACCESSORIES	31
BCM Chassis BCM-RFC/XX.....	31
Chassis front.....	31
Card Extender BCM-XTD.....	31
SCHEMATICS & BOARD LAYOUT	33
ACKNOWLEDGEMENTS.....	33

INITIAL INSPECTION, WARRANTY, ASSISTANCE

INITIAL INSPECTION

It is recommended that the shipment be inspected immediately upon delivery. If it is damaged in any way, contact Bergoz Instrumentation or your distributor. The content of the shipment should be compared to the items listed on the invoice. Any discrepancy should be notified to Bergoz Instrumentation or its distributor immediately. Unless promptly notified, Bergoz Instrumentation will not be responsible for such discrepancies.

WARRANTY

Bergoz Instrumentation warrants its beam monitors to operate within specifications under normal use for a period of 12 months from the date of shipment. Spares, repairs and replacement parts are warranted for 90 days. Products not manufactured by Bergoz Instrumentation are covered solely by the warranty of the original manufacturer. In exercising this warranty, Bergoz Instrumentation will repair, or at its option, replace any product returned to Bergoz Instrumentation or its distributor within the warranty period, provided that the warrantor's examination discloses that the product is defective due to workmanship or materials and that the defect has not been caused by misuse, neglect, accident or abnormal conditions or operations. Damages caused by ionizing radiations are specifically excluded from the warranty. Bergoz Instrumentation and its local distributors shall not be responsible for any consequential, incidental or special damages.

ASSISTANCE

Assistance in installation, use or calibration of Bergoz Instrumentation beam monitors is available from Bergoz Instrumentation, 01630 Saint Genis Pouilly, France. It is recommended to send a detailed description of the problem by e-mail.

SERVICE PROCEDURE

Products requiring maintenance should be returned to Bergoz Instrumentation or its distributor. Bergoz Instrumentation will repair or replace any product under warranty at no charge. The purchaser is only responsible for transportation charges.

For products in need of repair after the warranty period, the customer must provide a purchase order before repairs can be initiated. Bergoz Instrumentation can issue fixed price quotations for most repairs. However, depending on the damage, it may be necessary to return the equipment to Bergoz Instrumentation to assess the cost of repair.

RETURN PROCEDURE

All products returned for repair should include a detailed description of the defect or failure and name of the user. Contact Bergoz Instrumentation or your local distributor to determine where to return the product. Returns must be notified by e-mail prior to shipment.

Return should be made prepaid. Bergoz Instrumentation will not accept freight-collect shipment. Shipment should be made via United Parcel Service, DHL or Federal Express. Within Europe, the transportation service offered by the Post Offices "EMS" (Chronopost, Datapost, etc.) can be used. The delivery charges or customs clearance charges arising from the use of other carriers will be charged to the customer.

YOU JUST RECEIVED A BEAM CHARGE MONITOR

This manual applies only to BCM-RF revisions 204.2 and above

Its does NOT apply to either earlier BCM-RF revisions, BCM-IHR or BCM-CA.

The BCM-RF system includes:

Description Order code

- BCM-RF electronics module BCM-RF-E/XXXXMHz
 Where XXXX is operating frequency

Accessories

- 19" chassis with power supply BCM-RFC/XX,
 Where XX = number of BCM stations

Note : BCM-RFC/XX RF-shielded chassis is compatible with BCM-IHR, BCM-CA.

Hence, BCM-RF, BCM-IHR and BCM-CA can be mixed in same BCM-RFC/XX chassis.

Determine your BCM-RF operating mode (charge or current measurement)

- Sample & Hold mode (S&H) for bunch charge measurement.
- Track-Continuous mode (T-C) for CW or long macropulse beam current measurement.


BCM-RF-E is factory configured either for beam charge measurement or beam current measurement. The operating mode can be selected with mechanical switches or by USB remote control (see **OPERATING MODE SELECTION**).

The ex-factory mode of operation is defined in the Test report and in the Certificate of Calibration.

ACCESS to SWITCHES

To access switches, remove the module shield:

Note: It is possible to remote control BCM-RF-E via the USB port (*) without removing module shield.




(*) Available only from #006 BCM-RF serial number. Firmware revision is written on the microcontroller.

OPERATING MODE SELECTION


By using on-board switches ... (manual mode)

WARNING : the PIC microcontroller has the priority on switches control. To use the mechanical switches, please remove the microcontroller from its slot.

Track-Continuous mode (T-C)



Sample and Hold mode (S&H)



OPERATING MODE SELECTION (Cont'd)

By remote control via USB port


Available only from #006 BCM-RF-E serial number. Bergoz Instrumentation provides a Labview Graphical User Interface to control and read the BCM-RF-E. This software is available from our website. See GRAPHICAL USER INTERFACE section.

Name : BCM-RF-GUI-x.x

Extension : executable file (.EXE)

Requires : Labview Run Time Environment 7.1

Size : < 500 KB



QUICK CHECK

You can check immediately that your BCM is working. This is what you need:

- BCM-RF-E (Beam Charge Monitor RF): at least one BCM-RF-E electronics module in...
- BCM-RFC/xx chassis
- Turbo-ICT (Turbo-Integrating Current Transformer)
- 4-channel oscilloscope with 500 MHz bandwidth or higher
- Fast pulse generator <200 ps fwhm, <200 pC, <2 MHz repetition rate or with external trigger
- Pulse generator used as trigger: >2 V pos. edge, $T_w > 100$ ns, <2 MHz repetition rate
- Short (4...8 ns) BNC cables and SMA-BNC adapters.

Verify that this manual corresponds to your BCM version

The modular Electronics "Beam Charge Monitor, RF" (BCM-RF-E) version is the object of this manual. Other User's Manuals cover other versions: BCM-IHR-E and BCM-CA-E.

QUICK CHECK ...

Front panel

BCM-RF-E module



Signal View, Log signal envelope (1 Mohm or 50 ohm)

Signal View

Output signal, proportional to the input charge logarithm (1 Mohm)

Output View

Hold View, rising edge indicates time when log value is held (50 ohm, only in S&H mode)

Hold View

Front panel LED : blinks in S&H mode when an input signal or external trigger arrives.

USB connector type B : data readout and remote control

USB

Hold Delay front panel trimmer used to fine adjust the hold time.
Do not change if Turbo-ICT-CAL-FO option is installed!

Hold Delay

WARNING: Jumpers configuration & Potentiometers settings

Your BCM is in the "Ex-factory" configuration.


Jumper and timing adjustments (potentiometers) have been configured according to your order.
Do not change those settings until you are familiar with the Beam Charge Monitor.

QUICK CHECK

Charge measurement, select **Sample & Hold mode (S&H)**, to measure bunch charge

See OPERATING MODE SELECTION

SETUP



Connect the Turbo-ICT to "BCM Input" SMA connector on chassis rear panel.

Make sure your AC mains voltage corresponds to BCM-RFC chassis mains voltage.

Note: A label with the set AC mains voltage is affixed to the power supply unit front panel


Connect BCM-RFC chassis AC input to the mains.

- Apply a 50 ohm fast pulse through the Turbo-ICT aperture:
- Pulse charge: $< 200 \text{ pC}$.
- Pulse polarity: any.
- Pulse fwhm $< 200 \text{ ps}$.
- Pulse repetition rate: $< 2 \text{ MHz}$ (expected to be triggered by pulse generator below)
- Apply a 50 ohm pulse from generator to "Trigger in" SMA connector on BCM-RFC/XX chassis rear panel:
- Pulse width: $> 30 \text{ ns}$.
- Pulse polarity: rising edge.
- Pulse amplitude: $> 2.4 \text{ V}$ in 50 ohm
- Pulse repetition rate: $< 2 \text{ MHz}$ ($> 500 \text{ ns}$).

Apply same pulse to oscilloscope trigger input.

Charge measurement, Sample & Hold mode (cont'd)

Waveforms



Channel 4 (green) shows the input short pulse which goes thru the Turbo-ICT aperture.

Channel 2 (pink) shows “BCM Out”. This voltage is held constant for up to 100ms after the pulse.

Channel 3 (cyan) shows the “Hold View”. The rising edge is an indication when the hold time occurs. But it is too imprecise to set up trigger timing.

Channel 1 (orange) shows the “Signal View”. This is the logarithmic envelope of the signal from the Turbo-ICT. The voltage apex is proportional to the logarithm of the beam charge.


To adjust the hold time:

- Monitor the traces above with an oscilloscope. Check all waveforms.
- Adjust the hold time using the GUI or by adjusting the external trigger delay. Check that the rising edge of “Hold View” is close to the apex of “Signal View”.
- Connect a digital voltmeter to “BCM Out”.
- Fine adjust hold time to maximize the output voltage. The rising edge of “Hold View” might not coincide with the “Signal View”.
- Once installed in the accelerator, hold time has to be set up again to account for changes in the signal delay between Turbo-ICT and BCM-RF.

Charge measurement, Sample & Hold mode (cont'd)

ADC trigger waveforms

Observe ADC Trigger output "TRG.ADC.OUT" available from the DB9 connector pin 3 on the chassis rear panel.



BCM Output signal is ready to be used on the rising edge of TRG.ADC.OUT which occurs about 500 ns after the input pulse.

Your ADC can be triggered from DB9 pin 3 output signal TRG.ADC.OUT:


- First edge is 200 ns after the beam bunch.
- Second edge is 500 ns after the beam bunch.
- Polarity of edge is determined by ADC Trigger Edge Polarity switch.

Charge measurement, Sample & Hold mode (Cont'd)

Sample-to-sample noise

The BCM Output exhibits noise. The following test gives an idea of the noise, without the need for a sophisticated data acquisition system.

An oscilloscope takes a single sweep of Output View. Its time base is set to display –in a single sweep– the output corresponding to ~ 25 input pulses.




From a succession of 15 pC input pulses, the BCM-RF provides an output noise standard deviation of 2.5 mVrms. Applying the calibrated reverse "log" function from an ICT-Turbo1 to this output noise gives an equivalent input noise of 0.11 pCrms. One can see that the signal over noise ratio is 136, equivalent to 0.7 % resolution.

QUICK CHECK

Current measurement, select Track-Continuous mode (T-C), to measure CW beam or long macropulses average current

See OPERATING MODE SELECTION

SETUP



Connect the Turbo-ICT to "BCM Input" SMA connector on chassis rear panel.

Make sure your AC mains voltage corresponds to BCM-RFC chassis mains voltage.


Note: A label with the set AC mains voltage is affixed to the power supply unit front panel.

Connect BCM-RFC chassis AC input to the mains.

- Apply a 50 ohm sinewave signal into Turbo-ICT aperture.
- Sinewave centre frequency: BCM-RF-E operating frequency marked on BCM-RF-E label located on DIN41612M connector, or Certificate of Calibration.
- Output power: <-30 dBm.

Current measurement, Track-Continuous mode (T-C) (Cont'd)

Waveforms and reaction time




The signal Output View is proportional to the logarithm of the current from the generator.

The reaction time, measured from the input current change until the output value rises by 50%, is about 110 ns.

Current measurement, Track-Continuous mode (T-C) (Cont'd)


Noise



Output View contains ~2.7 mVrms of noise when -30 dBm signal is applied through the ICT-Turbo1 aperture.


The equivalent input noise taken from the calibrated reverse function is 1 uArms i.e. 0.7% resolution.

Current measurement, Track-Continuous mode (T-C) (Cont'd)



The noise (standard deviation) from -60 dBm signal applied through ICT-Turbo1 aperture is about 21 mVrms. The equivalent input noise taken from the calibrated reverse function is 0.14 uArms i.e. 0.7% resolution. Notice that the resolution remains the same although the signal decreases.

Current measurement, Track-Continuous mode (T-C) (Cont'd)





Without input signal, Output View is noisy because in the absence of an input signal the logarithmic amplifier gain increases.


When no signal is applied through the Turbo-ICT aperture, the noise on Output View is typically some 10 millivolts.

Current measurement , Track-Continuous mode (T-C) (Cont'd)

Wide-band output power spectrum



I/O AND SWITCHES



GENERAL DESCRIPTION

BCM-RF electronic modules (-E) embeds two modes of operation:

T-C mode: BCM-RF output provides the log envelope of the Turbo-ICT modulation in 10 MHz bandwidth. Suitable for CW or macropulse beam current measurements.

S&H mode: BCM-RF output holds the log envelope voltage apex of the Turbo-ICT modulation. The output is a DC level held up to 100ms or until the next bunch arrives. This mode allows repetition rates up to 2 MHz. Suitable for single bunch charge measurements.

This manual describes both modes.

System components

In a beam line or particle accelerator application, the BCM detects the beam signal with a non-destructive sensor:

- Turbo Integrating Current Transformer (Turbo-ICT).

Electronics are housed in a 3U-high, 19"-wide RF-shielded chassis, which can hold:



- BCM-RF-E.
- Power Supply.

The BCM output is a voltage up to +5V, proportional to the log of the beam charge.

In "S&H" mode, the output voltage level is held up to 100ms or until a next bunch arrives.

In "T-C" mode, the output voltage waveform is the log envelope of the Turbo-ICT modulation signal in 20 MHz bandwidth.

ARCHITECTURE



SPECIFICATIONS

Typical performance measured with Turbo2 option:

Charge measurement:

BCM-RF mode	S&H
Input charge	300 pC max (more without front-end amplifier)
Measurement single range	500 fC – 300 pC
Bunch repetition rate	Single bunch – 2 MHz
Output voltage	0 / 5V, log of the beam charge
Reaction time	500 ns to > 99% final value
Noise	50 fC or 1 % of charge (whichever is higher)
Non-linearity	3%
Time response	Hold till next bunch
Trigger	on-board or external (pos. edge, > 2.4V in 50 ohm, > 30 ns width)

Current measurement:

BCM-RF mode	T-C
Measurement single range	5 uA – 20 mA (more without front-end amplifier)
Beam RF	1 MHz – 350 MHz
Output voltage	0 / +5V, log of beam current
Risetime	< 70 ns
Reaction time	100 ns
Noise	0.1 uA or 0.2% of current (whichever is higher)
Non-linearity	3%
Time response	reports current variation to 10 MHz bandwidth

Input / output interfaces and others ...

Front-panel connectors (BNC) :	Signal View (for oscilloscope, 50 ohm or 1 Mohm) Output View (for oscilloscope, 1 Mohm) Hold View (for oscilloscope, 50 ohm)
Rear module connector :	DIN 41612-M / 24+8 male, with 1.0/2.4 coaxial inserts
Back-panel connectors (SMA) :	BCM Input, 50-ohm coaxial cable from Turbo-ICT BCM Output, for high-impedance readout
Trigger Input :	50 ohm, pos. edge, > 2.4V, > 30 ns width
Front-panel potentiometers :	Hold time delay adjustment
Power consumption :	+15V, 170 mA / -15V, 110 mA (Turbo-ICT connected)
Card size :	3U x 4F, i.e. Eurosize 100 x 160 mm, 20mm wide
Chassis size :	3U x 19"

Power supply and fuses

See Annex III: Delta Elektronika U-Series linear power supply data sheet

The mains voltage is factory set according to the label stuck to the front panel.

Please remove this label when you change the mains voltage selection.


Type :	5U 15-15 modular plug-in $\pm 15V$ linear power supply
Manufacturer :	Delta Elektronika, 4300A Zierikzee, The Netherlands
Output voltage :	$\pm 15V$, 200 mA
Mains voltage :	jumper selectable, 110, 220 Vac, 50-60 Hz tested at 90 Vac-50/60 Hz for 100 Vac Japanese mains voltage.
Mains voltage selector :	located under the power supply block
Card size :	3U x 10F, i.e. Eurosize 100 x 160 mm, 50mm wide
Back-panel connector :	The Power supply mains are wired to a IEC connector via EMI/RFI filter and fuse.

SENSITIVITY of the TURBO-ICT and BCM-RF

BCM Output gives a voltage proportional to the logarithm of the input bunch charge. The sensitivity is defined by the reverse transfer function of the complete signal chain from Turbo-ICT to BCM-RF-E.

The sensitivity of a given combination of Turbo-ICT, BCM-RF-E and connecting cable is provided in the Calibration Report. Each combination has its own transfer function.

This section describes the sensitivity calculation from the beam to the BCM-RF-E output:



Reverse transfer function (sensitivity):

$$Q_{in} = Q_{cal} * 10^{\left(\frac{BCM_{out}}{U_{cal}}\right)} \quad (\text{S - H mode})$$

$$I_{in} = I_{cal} * 10^{\left(\frac{BCM_{out}}{U_{cal}}\right)} \quad (\text{T - C mode})$$

U_{cal} and I_{cal} are measured in Volts and Q_{cal} is measured in pC. These calibration constants are factory measured and written on the Calibration Report.

Notice that losses in the cable connecting Turbo-ICT and BCM-RF were taken into account during calibration. Only in case the expected cable losses (see Calibration Report) differ from real cable losses, additional corrections have to be applied by the user. Real cable losses can be obtained, for example, by network analyser measurements. All losses from the Turbo-ICT output connector to the BCM-RF input connector must be taken into account.

Cable losses only affect Q_{cal} and I_{cal} . Hence, any correction is just an additional factor applied to the reverse transfer function:

$$Q_{in} = Q_{cal} * \frac{\text{Real Cable Losses}}{\text{Expected Cable Losses}} * 10^{\left(\frac{BCM_{out}}{U_{cal}}\right)} \quad (\text{S - H mode})$$

$$I_{in} = I_{cal} * \frac{\text{Real Cable Losses}}{\text{Expected Cable Losses}} * 10^{\left(\frac{BCM_{out}}{U_{cal}}\right)} \quad (\text{T - C mode})$$

As long as the applied correction remains small, overall performance is not affected. But in case larger corrections must be applied lower and upper limits could shift. Non-linearities will change, too. In such a case, Bergoz Instrumentation can provide further advice.

GRAPHICAL USER INTERFACE

Bergoz Instrumentation provides a demo GUI to control and read the BCM-RF via USB.

This software is an executable file developed on Labview 7.1 (will be updated to newer version).

It uses Windows generic USB to serial port driver. No external driver is required.


Operating systems supported:

In principle, any operating system that can run LABVIEW 7.1
or the corresponding run time environment, e.g. Win 2000, XP, Vista, 7

LABVIEW run time environments can be obtained from the National Instruments web site.

To control the BCM-RF from USB cable :

- Connect the USB cable from the BCM-RF to the PC.
- Windows recognizes the device and installs the generic serial to USB driver.
- In the devices manager, find the associated serial COM port number.
- On the GUI select the COM port number connected to the BCM-RF.
- Run the GUI.
- Data_In LED starts to blink which means that data comes from BCM-RF.
- Configure the BCM-RF to match your application.
- Store the customized configuration to the BCM-RF's microcontroller by clicking on :
"Configuration save to BCM-RF".



GRAPHICAL USER INTERFACE (Cont'd)

Table shows switches position for different measurement mode. User must set GUI switches to perform the measurement.

Measurement Mode	CAL-FO	Trigger Mode	Internal Clock	Hold Delay Control	Trigger ADC Out
S&H	Disable	On Board or External	Enable	Digital	Blinks on new trigger
T-C	Disable	Not Used	Not Used	Not Used	Not Used
CAL-FO	Enable	X	X	X	X

Measurement Mode :

- S&H – Sample and Hold mode measures the bunch charge. It works with single pass up to 2 MHz repetition rate. This mode uses a trigger and internal clock (must be activated).
- T-C – Track Continuous mode measures the average beam current in 10 MHz bandwidth. This mode is a free running process, it does not use any clock or trigger.

Trigger Mode :

- On-board – BCM-RF embeds an on-board trigger synchronized on the incoming signal from Turbo-ICT.
- External – BCM-RF can be triggered from an external TTL rising edge. “Trigger in” SMA connector is located on the rear chassis.

Internal clock :

- BCM-RF internal clock. Must be enabled for S&H mode.

Hold Delay :

- BCM-RF has an internal digital delay line USB controlled by the GUI. The user can adjust the hold delay from 0 ns to 255 ns.

CAL-FO (optional) :

- Online calibration check activation. Only available as option. This is a calibration generator installed on Turbo-ICT and controlled by fiber optics.

Trigger ADC Out :

- Indicator blinks when a new signal arrives. Only available in S&H mode

Averaging :

- User can change the sampling average number of points from 1 to 65535.

Reverse transfer function:

- Reverse function calculates the real bunch charge from the BCM-RF output voltage by using the given calibration constants. Per default the reverse function is not activated and the output voltage is plotted. When the reverse function is activated the calculated charge is plotted (pC). Currently the reverse function applies only to S-H mode.

Configuration saved to BCM-RF :

- BCM-RF configuration can be stored into microcontroller's EEPROM.

BCM-RF FIRMWARE

BCM-RF embeds a PIC 18F2458 from Microchip Technology Inc. This microcontroller has a 12 bits ADC and USB registers.

The firmware is written in C with MPLAB 8.9 and MPLAB C18 (C compiler) available from the Microchip website.

The code can be downloaded from www.bergoz.com on the Turbo-ICT / BCM-RF web page. Users can modify the code to match the BCM-RF to their own application.

File name : "BCM-RF-Firmware-x.x".

To program and debug the microcontroller, remove the cover shield and connect an ICD3 Microchip In-circuit debugger to the RJ11-R connector (see I/O AND SWITCHES section).

BCM-RF uses Communication Devices Class USB protocole in POLLING mode.

All characters are transmitted in hexadecimal.

A general frame used to write or read the microcontroller looks like this:

1 byte : Command	1 byte : Frame number of byte	1 byte : Read / write	2 bytes : Frame counter	1 byte : Separation character	2 bytes : Value	2 bytes : End characters
'A' .. 'Z' In Ascii	0 to 9	'.' write frame '?' read frame In Ascii	0 to 65535	'=' In Ascii	0 to 65535	/n/0 In Ascii

The analog BCM-RF output signal is periodically sampled by the microcontroller's 12 bits ADC. The sampled value is then automatically sent to the host via USB.

In S&H mode, when the microcontroller receives a new trigger from ADC_TRG_OUT signal, it sends automatically a "new trigger" frame to host.

Frames automatically sent by BCM-RF to host :

Command (ascii character)	Description	Example
A	BCM-RF's ADC sampled value	A0 : xxxx = 0xxx (1.5 byte)
!	New trigger frame (only available in S&H mode)	!0 : 0000 = 0001

MICROCONTROLLER FIRMWARE (Cont'd)

Write commands send by the host to BCM-RF :

Command Ascii character	Description	Example
D	Write on-board's digital delay line value.	D0 : 0000 = 00xx (1 byte)
E	Write configuration to microcontroller's EEPROM	M0 : 0000 = 0001
I	Set on-board's switches position	I0 : 0000 = 00xx (1 byte)
K	Activate on-board CAL-FO mode (optional)	K0 : 0000 = 0001
M	Activate microcontroller's reverse function algorithm.	M0 : 0000 = 0001
T	Send average number of points for ADC	T0 : 0000 = xxxx (2 bytes)
V	Write calibration constant Qcal (32 bits value) into microcontroller's EEPROM.	V1 : 0000 = xxxx (2 bytes high) V0 : 0000 = xxxx (2 bytes low)
W	Calibration constant Ucal (32 bits value) into microcontroller's EEPROM.	W1 : 0000 = xxxx (2 bytes high) W0 : 0000 = xxxx (2 bytes low)

Read commands send by the host to BCM-RF.

When the BCM-RF's microcontroller receives a read frame, it answers by a response frame which is similar to the write frame.

Command Ascii character	Description	Example
D	Read on-board's digital delay line value.	D0 ?
E	Read microcontroller's configuration stored in EEPROM.	E0 ?
I	Read on-board's switches position	I0 ?
K	Read if on-board's CAL-FO is activated (optional)	K0 ?
M	Read if microcontroller's reverse function algorithm is activated.	M0 ?
T	Read average number of points for ADC	T0 ?
V	Read calibration constant Qcal (32 bits value) into microcontroller's EEPROM.	V1 ? V0 ?
W	Read calibration constant Ucal (32 bits value) into microcontroller's EEPROM.	W1 ? W0 ?

CONNECTORS PINS ALLOCATION

BCM-RF Front panel BNC connectors					
RF-Chassis Rear SMA connectors					
DB9 female connector on BCM-RFC rear panel					
DIN41612M BCM-RF module rear connector					
INPUT SIGNALS					
BCM-RF Input	BCM Input	B8 *		SMA1	
* coaxial insert 1.0/2.3 type					
OUTPUT SIGNALS					
BCM –RF output	BCM Output	B11*		SMA2	
ADC Trigger TTL output pos/neg edge	TRG.ADC.OUT	C19	DB9,6		
BNC front-panel MONITORING					
Input signal after log. demodulation	Signal View				BNC 1
BCM-RF Output	Output View				BNC 2
Hold clock (on rising edge)	Hold View				BNC 3
EXTERNAL TRIGGER INPUT					
External trigger input 50 Ω, pos. edge > 2V	EXT.TRG.IN	B22*		SMA3	
POWER SUPPLY					
+ (8...15) V	+15V	C13			
- (8...15) V	-15V	C15			
Common	COM	C14			

BCM CABLE LAYOUT INSTALLATION

Unnecessary intermediate connectors should be avoided. When for practical reasons patch-panels must be used, the cables on either side of the patch-panel should be passed through tubular ferrite and nanocrystalline cores.

BCM chassis and modules should be kept away as much as possible from RF equipment, klystrons, cavities.

Connectors must be chosen carefully to match the cable used. Connectors manufacturers instructions must be followed meticulously. If cable layout is subcontracted, subcontractors must be informed of the extreme reliability expected from these cables. All cables with connectors must be checked before installation with a network analyzer, up to twice the operating frequency at least.

BCM modules must be installed in an RF-shielded chassis, BCM-RFC/XX or equivalent.

ACCESSORIES

BCM Chassis BCM-RFC/XX

The BCM-RFC/XX chassis is built around a 19" Schroff rackable RF chassis.


Dimensions of the bin: 3U x 84F

Schroff reference: Europac Lab HF/RF #20845-283

The BCM-RFC/XX is available equipped for 1 up to 16 BCM stations. XX being the number of stations.

BCM-RFC/XX with less than 16 stations are partially equipped BCM-RFC/16. As a result, all BCM chassis are field-upgradable to the full 16-station chassis.

Chassis rear view



Chassis front


Unequipped stations are masked with RF-shielded blank panels.

Card Extender BCM-XTD

The card extender allows access to the BCM-RF-E on-board switches while it is connected to the chassis, thus to the readout and control system.

ACCESSORIES (Cont'd)

BCM Chassis BCM-RFC/XX Outer dimensions



SCHEMATICS & BOARD LAYOUT

Schematics and board layouts of our instruments remain the exclusive property of Bergoz Instrumentation at all times. They are protected by the copyright laws.

Schematics and board layouts are not delivered with our instruments. They can be obtained at the specific request of the instrument's user.

A request should be sent by fax, worded in the following way:

To: Bergoz Instrumentation

From: User's name

Date:.....

I am a user of instrument type xxx-xxx serial nr. xxx,xxx,xxx,xxx, etc.

Please send me one copy of the corresponding schematics and board layout.

I will use it for the instrument's maintenance only.

I will make copies only for my own use.

I will inform others who need these schematics that they should request them from Bergoz Instrumentation.

Signed:

ACKNOWLEDGEMENTS

The fundamental principles of our BCM-RF-E module and Turbo-ICT were developed internally at Bergoz Instrumentation.

Our BCM-RF-E module was designed by Sebastien Artinian,. Some functions are borrowed from the earlier LR-BPM design by Alexander Kalinin.

Saint Genis Pouilly, December 2012