

# CSPAD Geometry Software (depricated)

 Unknown macro: 'html'

- Package PSCalib
  - Class GeometryAccess
  - Class GeometryObject
  - Class CalibFileFinder
  - Class CSPadCalibPars
  - Class CSPad2x2CalibPars
- Package CSPadPixCoords
  - Class Image2D
  - GlobalMethods
  - CSPadConfigPars
  - CSPad2x2ConfigPars
  - Class QuadParameters
  - Class PixCoords2x1
  - Class PixCoords2x1V2
  - Class PixCoordsCSPad2x2
  - Class PixCoordsCSPad2x2V2
  - Class PixCoordsQuad
  - Class PixCoordsCSPad
  - Class PixCoordsCSPadV2
- PSANA Modules
  - Module PixCoordsTest
  - Module CSPadImageProducer
  - Module CSPadInterpolImageProducer
  - Module CSPadImageGetTest
  - Module CSPadNDArrProducer
  - Module CSPad2x2NDArrProducer
  - Module CSPad2x2ImageProducer
- Interface description
  - Classes for CSPAD
  - Classes for CSPAD2x2
- References

 Unknown macro: 'html'

A few packages/classes dedicated to the CSPAD geometry are briefly described in this page with references to more detailed documentation.

## Package PSCalib

Package PSCalib provides an access to the calibration parameters.  
For complete reference select [Doxygen documentation](#).

### Class GeometryAccess

[Interface description in doxygen](#)

Class supports universal detector geometry description. More detailed documentation is available in the [Detector Geometry](#) page.

### Class GeometryObject

[Interface description in doxygen](#)

Class defines the elementary building block for hierarchical geometry description. More detailed documentation is available in the [Detector Geometry](#) page.

### Class CalibFileFinder

[Detailed description in Doxygen](#)

CalibFileFinder class finds the path to calibration file. Class constructor creates the object with elements of the path to the calibration file. When all input parameters are provided, the method findCalibFile(...) returns the path/name to the file with requested type of calibration parameters for indicated run number.

### Class CSPadCalibPars

[Detailed description in Doxygen](#)

CSPadCalibPars class loads/holds/provides access to the CSPad geometry calibration parameters. Class constructor creates an object which loads and holds the calibration parameters of all CSPad geometry types. Access methods allows to extract any of requested parameters.  
The CSPad geometry types are defined in the package `pdscalibdata` by the classes

- CalibParsCenterV1
- CalibParsCenterCorrV1
- CalibParsMargGapShiftV1
- CalibParsOffsetV1
- CalibParsOffsetCorrV1
- CalibParsRotationV1
- CalibParsTiltV1
- CalibParsQuadRotationV1
- CalibParsQuadTiltV1

CSPad position w.r.t. IP is defined by the classes:

- CsPadBeamVectorV1
- CsPadBeamIntersectV1

Alternative parameters for optical measurements of entire detector:

- CsPadCenterGlobalV1
- CsPadRotationGlobalV1

Calibration files location and structure is explained in [CsPad calibration in translator](#).

## Class CSPad2x2CalibPars

[Detailed description in Doxygen](#)

CSPad2x2CalibPars class loads/holds/provides access to the CSPad2x2 geometry calibration parameters. Class constructor creates an object which loads and holds the calibration parameters of all CSPad2x2 geometry types. Access methods allows to extract any of requested parameters.  
The CSPad2x2 geometry types are defined in the package `pdscalibdata` by the classes

- CsPad2x2CenterV1
- CsPad2x2TiltV1

## Package CSPadPixCoords

Package CSPadPixCoords pre-calculates the 2x1 section, quad, and CSPad pixel coordinates. CSPad image is also produced in combination of the pixel coordinates with event data.

For complete reference select [Doxygen documentation](#).

## Class Image2D

[Detailed description in Doxygen](#)

Image2D class creates the object with a pointer to the 2D array and its sizes. It is assumed that the 2D array contains image data. Methods of this class provide access to the 2D array and its transformed versions after up-down, left-right flips, transpose, and rotations by n\*90 degree.

## GlobalMethods

[Detailed description in Doxygen](#)

## CSPadConfigPars

[Detailed description in Doxygen](#)

## CSPad2x2ConfigPars

[Detailed description in Doxygen](#)

## Class QuadParameters

QuadParameters class object contains current parameters of the CSPad quad.

## Class PixCoords2x1

`PixCoords2x1` class defines the 2x1 section pixel coordinates in its local frame. Class constructor calculates the X, Y, and Z pixel coordinates for single 2x1 in its own (matrix style) frame; X coordinate is directed along ascending row index from 0 to 388 (from top to bottom for 0 rotation angle), Y coordinate is directed along descending (ONLINE specific...) column index from 184 to 0 (from left to right for 0 rotation angle). The `getPixCoorRot(...)` methods return these coordinates for rotated by N\*90 degree 2x1 in um(micrometer) or pixels. Rows and columns are defined like described in ONLINE modules: /reg/g/psdm/sw/external/lusi-xtc/2.12.0a/x86\_64-rhel5-gcc41-opt/pdsdata/cspad/ElementIterator.hh, Detector.hh  
No parameters needed in constructor; everything is defined through the known 2x1 chip geometry.

## Class `PixCoords2x1V2`

[Detailed description in Doxygen](#)

## Class `PixCoordsCSPad2x2`

[Detailed description in Doxygen](#)

## Class `PixCoordsCSPad2x2V2`

[Detailed description in Doxygen](#)

## Class `PixCoordsQuad`

`PixCoordsQuad` class defines the quad pixel coordinates in its local frame. The coordinate frame coincides with that of optical measurement, but in "matrix style" geometry:

- X axis directed along ascending row index (from top to bottom)
  - Y axis directed along ascending column index (from left to right)
- Pixel coordinates in quad are calculated using precise 2x1 section geometry from the class `PixCoords2x1` and the calibration parameters extracted by the `PSCalib` package.

## Class `PixCoordsCSPad`

`PixCoordsCSPad` class defines and provides access to the `CSPad` pixel coordinates. Class constructor fills the arrays of the `CSPad` pixel coordinates using the `PixCoordsQuad` class object and the calibration parameters extracted by the `PSCalib` package. Access methods return the `CSPad` pixel coordinate for indicated axis, quad, section, row, and column indexes. Pixel coordinates can be returned in um(micrometers), pix(pixels), or int(integer pixel index).

## Class `PixCoordsCSPadV2`

[Detailed description in Doxygen](#) - supports calibration type `center_global`

## PSANA Modules

All `psana` modules are described in [psana - Module Catalog](#).

## Module `PixCoordsTest`

See description in [psana - Module Catalog](#).

To execute this test module go to your release directory and run `psana` with appropriate configuration file:

```
cd <your-favorite-psana-release-directory>
psana -c CSPadPixCoords/data/psana_pixcoords_test.cfg
```

At the end of this test a few text files with images will be created:

`test_2x1.txt`, `test_q0.txt`, ..., `test_q3.txt`, and `test_cspad.txt`

To plot these images you may use the python script, for example:

```
./CSPadPixCoords/data/PlotCameraImageFromFile.py test_q0.txt
...
./CSPadPixCoords/data/PlotCameraImageFromFile.py test_2x1.txt
./CSPadPixCoords/data/PlotCameraImageFromFile.py test_cspad.txt
```

....

..

## Module `CSPadImageProducer`

See description in [psana - Module Catalog](#).

This module produces the uniform 2-D 1750x1750 image array for entire CSPad. Actual pixel coordinates are replaced by the nearest bin coordinates from this 2-D image array. Interpolation is not used in this case.

To test this module go to your release directory and run psana with appropriate configuration file:

```
cd <your-favorite-psana-release-directory>
psana -c CSPadPixCoords/data/psana_cspad_image_producer.cfg
```

This procedure demonstrates that module is doing something, but does not show any useful results. The next module `CSPadImageGetTest` is intended to demonstrate results.

## Module `CSPadInterpollImageProducer`

See description in [psana - Module Catalog](#).

This module is very similar to the `CSPadImageProducer`, but the amplitude in each bin of the uniform 2-D image array is interpolated using the 4 nearest sensor pixels. In this module we use the 4-node bi-linear interpolation algorithm.

## Module `CSPadImageGetTest`

See description in [psana - Module Catalog](#).

To execute this test go to your release directory and run psana with appropriate configuration file:

```
cd <your-favorite-psana-release-directory>
psana -c CSPadPixCoords/data/psana.cfg
```

At the end of this test the text file `cspad_image_get_test_ev000010.txt` with image will be created.

To plot this image you may use the python script:

```
./CSPadPixCoords/data/PlotCameraImageFromFile.py cspad_image_get_test_ev000010.txt
```

## Module `CSPadNDArrProducer`

[Detailed description in Doxygen](#)

## Module `CSPad2x2NDArrProducer`

[Detailed description in Doxygen](#)

## Module `CSPad2x2ImageProducer`

[Detailed description in Doxygen](#)

# Interface description

## Classes for CSPAD

- `CSPadCalibPars` - access to calibration parameters of types `center`, `tilt`, `rotation`, `offeset`, `center_global`, etc.
- `CSPadConfigPars` - access to the run-time configuration parameters.
- `PixCoords2x1V2` - access to pixel coordinates of the 2x1 section.
- `PixCoordsCSPad` - access to pixel coordinates of the CSPAD detector. It uses `PixCoordsCSPadV2` if the calibration file for `center_global` is available.
- `PixCoordsCSPadV2` - supports pixel coordinates of the CSPAD detector for calibration type `center_global`.
- `CSPadImageProducer` - PSANA module which produces image of the CSPAD detector and preserves it in the event (during run-time).
- `CSPadInterpollImageProducer` - PSANA module which produces interpolated image of the CSPAD detector and preserves it in the event (during run-time).

## Classes for CSPAD2x2

- `CSPad2x2CalibPars` - access to calibration parameters of types `center` and `tilt`.
- `PixCoords2x1V2` - access to pixel coordinates of the 2x1 section.

- [PixCoordsCSPad2x2V2](#) - access to pixel coordinates of the CSPAD2x2 detector.
- [CSPad2x2ImageProducer](#) - PSANA module which produces image of the CSPAD2x2 detector and preserves it in the event (during run-time).

## References

- [Detector Geometry](#)
- [CSPAD pixel coordinates and image producer in Python](#)
- [CSPAD image producer in Python](#)
- [CSPAD2x2 modules in Python \(deprecated\)](#)