

CSPAD Geometry and Alignment



Unknown macro: 'html'


- [Introduction](#)
- [2x1 Sensor Geometry](#)
- [Shield to sensor distance](#)
- [Optical measurement](#)
- [Quality Check Procedure](#)
- [Detector geometry model](#)
- [Alignment parameters from optical measurement](#)
- [Alignment of quads in the detector](#)
- [Calibration store](#)
 - [Calibration type](#)
 - [Archive and History](#)
- [Detector data access software](#)
- [References](#)



Unknown macro: 'html'

Introduction


Cornell SLAC Pixel Array Detector (CSPAD) is an imaging X-ray detector made of silicon sensors (2x1) covering about 20x20cm² surface, as shown in the plot:



Pixel coordinates in 2x1 sensor chip are known with sub-micrometer precision. Construction of the detector allows significant freedom in relative positions of 2x1 sensors. To get precise pixel positions in the detector the 2x1 sensor coordinates needs to be calibrated. In this note we describe geometry of the CSPAD detector, optical and quad alignment procedure, parameters, and software providing access to precise geometry information.

2x1 Sensor Geometry

The 2x1 sensor geometry was tested with microscopic measurement. Two slides from Chris Kenney's presentation shows the pixel sizes:




[The same slides in PDF format.](#)

Important 2x1 features:

- Number of rows x columns = 185 x 388. (In DAQ notation of rows and columns is interchanged...)
- Most of pixels have size 109.92 x 109.92 μm^2 .
- Gap between two ASICS is covered by the two rows of elongated pixels with size 109.92 x 274.8 μm^2 .
- Two versions of sensors have different dimensions between corners, so it is reasonable to define pixel coordinates w.r.t. the sensor center.

Shield to sensor distance




Chart of CXI Camera1 provided by Serge Guillet on 2017-06-12.



Optical measurement

Optical measurement is conducted by Gabriel Blaj. Detector or its quad is installed on microscope table and 3-d coordinates of all 2x1 sensor corners are measured with precision about 8 μm (RMS) in x-y plane. All corners in the measurement are numerated in arbitrary order. It is expected that numeration order should be the same for different measurements. This procedure depends on CSPAD construction;

- For CSPAD with movable quads (i.e. for CXI) optical measurement is done separately for each quad. The numeration of corners is shown in the plot:




The same plots in PDF format: [CSPAD quad metrology](#) and [CSPAD pixel layout in quads](#).


For each quad measurement is started from the point #1 which in assembled detector is closest to the beam. The 1-st point (x,y,z) coordinates are re-set to (0,0,0) in the beginning of measurements. At the end, it is checked that the 1-st point coordinates are reproduced within precision of measurement.



The order of points in optical measurement does not coincide with [numeration of 2x1 in DAQ](#), as shown in the plot (and in [PDF file](#)):



- For CSPAD with fixed quad geometry (i.e. for XPP) optical measurement is done for entire detector. The numeration of corners in this case is shown in the file [XPPMetrologyAnnotated.pdf](#) and in the plot:



The 1-st corner of the 3-rd quad (x,y,z) coordinates are re-set to (0,0,0) in the beginning of measurements. At the end, it is checked that the 1-st point coordinates are reproduced within precision of measurement.

Corner coordinates are measured in micrometers (μm) and are saved in the [xlsx](#) format table, also containing numeration of quads and points. Then, [xlsx](#) format table is converted to the text file format in order to use it in python script.

Example of tables for CXI:

- [Metrology in XLSX](#)
- [Metrology in TXT](#)

Example of tables for XPP:

- [Metrology in XLSX](#)
- [Metrology in TXT](#) Python script converts this table to the table with standard numeration of points in quads:
- [Metrology in standard TXT](#)

Then, text table with "standard" numeration of points in quads is feed to the python script which provides quality check of optical measurement and evaluates the alignment parameters for quads. In the beginning, this script changes the numeration of points adopted in optical measurement to numeration of 2x1 used in DAQ. Further, all calibration parameters are associated with numeration of 2x1 sensors and quads in DAQ.

Quality Check Procedure

For quality check of optical measurement we calculate

S1 - 1st short side length of 2x1

S2 - 2nd short side length of 2x1

L1 - 1st long side length of 2x1

L2 - 2nd long side length of 2x1

D1 - 1st diagonal of 2x1 between corners 1 and 3

D2 - 2nd diagonal of 2x1 between corners 2 and 4

dS and dL are the deviations of the 1st and 2nd corner along the short and long sides, respectively. The sign of all dS are chosen in order to provide correct sign for the tilt angle (the same direction for all 2x1 sensors).

<dS/L> - the tilt angle of 2x1 averaged over two sides in radians.

angle(deg) - the same angle in degrees.

dD = D1 - D2

d(dS) = dS1 - dS2

d(dL) = dL1 - dL2

dz3(um) - signed distance from 2x1 sensor plane and corner 3, where the 2x1 sensor plane contains the corner points p1, p2, and p4. This plane is defined by the vectors v21=p2-p1, v41=p4-p1, and their orthogonal vector

```
vort = [v21 x v41].
```

Scalar product with normalization defines the distance from point 3 to the 2x1 plane containing 3 other points:

```
dz3 = (v31 * vort) / |vort|.
```

Quality check parameters expected for perfect geometry:

S1=S2, L1=L2 - the 2x1 sides should have equal length and width,

D1=D2 - the 2x1 diagonals should be equal,

dS1 = dS2 ? (388/185)*dL1 = (388/185)*dL2 - tilt angle should provide consistent deviation for all corners,

dD=0, d(dS)=0, and d(dL)=0 - within precision of measurement.

dz3(um) = 0

Everything, excluding <dS/L> and angle(deg), are in micrometers.

Example of the table with quality check results:

| pair: D2 | S1 dD | S2 d(dS) | dS1 d(dL) | dS2 dz3(um) | L1 | L2 | dL1 | dL2 | <dS/L> | angle(deg) | D1 |
|-------------|----------|-------------|--------------|----------------|-------|-------|------|------|----------|------------|-------|
| Quad 0 | | | | | | | | | | | |
| pair: 0 | 20891 | 20913 | 200 | 222 | 43539 | 43541 | -102 | -100 | 0.00485 | 0.27766 | 48298 |
| 48297 | 1 | -22 | -2 | 2.981 | | | | | | | |
| pair: 1 | 20910 | 20894 | 293 | 277 | 43540 | 43535 | -127 | -132 | 0.00655 | 0.37506 | 48302 |
| 48289 | 13 | 16 | 5 | -23.986 | | | | | | | |
| pair: 2 | 20890 | 20906 | 99 | 83 | 43536 | 43536 | 42 | 42 | 0.00209 | 0.11976 | 48290 |
| 48293 | -3 | 16 | 0 | -3.034 | | | | | | | |
| pair: 3 | 20897 | 20895 | 131 | 133 | 43545 | 43543 | 65 | 63 | 0.00303 | 0.17369 | 48299 |
| 48297 | 2 | -2 | 2 | 6.003 | | | | | | | |
| pair: 4 | 20911 | 20896 | -30 | -45 | 43549 | 43547 | 17 | 15 | -0.00086 | -0.04934 | 48303 |
| 48306 | -3 | 15 | 2 | -5.994 | | | | | | | |
| pair: 5 | 20901 | 20898 | 10 | 7 | 43540 | 43544 | -8 | -4 | 0.00020 | 0.01119 | 48296 |
| 48299 | -3 | 3 | -4 | 9.993 | | | | | | | |
| pair: 6 | 20904 | 20903 | 104 | 105 | 43536 | 43540 | 55 | 59 | 0.00240 | 0.13752 | 48302 |
| 48290 | 12 | -1 | -4 | 52.002 | | | | | | | |
| pair: 7 | 20901 | 20901 | -7 | -7 | 43545 | 43543 | -3 | -5 | -0.00016 | -0.00921 | 48299 |
| 48301 | -2 | 0 | 2 | 14.001 | | | | | | | |
| Quad 1 | | | | | | | | | | | |
| pair: 0 | 20913 | 20914 | -343 | -342 | 43540 | 43550 | 165 | 175 | -0.00787 | -0.45066 | 48313 |
| 48303 | 10 | -1 | -10 | -24.002 | | | | | | | |
| pair: 1 | 20898 | 20901 | -145 | -142 | 43548 | 43551 | 62 | 65 | -0.00330 | -0.18880 | 48300 |
| 48309 | -9 | -3 | -3 | -23.005 | | | | | | | |
| pair: 2 | 20895 | 20903 | -151 | -159 | 43535 | 43532 | -74 | -77 | -0.00356 | -0.20400 | 48289 |
| 48291 | -2 | 8 | 3 | -17.995 | | | | | | | |

| | | | | | | | | | | | |
|---------|-------|-------|------|---------|-------|-------|------|------|----------|----------|-------|
| pair: 3 | 20872 | 20909 | -235 | -272 | 43341 | 43354 | -37 | -24 | -0.00585 | -0.33507 | 48201 |
| 48036 | 165 | 37 | -13 | -13.010 | | | | | | | |
| pair: 4 | 20940 | 20904 | -455 | -491 | 43527 | 43554 | 214 | 241 | -0.01086 | -0.62242 | 48309 |
| 48309 | 0 | 36 | -27 | 1.101 | | | | | | | |
| pair: 5 | 20910 | 20903 | -302 | -309 | 43546 | 43546 | 145 | 145 | -0.00702 | -0.40196 | 48304 |
| 48307 | -3 | 7 | 0 | 6.016 | | | | | | | |
| pair: 6 | 20901 | 20919 | -421 | -439 | 43529 | 43539 | -213 | -203 | -0.00988 | -0.56593 | 48296 |
| 48298 | -2 | 18 | -10 | -8.026 | | | | | | | |
| pair: 7 | 20907 | 20907 | -452 | -452 | 43548 | 43539 | -201 | -210 | -0.01038 | -0.59475 | 48315 |
| 48294 | 21 | 0 | 9 | -8.982 | | | | | | | |
| | | | | | | | | | | | |
| Quad 2 | | | | | | | | | | | |
| pair: 0 | 20914 | 20914 | -25 | -25 | 43536 | 43540 | 10 | 14 | -0.00057 | -0.03290 | 48300 |
| 48300 | 0 | 0 | -4 | -11.013 | | | | | | | |
| pair: 1 | 20901 | 20897 | 7 | 3 | 43546 | 43536 | -1 | -11 | 0.00011 | 0.00658 | 48293 |
| 48300 | -7 | 4 | 10 | 4.036 | | | | | | | |
| pair: 2 | 20899 | 20903 | -256 | -260 | 43533 | 43539 | -127 | -121 | -0.00593 | -0.33954 | 48293 |
| 48294 | -1 | 4 | -6 | -1.023 | | | | | | | |
| pair: 3 | 20912 | 20904 | -210 | -202 | 43540 | 43547 | -106 | -99 | -0.00473 | -0.27106 | 48300 |
| 48306 | -6 | -8 | -7 | 24.004 | | | | | | | |
| pair: 4 | 20910 | 20903 | -543 | -550 | 43535 | 43536 | 261 | 262 | -0.01255 | -0.71923 | 48298 |
| 48299 | -1 | 7 | -1 | 0.004 | | | | | | | |
| pair: 5 | 20904 | 20905 | -241 | -240 | 43538 | 43544 | 111 | 117 | -0.00552 | -0.31647 | 48298 |
| 48301 | -3 | -1 | -6 | -6.024 | | | | | | | |
| pair: 6 | 20903 | 20902 | 21 | 22 | 43539 | 43543 | 8 | 12 | 0.00049 | 0.02829 | 48298 |
| 48298 | 0 | -1 | -4 | 8.999 | | | | | | | |
| pair: 7 | 20902 | 20903 | 82 | 81 | 43546 | 43547 | 35 | 36 | 0.00187 | 0.10723 | 48300 |
| 48306 | -6 | 1 | -1 | 9.995 | | | | | | | |
| | | | | | | | | | | | |
| Quad 3 | | | | | | | | | | | |
| pair: 0 | 20902 | 20898 | -82 | -86 | 43536 | 43543 | 30 | 37 | -0.00193 | -0.11054 | 48289 |
| 48302 | -13 | 4 | -7 | 1.994 | | | | | | | |
| pair: 1 | 20900 | 20904 | 79 | 83 | 43548 | 43541 | -35 | -42 | 0.00186 | 0.10658 | 48301 |
| 48301 | 0 | -4 | 7 | -17.993 | | | | | | | |
| pair: 2 | 20912 | 20894 | 181 | 199 | 43536 | 43535 | 97 | 96 | 0.00436 | 0.25005 | 48298 |
| 48289 | 9 | -18 | 1 | 10.011 | | | | | | | |
| pair: 3 | 20912 | 20905 | 119 | 126 | 43539 | 43538 | 57 | 56 | 0.00281 | 0.16121 | 48296 |
| 48301 | -5 | -7 | 1 | -16.000 | | | | | | | |
| pair: 4 | 20894 | 20912 | -454 | -436 | 43534 | 43545 | 212 | 223 | -0.01022 | -0.58560 | 48303 |
| 48296 | 7 | -18 | -11 | 2.023 | | | | | | | |
| pair: 5 | 20906 | 20919 | -336 | -323 | 43527 | 43535 | 155 | 163 | -0.00757 | -0.43369 | 48295 |
| 48294 | 1 | -13 | -8 | 5.993 | | | | | | | |
| pair: 6 | 20902 | 20905 | -203 | -206 | 43537 | 43525 | -89 | -101 | -0.00470 | -0.26916 | 48293 |
| 48287 | 6 | 3 | 12 | 2.981 | | | | | | | |
| pair: 7 | 20900 | 20897 | -140 | -137 | 43539 | 43544 | -68 | -63 | -0.00318 | -0.18225 | 48298 |
| 48296 | 2 | -3 | -5 | 29.997 | | | | | | | |

This quality check works well to catch significant typos in input table. In case of obvious typos input table can be corrected. When the quality check is passed successfully the alignment parameters are saved and deployed under the `calib`.

Detector geometry model

Since 2014 we support universal detector geometry software which is documented in the [Detector Geometry](#) page and in [CSPAD-geometry-parameters.pdf](#)


Alignment parameters from optical measurement

From optical measurement we extract coordinates of the center of each 2x1 sensor and its tilt angle. The center coordinates are evaluated as an averaged over 4 corners measurements for each axis.

The `tilt` parameters are used along with `rotation` to completely define orientation of 2x1 in quad (for CXI) or in detector (for XPP).

Alignment of quads in the detector

For CSPad with fixed quad geometry (like in XPP) optical measurement of entire detector (should) produces complete information for geometry alignment. For CSPad with moveable quads (like in CXI) quads relative position needs to be adjusted through the alignment parameters for quads. It is usually done using typical images with diffraction rings, wires or other shading objects:



Although few algorithms of automatic quad alignment were tried, we did not find good generic way for automated quad tuning. Currently, the quad tuning parameters in `marg_gap_shift` and `offset_corr` are adjusted manually for runs with specific images.

Calibration store

The official space for CSPAD alignment parameters is

`/reg/d/psdm/<INSTRUMENT>/<experiment>/calib/CsPad::Calib<VERSION>/<CSPad-name>/<type>/<run-range>.data`

For example:

```
/reg/d/psdm/CXI/cxi80410/calib/CsPad::CalibV1/CxiDs1.0:Cspad.0/geometry/1142-end.data
```

The file name consists of the run range followed by the `.data` extension, for example, `0-end.data`, `11-end.data`, `47-52.data`, etc.

Calibration type

Detector geometry calibration information is located in a single file of type

- `geometry` - contains hierarchical description of all detector components; for example for CSPAD, sensors' location and rotation in the quads, quads - in the detector, detector - in the setup, etc.

Archive and History

Optical measurement and other alignment files can be found in

- `/reg/g/psdm/detector/alignment/cspad/`
- [Geometry History](#)

Detector data access software

- [Detector data access software](#)

References

- [CSPAD Geometry and Alignment - Deprecated](#) - old version of this page
- [Detector Geometry](#) - confluence page
- [CSPAD in DAQ](#) - schematic description of CSPAD geometry available in DAQ.
- [CSPAD quad metrology](#) - slides for CXI type CSPAD quads
- [CSPad pixel layout in quads](#) - pdf file with numeration of ASICs in the CSPAD quads
- [XPPMetrologyAnnotated.pdf](#) - order of measurements of XPP camera.
- [Geometry History](#) - page with references to calibration files.
- [Detector data access software](#) - auto-generated documentation of the Detector package.