

Note on XCS code from Marcin

Content

- Content
- Get data from xtc files
 - Scripts for pre-processing in pyana
 - Content of pyana.cfg
 - import_data.py
 - Essential code
 - Sequence of operations on pixels
- Evaluate correlators
 - Scripts
 - Content of run_correlator.py
 - Essential code of xcs_correlator.py

References

Get data from xtc files

Scripts for pre-processing in pyana

Pyana configuration file and module:

- /reg/neh/home1/sikorski/xcs_pyana_current/pyana.cfg
- /reg/neh/home1/sikorski/xcs_pyana_current/xcs_timepix_pkg/src/import_data.py

Main features of import_data.py:

splits image for 10x10 parts shaped as (130,134), and saves them in sub-directory like

/reg/neh/home1/sikorski/xcs_pyana_current/e167-r0015-s00-c00/2013-04-03-10-39-22-734268/

Content of pyana.cfg

on 2013-10-14

```
[pyana]
modules = xcs_timepix_pkg.import_data
[xcs_timepix_pkg.import_data]
detector = 1      <== 1 = PI, 2 = TimePix
dark_fina = /reg/neh/home1/sikorski/xcs_pyana_current/e167-r0020-s00-c00/2013-04-03-09-43-06-939033/e167-r0020-
s00-c00_dark_img.txt
flag_bin = 1
bin_x = 1
bin_y = 1
LLD = 20.0
s_x = 10
s_y = 10
saxs_interval = 100
not_zero = 1.0E-9
r1 = 0
c1 = 0
r2 = 1300
c2 = 1340
saturation = 65535
```

import_data.py

Essential code

```

def __init__ ( self, ... )
    self.frame_size = ... (1300,1340)
    self.step_x = ... 130
    self.step_y = ... 134
    self.LLD      = ... 20.
    self.not_zero = 1.0E-9
    self.saturation = 65535
    self.dark = array with pedestals

    self.avg_img = np.zeros(self.frame_size)
    self.saturated_pixels = np.ones((1300,1340))

def beginjob( self, evt, env ) :
    ...
    for i in range(0,self.frame_size[0],self.step_x):
        for j in range(0,self.frame_size[1],self.step_y):
            fina = path + '/' +fina_gen + '_data_' + str(i) + '_' +str(j) + '.npy'
            self.f_id.append(open(fina, 'w'))

def event( self, evt, env ) :
    ...
    PI = evt.get(xtc.TypeId.Type.Id_PrincetonFrame)
    frame = np.array(PI.data()[self.R1:self.R2,self.C1:self.C2], dtype = np.float64)
    a = np.where(frame >= self.saturation)
    self.saturated_pixels[a] = 0
    np.subtract(frame, self.dark, frame)
    frame[frame < self.LLD] = self.not_zero

    np.add(self.avg_img, frame , self.avg_img)

    counter = 0
    for i in range(0,self.frame_size[0],self.step_x):
        for j in range(0,self.frame_size[1],self.step_y):
            np.save(self.f_id[counter], frame[i:i+self.step_x,j:j+self.step_y])
            counter += 1

def endjob( self, evt, env ) :
    ...
    counter = 0
    self.avg_img = self.avg_img/self.calibcycle_counter
    for i in range(0,self.frame_size[0],self.step_x):
        for j in range(0,self.frame_size[1],self.step_y):
            self.f_id[counter].close()
            counter += 1

    np.savetxt(self.f_id_sat, self.saturated_pixels,fmt='%-15f')
    np.savetxt(self.f_id_avg,self.avg_img, fmt='%-15f')

```

Sequence of operations on pixels

1. Get data frame in np.float64: frame = np.array(PI.data()[self.R1:self.R2,self.C1:self.C2], dtype = np.float64)
2. Set saturated pixels with intensity >=65535 to 0 a = np.where(frame >= self.saturation); self.saturated_pixels[a] = 0
3. Subtract pedestals: np.subtract(frame, self.dark, frame)
4. Apply LLD: frame[frame < self.LLD] = self.not_zero
5. Accumulate frame in avg_img: np.add(self.avg_img, frame , self.avg_img)
6. Accumulate split frame in files:
np.savetxt(self.f_id[counter], frame[i:i+self.step_x,j:j+self.step_y]) #, fmt='%-20.10f')

Evaluate correlators

Scripts

```
/reg/neh/home1/sikorski/scripts/run_correlator.py  
/reg/neh/home1/sikorski/scripts/correlator/xcs_correlator.py
```

Use split image files from

- in 2013-04-03 test: /reg/neh/home1/sikorski/xcs_pyana_current/e167-r0015-s00-c00/2013-04-03-10-39-22-734268/
- in 2013-08-23 test: /reg/neh/home1/sikorski/xcs_pyana_current/e167-r0011-s00-c00/2013-08-23-11-30-43-581071/

Content of run_correlator.py

on 2013-10-14

```
import sys  
sys.path.append("/reg/neh/home/sikorski/scripts/correlator")  
import xcs_correlator as xcs_corr  
#####  
path = r"/reg/neh/home1/sikorski/xcs_pyana_current/e167-r0011-s00-c00/2013-08-23-11-30-43-581071/  
path_mask = r"/reg/neh/home/sikorski/scripts"  
fina_mask = r"test.mat.npz"  
ts_fina = r"e167-r0011-s00-c00_time_stamps.txt"  
monitor_number = -1  
threshold_low = -0.002  
threshold_high = 0.1E9  
LLD = 20  
tau = [1,2,3,4,5,6,7,8,9,10,12,14,16,20]  
##### | Define parameters in xcs_correlator.py  
c = xcs_corr.img_correlator(path) | self.path = path; self.notzero = 1; self.results = {};  
c.set_time_stamps_fina(ts_fina) | self.time_stamps_fina = "e167-r0011-s00-c00_time_stamps.txt"  
c.initialize_correlator() | defines many self.*  
c.set_monitor(monitor_number) | self.monitor_name = 'no normalization'; self.monitor = np.ones_like  
(monitor[:,0]) = all 1  
c.set_time_stamps_fina(ts_fina) | !!! repetition...  
c.set_LLD(LLD) | self.LLD = 20  
c.set_tau(tau) | self.tau = [1,2,3,4,5,6,7,8,9,10,12,14,16,20]  
c.set_mask_folder(path_mask) | self.mask_folder = "/reg/neh/home/sikorski/scripts"  
c.set_mask_file(fina_mask) | self.mask_fina = "test.mat.npz"  
c.set_mask() | loads a few files in dictionary self.mask = {}  
c.set_threshold(threshold_low) | self.threshold = -0.002 (def=0)  
c.set_threshold_up(threshold_high) | self.threshold_up = 0.1E9 (def=1E4)  
c.submit_correlator_jobs() | calculates correlators  
c.avg_g2()  
c.harvest_SAXS_data()  
c.plot_g2()  
c.plot_fit_results(777)  
c.plot_SAXS()  
plt.show()
```

Essential code of xcs_correlator.py

on 2013-10-14

```
class img_correlator:  
    def __init__(self, path):  
        self.path = path  
        self.notzero = 1  
        self.results = {}  
        #self.LLD = 20 # it was 200 in the past  
  
    def submit_correlator_jobs(self):  
        filenames = os.listdir(self.path)  
        for filename in filenames:  
  
            self.correlator_2D(filename)
```

```

def correlator_2D(self, fina):

    f_id = open(os.path.join(self.path, fina), 'r')
    data = np.zeros((self.ncols, self.nrows, self.nframes)) # (130,134,500)

    for i in range(self.nframes):
        data[:, :, i] = np.load(f_id)
        data[:, :, i] = data[:, :, i]/ self.monitor[i]
        self.log('Frame ' + str(i) + ' was normalized by ' + str(self.monitor[i]))
        data[:, :, i][data[:, :, i]<(self.LLD)] = self.notzero

    f_id.close()
    s = data.shape
    self.tau = ... tau = [1,2,3,4,5,6,7,8,9,10,12,14,16,20]
    ...

    IP = np.zeros(s,np.float64)
    IF = np.zeros_like(IP)
    G2 = np.zeros_like(IP)
    counter = np.zeros_like(self.tau)

    for t in range(len(self.tau)):
        delta = self.tau[t]
        for i in range(data.shape[2]-delta):

            if (self.monitor[i] >= self.threshold) and (self.monitor[i + delta] >= self.threshold) and (self.monitor[i] <= self.threshold_up) and (self.monitor[i + delta] <= self.threshold_up):

                np.add(G2[:, :, t],np.multiply(data[:, :, i], data[:, :, i+delta]),G2[:, :, t])
                np.add(IP[:, :, t],data[:, :, i],IP[:, :, t])
                np.add(IF[:, :, t],data[:, :, i+delta],IF[:, :, t])
                counter[t] += 1

        ...

        to_be_deleted = []
        for i in range(0,len(self.tau)):
            if counter[i] == 0:
                to_be_deleted.append(i)
        self.tau = np.delete(self.tau,to_be_deleted)
        print "tbd = ", to_be_deleted
        counter = np.delete(counter, to_be_deleted)
        G2 = np.delete(G2,to_be_deleted, axis = 2)
        IP = np.delete(IP,to_be_deleted, axis = 2)
        IF = np.delete(IF,to_be_deleted, axis = 2)
        for i in range(0,len(self.tau)):
            G2[:, :, i] = G2[:, :, i]/counter[i]
            IP[:, :, i] = IP[:, :, i]/counter[i]
            IF[:, :, i] = IF[:, :, i]/counter[i]
            IP[:, :, i][IP[:, :, i] < self.notzero] = self.notzero
            IF[:, :, i][IF[:, :, i] < self.notzero] = self.notzero
        ...

        f_id = open(path_G2, 'w')
        for i in range(G2.shape[2]):
            np.save(f_id, G2[:, :, i])
        f_id.close()

        f_id = open(path_IP, 'w')
        for i in range(IP.shape[2]):
            np.save(f_id, IP[:, :, i])
        f_id.close()

        f_id = open(path_IF, 'w')
        for i in range(IF.shape[2]):
            np.save(f_id, IF[:, :, i])
        f_id.close()

```

References

- [Algorithms for Time Correlation Experiments](#)
- [Weekly progress of the IDPE project for TCE](#)