

# Fetching data from the database

- Add package (if necessary)
- Simple test
- Method detectors
- Method experiment\_runs
- Method run\_attributes
- Method calibration\_runs
- Igor's script to find data type

This note is composed from [Igor's e-mail from 2013-08-14](#) and further e-mail exchange.

## Add package (if necessary)

```
addpkg RegDB V00-01-16  
scons
```

## Simple test

All methods can be tested in a single command:

```
% python RegDB/src/experiment_info.py
```

## Method detectors

From Python code:

```
from RegDB import experiment_info  
detectors = experiment_info.detectors('XPP','xppa4513',1)  
print detectors  
  
['BldEb-0|NoDevice-0',  
'EpicsArch-0|NoDevice-0',  
'NoDetector-0|Evr-0',  
'XppEndstation-0|Opal1000-2',  
'XppGon-0|Cspad-0',  
'XppSb4Pim-1|Tm6740-1']
```

## Method experiment\_runs

To get info which runs (and WHAT kind of runs) exist for the specified experiment do:

```

runs = experiment_info.experiment_runs('XPP','xppa4513')
print runs

[{'begin_time': 1375417636042155759L,
 'begin_time_unix': 1375417636,
 'end_time': 1375417646535192694L,
 'end_time_unix': 1375417646,
 'exper_id': 329L,
 'id': 69762L,
 'num': 1L,
 'type': 'DATA'},
 {'begin_time': 1375437784068608263L,
 'begin_time_unix': 1375437784,
 'end_time': 1375437816700510685L,
 'end_time_unix': 1375437816,
 'exper_id': 329L,
 'id': 69763L,
 'num': 2L,
 'type': 'DATA'},
 ...
 {'begin_time': 1375806447772306081L,
 'begin_time_unix': 1375806447,
 'end_time': 1375806491980944918L,
 'end_time_unix': 1375806491,
 'exper_id': 329L,
 'id': 70074L,
 'num': 189L,
 'type': 'DATA'}]

```

The run number is given by key:

```
'num': 189L
```

The last key of the dictionary:

```
'type': 'DATA'
```

That is for normal data runs. A guess the idea is that people were able to change that to specify:

```
'type': 'CALIB'
```

This is already supported by I don't think we have any experiment where they're using it.

## Method `run_attributes`

To access tags associated with run use method `run_attributes(instr_name, exper_name, runnum)`:

```

instr_name = 'CXI'
exper_name = 'cxic0213'
runnum     = 215

for attr in experiment_info.run_attributes(instr_name, exper_name, runnum):
    print 'class:',attr['class']
    print 'name:',attr['name']
    print 'type (of the value):',attr['type']
    print 'value (optional):',attr['val']
    print 'description (optional):',attr['descr']

```



The value of an attribute is optional. If no value is set then you'll see None. The default description is the empty string

## Method calibration\_runs

Implemented since RegDB V00-01-18.

```
runs = calibration_runs(instr_name, exper_name)
for run in runs:
    print 'comment for this run:',run['comment']
    for calibtype in run['calibrations']:
        print calibtype
```

The function returns the dictionary of pairs:

```
{runnum:dict_recs}
```

where dict\_recs is a dictionary of records:

```
{'comment':<string-of-comment>, 'calibrations':[<list-of-calibrations>]}
```

and [<list-of-calibrations>] consists of calibration types, for which this file can be applied.

For example:

```
{runnum:dict_recs}
{
    123 : {
        'calibrations' : ['dark', 'flat', ...] ,
        'comment'      : 'Phillip marked this as DARK and FLAT filed run'
    },
    200 : {
        ...
    }
}
```

## Igor's script to find data type

```
#!/usr/bin/env python
import sys
import MySQLdb as db
import RegDB.experiment_info as expinfo
if len(sys.argv) != 2:
    print "usage: <text-to-search>"
    sys.exit(1)
text2search = sys.argv[1]
conn = db.Connect(host='psdb',db='regdb',user='regdb_reader')
cursor = conn.cursor()
cursor.execute('SELECT i.name AS `instr`, e.name as `exper` FROM instrument `i`, experiment `e` WHERE i.id=e.instr_id')
for (instr,exper) in cursor.fetchall():
    runs = expinfo.experiment_runs(instr,exper)
    if len(runs):
        print "%s:%s (%d runs)" % (instr,exper,len(runs))
        first = runs[0]['num']
        last = runs[-1]['num']
        for runnum in range(first,last+1):
            try:
                for det in expinfo.detectors(instr,exper,runnum):
                    idx = det.find(text2search)
                    if idx != -1:
                        print "%s:%s:%d %s" % (instr,exper,runnum,det)
            except ValueError:
                pass
```