

psana - Module Examples



Unknown macro: 'html'

- [About](#)
- [Auxiliary scripts](#)
- [Examples for package cspad_mod](#)
 - [Example for module cspad_mod::CsPad2x2Pedestals](#)
- [Examples for package CSPadPixCoords](#)
 - [Example for module CSPadPixCoords::CSPadImageProducer](#)
 - [Example for module CSPadPixCoords::CSPadNDArrProducer](#)
 - [Example for module CSPadPixCoords::CSPad2x2NDArrProducer](#)
 - [Example for module CSPadPixCoords::CSPad2x2NDArrReshape](#)
 - [Example for module CSPadPixCoords::CSPad2x2ImageProducer](#)
- [Example for package ImgPixSpectra](#)
 - [Example for module ImgPixSpectra::CSPadPixSpectra](#)
 - [Example for module ImgPixSpectra::CSPad2x2PixSpectra](#)
- [Examples for package ImgAlgos](#)
 - [Example for module ImgAlgos::Tahometer](#)
 - [Example for module ImgAlgos::EpixNDArrProducer](#)
 - [Example for module ImgAlgos::PnccdImageProducer](#)
 - [Example for module ImgAlgos::CameraImageProducer](#)
 - [Example for module ImgAlgos::PrincetonImageProducer](#)
 - [Example for module ImgAlgos::AcqirisArrProducer](#)
 - [Example for module ImgAlgos::AcqirisAverage](#)
 - [Example for module ImgAlgos::AcqirisCalib](#)
 - [Example for module ImgAlgos::AcqirisCFD](#)
 - [Example for combination of Acqiris modules](#)
 - [Example for module ImgAlgos::NDArrImageProducer](#)
 - [Example for module ImgAlgos::NDArrAverage](#)
 - [Example for module ImgAlgos::NDArrCalib](#)
 - [Example for module ImgAlgos::NDArrDropletFinder](#)
 - [Example for module ImgAlgos::PixCoordsProducer](#)
 - [Example for module ImgAlgos::ImgAverage](#)
 - [Example for module ImgAlgos::ImgMaskEvaluation](#)
 - [Masks: noisy, saturated, and combined:](#)
 - [Example for module ImgAlgos::ImgCalib](#)
 - [Example for module ImgAlgos::ImgRadialCorrection](#)
 - [Example for module ImgAlgos::ImgPeakFinder](#)
 - [Example for module ImgAlgos::ImgPeakFilter](#)
 - [Example for module ImgAlgos::ImgPeakFinderAB](#)
 - [Example for module ImgAlgos::ImgHitFinder](#)
 - [Example for module ImgAlgos::ImgSpectra](#)
 - [Example for module ImgAlgos::ImgSpectraProc](#)
 - [Example for module ImgAlgos::ImgSaveInFile](#)
 - [Example for module ImgPeakFinder and ImgPeakFilter for CSPad](#)
 - [Example for module ImgAlgos::CSPadArrAverage](#)
 - [Example for module ImgAlgos::CSPadBkgdSubtract](#)
 - [Example for Module ImgAlgos::CSPadMaskApply](#)
 - [Example for module ImgAlgos::CSPadArrNoise](#)
 - [Example for Module ImgAlgos::CSPadArrPeakFinder](#)
 - [Example for module ImgAlgos::CSPadArrPeakAnalysis](#)
 - [Example for TimeStampFilter and XtcOutputModule](#)
 - [Example for module ImgAlgos::UsdUsbEncoderFilter](#)
- [Examples for Package pyimgalogs](#)
 - [Example of configuration file for CSPAD](#)
 - [Example of configuration file for CSPAD2x2](#)
- [References](#)



Unknown macro: 'html'

About

This page provides examples for selected modules from [psana - Module Catalog](#).

Auxiliary scripts

A few python scripts in `ImgPixSpectra/data/` show how to process/plot the spectral array stored in the file.

- `PlotSpectralArrayFromFile.py` - allows to plot content of the spectral array as a 2-d plot.

- `SpectralArray.py` - provides access to the spectral array stored in the file. The class `SpectralArray` defined in this script is used in the `PlotSpectralArrayFromFile.py`.
- `MergeArrays.py` - sums the arrays from different files defined in the list and saves resulting array in a single file with the same shape. In this script the list of files is hardwired in the `get_list_of_input_file_names()` method. The output file name, `out_fname`, is also hardwired in the call to `spectra_merging(out_fname)`.

A few auxiliary scripts for example are located in the directory `ImgAlgos/data`:

- `psana.cfg` - set non-default parameters to run `psana` for `ImgAlgos::ImgPeakFinder` and `ImgAlgos::ImgPeakFilter`. The `psana` running this script saves images and peaks for event 115 in text files.
- `PlotCameraImageFromFile.py` - Plots image and spectrum for image saved in file.
- `PlotCameraImageAndPeaks.py` - Plots image with found peaks and spectrum.

Examples for package `cspad_mod`

Example for module `cspad_mod::CsPad2x2Pedestals`

Configuration file for pedestals calibration of two CSPAD2x2 simultaneously:

```
[psana]
files = /reg/d/psdm/xpp/xpptut13/xtc/e308-r0070-s02-c00.xtc \
        /reg/d/psdm/xpp/xpptut13/xtc/e308-r0070-s03-c00.xtc

modules = cspad_mod.CsPad2x2Pedestals:0 \
          cspad_mod.CsPad2x2Pedestals:1

#skip-events = 100
#events = 1111

[cspad_mod.CsPad2x2Pedestals:0]
source = DetInfo(XppGon.0:Cspad2x2.0)
output = pedestals-ave-xppi0513-r0070-Cspad2x2.0.dat
noise = pedestals-rms-xppi0513-r0070-Cspad2x2.0.dat

[cspad_mod.CsPad2x2Pedestals:1]
source = DetInfo(XppGon.0:Cspad2x2.1)
output = pedestals-ave-xppi0513-r0070-Cspad2x2.1.dat
noise = pedestals-rms-xppi0513-r0070-Cspad2x2.1.dat
```

Command to run this script:

```
psana -c psana-xppi0513-r0070-cspad2x2-pedestals.cfg
```

The xtc file name(s) may be passed as an argument:

```
psana -c psana-xppi0513-r0070-cspad2x2-pedestals.cfg /reg/d/psdm/xpp/xpptut13/xtc/e308-r0070-*.xtc
```

Output files contain results, which can be plotted for average values:

-
and rms values:
-

Examples for package `CSPadPixCoords`

Example for module `CSPadPixCoords::CSPadImageProducer`

How to write the CSPad image in text file:

```

[psana]
files      = /reg/d/psdm/<instrument>/<experiment>/xtc/<file-name-1>.xtc
events     = 5

modules    = cspad_mod.CsPadCalib CSPadPixCoords.CSPadImageProducer ImgAlgos.ImgSaveInFile

[cspad_mod.CsPadCalib]
inputKey   =
outputKey  = calibrated
doPedestals = yes
doPixelStatus = no
doCommonMode = no

#[CSPadPixCoords.CSPadInterpolImageProducer]
[CSPadPixCoords.CSPadImageProducer]
calibDir   = /reg/d/psdm/<instrument>/<experiment>/calib
typeGroupName = CsPad::CalibV1
source     = CxiDs1.0:Cspad.0
key        = calibrated
imgkey     = Image2D
tiltIsApplied = true

[ImgAlgos.ImgSaveInFile]
source     = CxiDs1.0:Cspad.0
#eventSave = 1
saveAll    = true

```

See [Package CSPadPixCoords](#)

Example for module CSPadPixCoords::CSPadNDArrProducer

See [Module CSPadPixCoords::CSPadNDArrProducer](#)

Example of the module CSPadNDArrProducer working in sequence with ImgAlgos::NDArrAverage is shown in [Example for module ImgAlgos::NDArrAverage](#)

Example for module CSPadPixCoords::CSPad2x2NDArrProducer

```
# Run this script:
# psana -c psana-meca1113-r0045-cspad-cspad2x2-dark-hotpix.cfg
```

```
[psana]
# Default calibration directory:
# calib-dir = /reg/d/psdm/mec/meca1113/calib

files = exp=meca1113:run=45
events = 400
#skip-events = 0

modules = CSPadPixCoords.CSPad2x2NDArrProducer:1 \
          CSPadPixCoords.CSPad2x2NDArrProducer:2 \
          ImgAlgos.NDArrAverage:1 \
          ImgAlgos.NDArrAverage:2
```

```
[CSPadPixCoords.CSPad2x2NDArrProducer:1]
source = MecTargetChamber.0:Cspad2x2.1
inkey =
outkey = cspad2x2.1_ndarr
outtype = int16
print_bits = 3
```

```
[CSPadPixCoords.CSPad2x2NDArrProducer:2]
source = MecTargetChamber.0:Cspad2x2.2
inkey =
outkey = cspad2x2.2_ndarr
outtype = int16
print_bits = 3
```

```
[ImgAlgos.NDArrAverage:1]
source = MecTargetChamber.0:Cspad2x2.1
key = cspad2x2.1_ndarr
avefile = cspad2x2.1-ave
rmsfile = cspad2x2.1-rms
maskfile = cspad2x2.1-msk
hotpixfile = cspad2x2.1-hot
#evts_stage1 = 100
#gate_width1 = 100.
thr_rms_ADU = 10
thr_min_ADU = 4
thr_max_ADU = 10000
print_bits = 29
```

```
[ImgAlgos.NDArrAverage:2]
source = MecTargetChamber.0:Cspad2x2.2
key = cspad2x2.2_ndarr
avefile = cspad2x2.2-ave
rmsfile = cspad2x2.2-rms
maskfile = cspad2x2.2-msk
hotpixfile = cspad2x2.2-hot
#evts_stage1 = 100
#gate_width1 = 100.
thr_rms_ADU = 10
thr_min_ADU = 4
thr_max_ADU = 10000
print_bits = 29
```

This script makes ndarrays for two cspad2x2 detectors and use them evaluate average, rms, mask and hot pixel arrays.

See [Module CSPadPixCoords::CSPad2x2NDArrProducer](#)

Example for module CSPadPixCoords::CSPad2x2NDArrReshape

Example of configuration file

```
[psana]
#files = exp=meca1113:run=376
#events = 10
##skip-events = 0

modules = cspad_mod.CsPadCalib \
          CSPadPixCoords.CSPad2x2NDArrProducer:clb \
          CSPadPixCoords.CSPad2x2NDArrProducer:raw \
          CSPadPixCoords.CSPad2x2NDArrReshape \
          ImgAlgos.NDArrAverage:clb \
          ImgAlgos.NDArrAverage:raw
#      EventKeys

[cspad_mod.CsPadCalib]
inputKey      =
outputKey     = clb_data
doPedestals   = yes
doPixelStatus = yes
doCommonMode  = yes

[CSPadPixCoords.CSPad2x2NDArrProducer:clb]
source        = MecTargetChamber.0:Cspad2x2.1
inkey         = clb_data
outkey        = cspad2x2.1_clb:as_data
outtype       = int16
print_bits    = 5

[CSPadPixCoords.CSPad2x2NDArrProducer:raw]
source        = MecTargetChamber.0:Cspad2x2.1
inkey         =
outkey        = cspad2x2.1_raw:as_data
outtype       = int16
print_bits    = 5

[CSPadPixCoords.CSPad2x2NDArrReshape]
source        = MecTargetChamber.0:Cspad2x2.1
keys_in       = cspad2x2.1_raw:as_data cspad2x2.1_clb:as_data
print_bits    = 255

[ImgAlgos.NDArrAverage:clb]
source        = MecTargetChamber.0:Cspad2x2.1
key           = cspad2x2.1_clb
avefile       = arr-ave-clb
rmsfile       = arr-rms-clb
print_bits    = 255

[ImgAlgos.NDArrAverage:raw]
source        = MecTargetChamber.0:Cspad2x2.1
key           = cspad2x2.1_raw
avefile       = arr-ave-raw
rmsfile       = arr-rms-raw
print_bits    = 255
```

This script produces raw and calibrated ndarrays shaped as data (185,388,2), then pass these arrays as

```
keys_in      = cspad2x2.1_raw:as_data cspad2x2.1_clb:as_data
```

to the CSPadPixCoords.CSPad2x2NDArrReshape module, which produces two reshaped (2,185,388) arrays and saves them in the event store with default keys

```
cspad2x2.1_raw, cspad2x2.1_clb (suffixes are dropped)
```

then both re-shaped arrays are averaged and saved in files on disk.

Example of associated python script:

```
#!/usr/bin/env python

import sys
import numpy as np
import matplotlib.pyplot as plt
import psana

psana.setConfigFile('psana-meca1113-r0376-cspad2x2-CSPad2x2NDArrReshape.cfg')

dsname = 'exp=meca1113:run=376'
print "Data source: %s" % dsname
ds = psana.DataSource(dsname)

#-----
print "Initializing Matplotlib Plotter"
fig = plt.figure(figsize=(10,5), dpi=80, facecolor='w', edgecolor='w', frameon=True)
plt.ion()
plt.show()

evnum = 0
evnum_max = 50

for evt in ds.events() :

    evtid = evt.get(psana.EventId)

    evnum += 1
    if evnum > evnum_max : break

    img_as_data = evt.get(psana.ndarray_int16_3, psana.Source('DetInfo(MecTargetChamber.0:Cspad2x2.1)'),
    'cspad2x2.1_clb:as_data')
    img_resdpd = evt.get(psana.ndarray_int16_3, psana.Source('DetInfo(MecTargetChamber.0:Cspad2x2.1)'),
    'cspad2x2.1_clb')

    print 'img_as_data.shape = ', img_as_data.shape
    print 'img_resdpd.shape = ', img_resdpd.shape

    img_as_data.shape = (2*185,388)
    img_resdpd.shape = (2*185,388)

    ax1 = fig.add_axes([0.05, 0.06, 0.44, 0.87])
    ax2 = fig.add_axes([0.55, 0.06, 0.44, 0.87])

    fig.canvas.set_window_title('Image from arrays "as-data" and "reshaped"')

    imsh1 = ax1.imshow(img_as_data, interpolation='nearest', aspect='auto', origin='upper') # ,
    extent=img_range)
    imsh2 = ax2.imshow(img_resdpd, interpolation='nearest', aspect='auto', origin='upper') # ,
    extent=img_range)

    plt.title("Event: %d.%d" % evtid.time())
    plt.draw()
    plt.clf()
```

This script retrieves from the event store calibrated ndarrays shaped as data (185,388,2) and re-shaped (2,185,388), then re-shape both to 2-d arrays (2*185,388) and plot them. Re-shaped array looks as recognizable (without spaces between ASICs) cspad2x2 image. Array shapes "as data" looks as an overlap of two images.

See [Module CSPadPixCoords::CSPad2x2NDArrReshape](#)

JIRA issue: [PSAS-45](#) - Getting issue details...

STATUS

Example for module CSPadPixCoords::CSPad2x2ImageProducer

See [Module CSPadPixCoords::CSPad2x2ImageProducer](#)

Example of the configuration script for psana (cspad2x2-test.cfg):

```
[psana]
files      = /reg/d/psdm/<instrument>/<experiment>/xtc/<file-name-1>.xtc
events     = 5
modules    = CSPadPixCoords.CSPad2x2ImageProducer ImgAlgos.ImgSaveInFile

[CSPadPixCoords.CSPad2x2ImageProducer]
source     = DetInfo(:Cspad2x2)
inkey      =
outimgkey  = Image
tiltIsApplied = true
print_bits = 15

[ImgAlgos.ImgSaveInFile]
source     = DetInfo(:Cspad2x2)
key        = Image
fname      = cspad2x2
saveAll    = true
#eventSave = 5
```

Command to run:

```
psana -c cspad2x2-test.cfg
```

One of the saved files cspad2x2-<run>-<timestamp>.txt is plotted as an image by the command

```
./PlotCameraImageFromFile.py cspad2x2-<run>-<timestamp>.txt 0 1200
```

-

Configuration file for cspad2x2 with pedestal subtraction



Access to the CSPad2x2 aligned geometry is added on 2013-02-13
and is available for offline release > ana-0.9.1.

Example of <psana-config-file.cfg>:

```

[psana]
files      = /reg/d/psdm/mec/mec73313/xtc/e268-r0180-s02-c00.xtc
#calib-dir = ./calib

modules    = cspad_mod.CsPadCalib CSPadPixCoords.CSPad2x2ImageProducer ImgAlgos.ImgSaveInFile
events     = 5

[cspad_mod.CsPadCalib]
source     = DetInfo(MecTargetChamber.0:Cspad2x2.3)
inputKey   =
outputKey  = calibrated_arr
doPedestals = yes
doPixelStatus = no
doCommonMode = no

[CSPadPixCoords.CSPad2x2ImageProducer]
calibDir   = /reg/d/psdm/mec/mec73313/calib
typeGroupName = CsPad2x2::CalibV1
#source    = DetInfo(MecTargetChamber.0:Cspad2x2.3)
source     = :Cspad2x2.3
inkey      = calibrated_arr
outimgkey  = Image
tiltIsApplied = true
print_bits = 15

[ImgAlgos.ImgSaveInFile]
source     = DetInfo(MecTargetChamber.0:Cspad2x2.3)
key        = Image
fname      = cspad2x2.3
saveAll    = true
print_bits = 3
#eventSave = 5

```

Example of psana configuration file to get cspad2x2 images for two detectors and save them in files, one in txt, another in tiff formats:

```

[psana]
files      = /reg/d/psdm/xpp/xpptut13/xtc/e308-r0008-s02-c00.xtc \
             /reg/d/psdm/xpp/xpptut13/xtc/e308-r0008-s03-c00.xtc

#modules    = cspad_mod.CsPad2x2Pedestals

#calib-dir  = ./calib
calib-dir   = /reg/d/psdm/xpp/xpptut13/xtc/calib

modules     = cspad_mod.CsPadCalib:0 \
              cspad_mod.CsPadCalib:1 \
              CSPadPixCoords.CSPad2x2ImageProducer:0 \
              CSPadPixCoords.CSPad2x2ImageProducer:1 \
              ImgAlgos.ImgSaveInFile:0 \
              ImgAlgos.ImgSaveInFile:1

events      = 5

[cspad_mod.CsPadCalib:0]
source      = DetInfo(XppGon.0:Cspad2x2.0)
inputKey    =
outputKey   = calibrated_arr0
doPedestals = yes
doPixelStatus = no
doCommonMode = yes

[cspad_mod.CsPadCalib:1]
source      = DetInfo(XppGon.0:Cspad2x2.1)
inputKey    =
outputKey   = calibrated_arr1
doPedestals = yes

```



```

doPixelStatus = no
doCommonMode = yes

[CSPadPixCoords.CSPad2x2ImageProducer:0]
calibDir      = /reg/d/psdm/xpp/xpptut13/xtc/calib
typeGroupName = CsPad2x2::CalibV1
source        = DetInfo(XppGon.0:Cspad2x2.0)
inkey         = calibrated_arr0
outimgkey     = Image
tiltIsApplied = false
useWidePixCenter = false
print_bits    = 15

[CSPadPixCoords.CSPad2x2ImageProducer:1]
calibDir      = /reg/d/psdm/xpp/xpptut13/xtc/calib
typeGroupName = CsPad2x2::CalibV1
source        = DetInfo(XppGon.0:Cspad2x2.1)
inkey         = calibrated_arr1
outimgkey     = Image
tiltIsApplied = false
useWidePixCenter = false
print_bits    = 15

[ImgAlgos.ImgSaveInFile:0]
source        = DetInfo(:Cspad2x2.0)
key           = Image
fname         = cspad2x2.0
ftype         = txt
#ftype        = tiff
saveAll       = true
print_bits    = 3
#eventSave    = 5

[ImgAlgos.ImgSaveInFile:1]
source        = DetInfo(:Cspad2x2.1)
key           = Image
fname         = cspad2x2.1
#ftype        = txt
ftype         = tiff
saveAll       = true
print_bits    = 3
#eventSave    = 5

```

Example for package ImgPixSpectra

See [Package ImgPixSpectra](#)

Modules:

- `ImgPixSpectra::CSPadPixSpectra`
- `ImgPixSpectra::CSPad2x2PixSpectra`
- `ImgPixSpectra::CameraPixSpectra`

Example for module `ImgPixSpectra::CSPadPixSpectra`

See module description in [Module `ImgPixSpectra::CSPadPixSpectra`](#)

Configuration file `psana-cxib2313-r0114-cspad-pix-spectra.cfg`:

```
# Command to run this script:
# psana -c psana-cxib2313-r0114-cspad-pix-spectra.cfg
#
# Other useful commands:
# psana -n 5 -m PrintSeparator,PrintEventId,psana_examples.DumpCsPad /reg/d/psdm/cxi/cxib2313/xtc/e336-r0114*
# psana -n 5 -m EventKeys exp=cxib2313:run=114:xtc

[psana]
files      = exp=cxib2313:run=114:xtc
#calib-dir = ./calib
skip-events = 0
events     = 100
modules    = cspad_mod.CsPadCalib ImgPixSpectra.CSPadPixSpectra

[cspad_mod.CsPadCalib]
source     = DetInfo(CxiDs1.0:Cspad.0)
inputKey   =
outputKey  = calibrated
doPedestals = yes
doPixelStatus = no
doCommonMode = yes

[ImgPixSpectra.CSPadPixSpectra]
source     = CxiDs1.0:Cspad.0
inputKey   = calibrated
amin       = -20.
amax       = 20.
nbins      = 10
arr_fname  = cspad_spectral_array.txt
```

where module `cspad_mod.CsPadCalib` subtracts pedestals, apply common mode correction, and save CSPAD array in the event store with key "calibrated", which is used in the next module `ImgPixSpectra.CSPadPixSpectra`.

To run this script use command

```
psana -c psana-cxib2313-r0114-cspad-pix-spectra.cfg
```

which produces two files:

- `cspad_spectral_array.txt` – array of 10-bin amplitude spectra for all pixels
- `cspad_spectral_array.txt.sha`– file with metadata for array shape:

```
NPIXELS  2296960
NBINS    10
AMIN     -20
AMAX     20
NEVENTS  100
ARRFNAME cspad_spectral_array.txt
```

Plot for content of the file `cspad_spectral_array.txt`:

.

Example for module `ImgPixSpectra::CSPad2x2PixSpectra`

See module description in [Module `ImgPixSpectra::CSPad2x2PixSpectra`](#)
 Configuration file example for `CSPad2x2PixSpectra`:

```
[psana]
files      = /reg/d/psdm/<instrument>/<experiment>/<file-name>.xtc
modules    = ImgPixSpectra.CSPad2x2PixSpectra
```

```
[ImgPixSpectra.CSPad2x2PixSpectra]
source     = CxiSc1.0:Cspad2x2.0
amin       = 500.
amax       = 1000.
nbins      = 100
arr_fname  = cspad2x2-pix-spectra.txt
```

To get images from saved file one may execute the auxiliary script:

```
ImgPixSpectra/data/PlotSpectralArrayFromFile.py cspad2x2-pix-spectra.txt
```

generates image for limited range of pixels for CSPad, CSPad2x2, or Camera, respectively:

```
--
```

Examples for package ImgAlgos

See [Package ImgAlgos](#)

Example for module ImgAlgos::Tahometer

See [Module ImgAlgos::Tahometer](#)

Example of the psana configuration file:

```
[psana]
files      = /reg/d/psdm/<INSTRUMENT>/<experiment>/xtc/<file-name-1>.xtc\
            /reg/d/psdm/<INSTRUMENT>/<experiment>/xtc/<file-name-2>.xtc

modules    = ImgAlgos.Tahometer

[ImgAlgos.Tahometer]
dn         = 10
print_bits = 7
```

Example for module ImgAlgos::EpixNDArrProducer

See [Module ImgAlgos::EpixNDArrProducer](#)

Example of the psana configuration file:

```
[psana]
files      = exp=xppi0414:run=94
events     = 100
modules    = ImgAlgos.EpixNDArrProducer \
            ImgAlgos.NDArrAverage

[ImgAlgos.EpixNDArrProducer]
source     = DetInfo(:Epix10k)
key_in     =
key_out    = epix-nda
outtype    = float
print_bits = 255

[ImgAlgos.NDArrAverage]
source     = DetInfo(:Epix10k)
key        = epix-nda
avefile    = epix-ave
print_bits = 29
```

EpixNDArrProducer gets Epix10k ndarray from data and saves it in the events store with possible type conversion. Then NDArrAverage averages this ndarray over requested number of events (100) and saves it in file.

Example for module ImgAlgos::PnccdImageProducer

See [Module ImgAlgos::PnccdImageProducer](#)

Example of the psana configuration file:

```
[psana]
files      = /reg/d/psdm/<INSTRUMENT>/<experiment>/xte/<file-name-1>.xte\
             /reg/d/psdm/<INSTRUMENT>/<experiment>/xte/<file-name-2>.xte
#skip-events = 100
events      = 5
modules = ImgAlgos.PnccdImageProducer ImgAlgos.ImgSaveInFile

[ImgAlgos.PnccdImageProducer]
source      = DetInfo(:pnCCD)
inkey       =
outimgkey   = imgpnccd
print_bits  = 1

[ImgAlgos.ImgSaveInFile]
source      = DetInfo(:pnCCD)
key         = imgpnccd
fname       = pnccd-img-ev
saveAll     = true
#eventSave  = 82
print_bits  = 1
```

This script saves text files with images like pnccd-img-ev-<run-date-time.nsec>.txt, which can be presented as:

.

Advanced example for [PnccdImageProducer](#):

- get pnccd ndarray from data,
- calibrate ndarray (subtract pedestals, common mode, remove pixels with bad status),
- produce image with two gaps from calibrated ndarray,
- average image for 10 events:

```
[psana]
#calib-dir = /reg/d/psdm/SXR/sxrb5914/calib
files = exp=sxrb5914:run=245
events = 10
```

```
modules = ImgAlgos.Tahometer \
          ImgAlgos.PncddNDArrProducer \
          ImgAlgos.NDArrCalib \
          ImgAlgos.PncddImageProducer \
          ImgAlgos.NDArrAverage
```

```
[ImgAlgos.Tahometer]
dn      = 100
print_bits = 7
```

```
[ImgAlgos.PncddNDArrProducer]
source  = DetInfo(Camp.0:pnCCD.1)
key_in  =
key_out = pncdd-ndarr
outtype = asdata
print_bits = 0
```

```
[ImgAlgos.NDArrCalib]
source = DetInfo(Camp.0:pnCCD.1)
key_in = pncdd-ndarr
key_out = calibrated
do_peds = yes
do_cmod = yes
do_stat = yes
do_mask = no
do_bkgd = no
do_gain = no
do_nrms = no
do_thre = no
#fname_mask = pncdd-test-mask.txt
#fname_bkgd = pncdd-test-bkgd.txt
masked_value = 0
threshold_nrms = 4.0
threshold = 100
below_thre_value = 0
bkgd_ind_min = 10000
bkgd_ind_max = 10200
bkgd_ind_inc = 1
print_bits = 1
```

```
[ImgAlgos.PncddImageProducer]
source = DetInfo(Camp.0:pnCCD.1)
inkey = calibrated
outimgkey = pncdd-img
gap_rows = 0
gap_cols = 16
gap_value = 0
print_bits = 1
```

```
[ImgAlgos.NDArrAverage]
source = DetInfo(Camp.0:pnCCD.1)
key = pncdd-img
avefile = pncdd-ave
rmsfile = pncdd-rms
#maskfile = pncdd-msk
#hotpixfile = pncdd-hot
thr_rms_ADU = 160
thr_min_ADU = 2
thr_max_ADU = 10000
print_bits = 29
```

Example for module `ImgAlgos::CameraImageProducer`

See [Module `ImgAlgos::CameraImageProducer`](#)

Example of the psana configuration file:

```
[psana]
files      = /reg/d/psdm/<INSTRUMENT>/<experiment>/xtc/<file-name-1>.xtc
modules    = ImgAlgos.CameraImageProducer  ImgAlgos.ImgSaveInFile
events     = 5

[ImgAlgos.CameraImageProducer]
source     = DetInfo(:Opal1000)
key_in     =
key_out    = img
subtract_offset = true
print_bits = 15

[ImgAlgos.ImgSaveInFile]
source     = DetInfo(:Opal1000)
key        = img
fname      = img-from-my-experiment
saveAll    = true
#eventSave = 1
```

This script saves text files with images like `img-from-my-experiment-<run-date-time.nsec>.txt`, which can be drawn by the python script

```
./ImgAlgos/data/PlotCameraImageFromFile.py <filename>.txt <Amin> <Amax>
```

Example for module `ImgAlgos::PrincetonImageProducer`

See [Module `ImgAlgos::PrincetonImageProducer`](#)

Example of the psana configuration file:

```
[psana]
files      = /reg/d/psdm/<INSTRUMENT>/<experiment>/xtc/<file-name-1>.xtc\
            /reg/d/psdm/<INSTRUMENT>/<experiment>/xtc/<file-name-2>.xtc
modules    = ImgAlgos.PrincetonImageProducer \
            ImgAlgos.ImgSaveInFile
events     = 3

[ImgAlgos.PrincetonImageProducer]
source     = DetInfo(:Princeton)
key_in     =
key_out    = img
subtract_offset = true
print_bits = 31

[ImgAlgos.ImgSaveInFile]
source     = DetInfo(:Princeton)
key        = img
fname      = img-xcs
saveAll    = true
print_bits = 31
```

Example for module `ImgAlgos::AcqirisArrProducer`

See description of parameters in [Module `ImgAlgos::AcqirisArrProducer`](#)

Example of the psana configuration file `psana-am001509-r0125-acqiris.cfg`:

```

# Command to run this script:
# psana -c psana-amo01509-r0125-acqiris.cfg;
#
# Useful commands:
# psana -n 5 -m EventKeys exp=amo01509:run=125:xtc > test-acqiris-file.txt;
# psana -n 5 -m psana_examples.DumpAcqiris exp=amo01509:run=125:xtc > test-acqiris-file.txt;    (C++ version)
# psana -n 1 -m psana_examples.dump_acqiris exp=amo01509:run=125:xtc > test-acqiris-file.txt;    (python version)

[psana]
#files = /reg/d/psdm/AM0/amo01509/xtc/e8-r0125-s00-c00.xtc /reg/d/psdm/AM0/amo01509/xtc/e8-r0125-s01-c00.xtc
files = exp=amo01509:run=125:xtc

modules = ImgAlgos.AcqirisArrProducer ImgAlgos.ImgSaveInFile:wf ImgAlgos.ImgSaveInFile:wt

skip-events = 0
events      = 100

[ImgAlgos.AcqirisArrProducer]
source      = AmoETOFF.0:Acqiris.0
key_in      =
key_wform   = acqiris_wform
key_wtime   = acqiris_wtime
fname_prefix = acq
print_bits  = 11

[ImgAlgos.ImgSaveInFile:wf]
source      = AmoETOFF.0:Acqiris.0
key         = acqiris_wform
fname       = acq-AmoETOFF-wform
ftype       = txt
#ftype      = tiff
#saveAll    = true
print_bits  = 3
eventSave   = 5

[ImgAlgos.ImgSaveInFile:wt]
source      = AmoETOFF.0:Acqiris.0
key         = acqiris_wtime
fname       = acq-AmoETOFF-wtime
ftype       = txt
#ftype      = tiff
#saveAll    = true
print_bits  = 3
eventSave   = 5

```

This script with psana does a few things:

- module AcqirisArrProducer gets Acqiris data objects from event store, adjusts trigger time corrections, and saves them back in the event store as uniform ndarrays<double,2> objects for waveforms and times
- two instances of the module ImgSaveInFile save arrays of waveforms and wave-times for locally counted event 5.

This script saves 3 text files:

- acq-amo01509-r0125.txt -- with configuration parameters:

```

Acqiris::ConfigV1:
  nbrBanks=1 channelMask=69905 nbrChannels=5 nbrConvertersPerChannel=4
  horiz: sampInterval=2.5e-10 delayTime=0 nbrSegments=1 nbrSamples=10000
  vert(0): fullScale=0.1 slope=1.52588e-06 offset=0 coupling=3 bandwidth=0
  vert(1): fullScale=0.1 slope=1.52588e-06 offset=0 coupling=3 bandwidth=0
  vert(2): fullScale=0.1 slope=1.52588e-06 offset=0 coupling=3 bandwidth=0
  vert(3): fullScale=0.1 slope=1.52588e-06 offset=0 coupling=3 bandwidth=0
  vert(4): fullScale=0.1 slope=1.52588e-06 offset=0 coupling=3 bandwidth=0

```

- acq-AmoETOFF-wform-r0125-20091018-182512.194787218.txt – with waveforms for 5th event

- acq-AmoETOF-wtime-r0125-20091018-182512.194787218.txt – with wave-times for 5th event

Arrays from these files can be plotted for all channels, shown by different colors:

Example for module `ImgAlgos::AcqirisAverage`

See description of parameters in [Module `ImgAlgos::AcqirisAverage`](#)

Configuration file `psana-amo01509-r0125-acqiris-average.cfg`:

```
# Command to run this script:
# psana -c psana-amo01509-r0125-acqiris-average.cfg;
#
# Useful commands:
# psana -n 5 -m EventKeys exp=amo01509:run=125:xtc > test-acqiris-file.txt;
# psana -n 5 -m psana_examples.DumpAcqiris exp=amo01509:run=125:xtc > test-acqiris-file.txt;    (C++ version)
# psana -n 1 -m psana_examples.dump_acqiris exp=amo01509:run=125:xtc > test-acqiris-file.txt;    (python version)

[psana]
#files = /reg/d/psdm/AMO/amo01509/xtc/e8-r0125-s00-c00.xtc /reg/d/psdm/AMO/amo01509/xtc/e8-r0125-s01-c00.xtc
files = exp=amo01509:run=125:xtc

modules = ImgAlgos.AcqirisArrProducer ImgAlgos.AcqirisAverage

skip-events = 0
events      = 1000

[ImgAlgos.AcqirisArrProducer]
source      = AmoETOF.0:Acqiris.0
key_in      =
key_wform   = acqiris_wform
key_wtime   = acqiris_wtime
fname_prefix = acq
print_bits  = 3

[ImgAlgos.AcqirisAverage]
source      = AmoETOF.0:Acqiris.0
key_in      = acqiris_wform
key_ave     = acqiris_average
fname_ave_prefix = acq
thresholds  = -0.005 -0.005 -0.005 -0.005 -0.005
is_positive_signal = no
do_inverse_selection = yes
#skip_events = 0
#proc_events = 100
print_bits  = 255
```

Psana with this script runs over 1000 events apply threshold-based selection algorithm and produces files:

- `acq-amo01509-r0125-config.txt` -- with Acqiris configuration parameters:

```
Acqiris::ConfigV1:
  nbrBanks=1 channelMask=69905 nbrChannels=5 nbrConvertersPerChannel=4
  horiz: sampInterval=2.5e-10 delayTime=0 nbrSegments=1 nbrSamples=10000
  vert(0): fullScale=0.1 slope=1.52588e-06 offset=0 coupling=3 bandwidth=0
  vert(1): fullScale=0.1 slope=1.52588e-06 offset=0 coupling=3 bandwidth=0
  vert(2): fullScale=0.1 slope=1.52588e-06 offset=0 coupling=3 bandwidth=0
  vert(3): fullScale=0.1 slope=1.52588e-06 offset=0 coupling=3 bandwidth=0
  vert(4): fullScale=0.1 slope=1.52588e-06 offset=0 coupling=3 bandwidth=0
```


- acq-amo01509-r0125-ave-wfs.txt with averaged array of waveforms, which were below threshold (averaging for baseline level):

Example for module `ImgAlgos::AcqirisCalib`

See description of parameters in [Module `ImgAlgos::AcqirisCalib`](#)

Configuration file `psana-am01509-r0125-acqiris-calib.cfg`:

```
# Command to run this script:
# psana -c psana-am01509-r0125-acqiris-calib.cfg;
#
# Useful commands:
# psana -n 5 -m EventKeys exp=amo01509:run=125:xtc > test-acqiris-file.txt;
# psana -n 5 -m psana_examples.DumpAcqiris exp=amo01509:run=125:xtc > test-acqiris-file.txt;    (C++ version)
# psana -n 1 -m psana_examples.dump_acqiris exp=amo01509:run=125:xtc > test-acqiris-file.txt;    (python version)

[psana]
#files = /reg/d/psdm/AM0/amo01509/xtc/e8-r0125-s00-c00.xtc /reg/d/psdm/AM0/amo01509/xtc/e8-r0125-s01-c00.xtc
files = exp=amo01509:run=125:xtc

modules = ImgAlgos.AcqirisArrProducer ImgAlgos.AcqirisCalib ImgAlgos.ImgSaveInFile:wf_raw ImgAlgos.
ImgSaveInFile:wf_calib

skip-events = 0
events      = 10

[ImgAlgos.AcqirisArrProducer]
source      = AmoETOFF.0:Acqiris.0
key_in      =
key_wform   = acqiris_wform
key_wtime   = acqiris_wtime
fname_prefix = acq
print_bits  = 1

[ImgAlgos.AcqirisCalib]
source      = AmoETOFF.0:Acqiris.0
key_in      = acqiris_wform
key_out     = wf-calibrated
fname_base_line = acq-am01509-r0125-ave-wfs.txt
#skip_events = 0
#proc_events = 100
print_bits  = 255

[ImgAlgos.ImgSaveInFile:wf_raw]
source      = AmoETOFF.0:Acqiris.0
key         = acqiris_wform
fname       = acq-AmoETOFF-wform-raw
ftype       = txt
#ftype      = tiff
#saveAll    = true
print_bits  = 3
eventSave   = 5

[ImgAlgos.ImgSaveInFile:wf_calib]
source      = AmoETOFF.0:Acqiris.0
key         = wf-calibrated
fname       = acq-AmoETOFF-wform-calibrated
ftype       = txt
#ftype      = tiff
#saveAll    = true
print_bits  = 3
eventSave   = 5
```

In this script the base-line level for all waveforms is loaded from file and is subtracted in module `ImgAlgos.AcqirisCalib`.

For example, the raw and calibrated waveforms for event 5 were saved in the files:

- `acq-AmoETOF-wform-raw-r0125-e000000005-20091018-182512.194787218.txt`
- `acq-AmoETOF-wform-calibrated-r0125-e000000005-20091018-182512.194787218.txt`

which content is presented on plots:

Example for module `ImgAlgos::AcqirisCFD`

See description of parameters in [Module `ImgAlgos::AcqirisCFD`](#)

Configuration file `psana_cfd.cfg`:

```
# Command to run this script:
# psana -c psana_cfd.cfg;
#
[psana]
modules = ImgAlgos.AcqirisArrProducer ImgAlgos.AcqirisCalib ImgAlgos.AcqirisCFD
files = /reg/d/psdm/AMO/amo01509/xtc/e8-r0125-s00-c00.xtc
[ImgAlgos.AcqirisArrProducer]
source = AmoETOF.0:Acqiris.0
key_in =
key_wform = acqiris_wform
key_wtime = acqiris_wtime
fname_prefix = acq
print_bits = 0
[ImgAlgos.AcqirisCalib]
source = AmoETOF.0:Acqiris.0
key_in = acqiris_wform
key_out = wf-calibrated
fname_base_line = acq--r0000-ave-wfs.txt
#skip_events = 0
#proc_events = 100
print_bits = 0
[ImgAlgos.AcqirisCFD]
source = AmoETOF.0:Acqiris.0
key_wform = wf-calibrated
baselines = 0.0 0.0 0.0 0.0 0.0
thresholds = -0.005 -0.005 -0.005 -0.005 -0.005
fractions = 0.5 0.5 0.5 0.5 0.5
deadtimes = 20.0e-9 20.0e-9 20.0e-9 20.0e-9 20.0e-9
leading_edges = 1 1 1 1 1
```

This script analyzes an AMO run where 5 acqiris channels were in use. It uses an `AcqirisArrProducer` to producer the list of waveforms/times for all channels. It then uses an `AcqirisCalib` module to do a baseline subtraction using a previously generated file of baselines `act--r0000-ave-wfs.txt`. Finally, a constant fraction discriminator algorithm is run on all the waveforms with user specified parameters. The edges calculated by `AcqirisCFD` are saved to the event as one `ndarray<double,1>` per channel, each with a (default) key "`acqiris_edges_N`" where N is the channel number. Channels where no edges were found are not saved to the event.

Example for combination of Acqiris modules

See description of parameters in [Modules `ImgAlgos::AcqirisArrProducer`, `AcqirisAverage`, and `AcqirisCalib`](#)

Configuration file `psana-amo01509-r0125-acqiris-comb.cfg`:

```

# Command to run this script:
# psana -c psana-amo01509-r0125-acqiris-comb.cfg;

[psana]
files = exp=amo01509:run=125:xtc

modules = ImgAlgos.AcqirisArrProducer ImgAlgos.AcqirisAverage:bl ImgAlgos.AcqirisCalib ImgAlgos.AcqirisAverage:
signal ImgAlgos.Tahometer

#skip-events = 0
events      = 2010

[ImgAlgos.AcqirisArrProducer]
source      = AmoETOF.0:Acqiris.0
key_in      =
key_wform   = acqiris_wform
key_wtime   = acqiris_wtime
fname_prefix = acq
print_bits  = 7

[ImgAlgos.AcqirisAverage:bl]
source      = AmoETOF.0:Acqiris.0
key_in      = acqiris_wform
#key_ave    =
fname_ave_prefix = acq-bline
thresholds  = -0.005 -0.005 -0.005 -0.005 -0.005
is_positive_signal = no
do_inverse_selection = yes
skip_events = 0
proc_events = 1000
print_bits  = 31

[ImgAlgos.AcqirisCalib]
source      = AmoETOF.0:Acqiris.0
key_in      = acqiris_wform
key_out     = wf-calibrated
fname_base_line = acq-bline
skip_events = 1001
proc_events = 1000
print_bits  = 47

[ImgAlgos.AcqirisAverage:sigal]
source      = AmoETOF.0:Acqiris.0
key_in      = wf-calibrated
#key_ave    =
fname_ave_prefix = acq-signal
thresholds  = -0.01 -0.01 -0.01 -0.01 -0.01
is_positive_signal = no
do_inverse_selection = no
skip_events = 1001
proc_events = 1000
print_bits  = 31

[ImgAlgos.Tahometer]
print_bits  = 7

```

This script works with psana as follows:

- for the 1st 1000 events averages waveforms below threshold and saves results in the file acq-bline-amo01509-r0125-ave-wfs.txt;
- for the next 1000 events subtracts baseline level and averages waveforms above thresholds and saves results in the file acq-signal-amo01509-r0125-ave-wfs.txt.

This script produces three files:

- acq-amo01509-r0125-config.txt - with Acqris configuration parameters:

```
Acqiris::ConfigV1:
  nbrBanks=1 channelMask=69905 nbrChannels=5 nbrConvertersPerChannel=4
  horiz: sampInterval=2.5e-10 delayTime=0 nbrSegments=1 nbrSamples=10000
  vert(0): fullScale=0.1 slope=1.52588e-06 offset=0 coupling=3 bandwidth=0
  vert(1): fullScale=0.1 slope=1.52588e-06 offset=0 coupling=3 bandwidth=0
  vert(2): fullScale=0.1 slope=1.52588e-06 offset=0 coupling=3 bandwidth=0
  vert(3): fullScale=0.1 slope=1.52588e-06 offset=0 coupling=3 bandwidth=0
  vert(4): fullScale=0.1 slope=1.52588e-06 offset=0 coupling=3 bandwidth=0
```

- acq-blinc-amo01509-r0125-ave-wfs.txt - contains baseline averaged arrays, which can be presented by full scale and zoomed plots:
- acq-signal-amo01509-r0125-ave-wfs.txt- contains signal averaged arrays, which can be presented by full scale and zoomed plots:

Example for module `ImgAlgos::NDArrImageProducer`

See [Module `ImgAlgos::NDArrImageProducer`](#)

Module `ImgAlgos.NDArrImageProducer` produces image from any detector-associated ndarray

Possible chain of modules:

- `<Package>.<Detector>NDArrProducer`
- `ImgAlgos.NDArrCalib`
- `ImgAlgos.NDArrImageProducer`
- `ImgAlgos.NDArrAverage`

Example of the configuration file for cspad

```

[psana]
files = exp=cxi0114:run=227
events = 100

modules = ImgAlgos.Tahometer \
          CSPadPixCoords.CSPadNDArrProducer \
          ImgAlgos.NDArrCalib \
          ImgAlgos.NDArrImageProducer \
          ImgAlgos.NDArrAverage

[CSPadPixCoords.CSPadNDArrProducer]
source      = DetInfo(CxiDs1.0:Cspad.0)
inkey       =
outkey      = cspad_ndarr
outtype     = float
is_fullsize = yes
print_bits  = 3

[ImgAlgos.NDArrCalib]
source      = DetInfo(CxiDs1.0:Cspad.0)
key_in     = cspad_ndarr
key_out    = calibrated
do_peds    = yes
do_cmod    = yes
do_stat    = yes
do_mask    = no
do_bkgd    = no
do_gain    = no
do_nrms    = no
do_thre    = no
fname_mask =
fname_bkgd =
masked_value = -10
threshold_nrms = 4
threshold   = 7
below_thre_value = 0
print_bits  = 3

[ImgAlgos.NDArrImageProducer]
source      = CxiDs1.0:Cspad.0
key_in     = calibrated
key_out    = cspad_img
#type_out  = asinp
#type_out  = float
#x0_off_pix = 50
#y0_off_pix = 50
print_bits = 255
#calibdir  = /reg/neh/home/dubrovin/LCLS/CSPadAlignment-v01/calib-cxi-ds1-2014-03-19/calib

[ImgAlgos.NDArrAverage]
source      = DetInfo(CxiDs1.0:Cspad.0)
key         = cspad_img
avefile     = cspad-img-ave
rmsfile     = cspad-img-rms
#maskfile   = cspad-img-msk
#hotpixfile = cspad-img-hot
thr_rms_ADU = 0
thr_min_ADU = 4
thr_max_ADU = 65000
print_bits  = 29

```

Example of the configuration file for cspad2x2

```

[psana]
# calib-dir = /reg/d/psdm/mec/meca1113/calib
#calib-dir = /reg/neh/home/dubrovin/LCLS/CSPad2x2Alignment/calib-cspad2x2-01-2013-02-13/calib

files = exp=meca1113:run=376
events = 10
#skip-events = 0

modules = CSPadPixCoords.CSPad2x2NDArrProducer \
          ImgAlgos.NDArrImageProducer \
          ImgAlgos.NDArrAverage

[CSPadPixCoords.CSPad2x2NDArrProducer]
source      = MecTargetChamber.0:Cspad2x2.1
inkey       =
outkey      = cspad2x2_ndarr
outtype     = int16
print_bits  = 5

[ImgAlgos.NDArrImageProducer]
source      = MecTargetChamber.0:Cspad2x2.1
key_in      = cspad2x2_ndarr
key_out     = cspad2x2_img
#type_out   = asinp
#type_out   = float
#x0_off_pix = 50
#y0_off_pix = 50
print_bits  = 255
#oname      = CSPAD2X2:V1
#oindex     = 0
#pix_scale_size_um = 218.
#calibdir   = /reg/neh/home/dubrovin/LCLS/CSPad2x2Alignment/calib-cspad2x2-01-2013-02-13/calib
#calibgroup = CsPad2x2::CalibV1

[ImgAlgos.NDArrAverage]
source      = MecTargetChamber.0:Cspad2x2.1
#key        = cspad2x2_ndarr
key         = cspad2x2_img
avefile     = arr-ave
rmsfile     = arr-rms
maskfile    = arr-msk
hotpixfile  = arr-hot
thr_rms_ADU = 10
#thr_min_ADU = 2
#thr_max_ADU = 20000
print_bits  = 255

```

Example for module `ImgAlgos::NDArrAverage`

- See [Module `ImgAlgos::NDArrAverage`](#) and [Module `CSPadPixCoords::CSPadNDArrProducer`](#)
- The `NDArrAverage` module in combination with `CSPadNDArrProducer` (or any other device `NDArrProducer`) can be used for evaluation of averaged pedestals or background using dedicated runs. Typical configuration file may looks like this:

```
# Run this script:
# psana -c psana-meca1113-r0045-cspad-cspad-dark-hotpix.cfg

[psana]
files = exp=meca1113:run=45
events = 400
#skip-events = 0

modules = CSPadPixCoords.CSPadNDArrProducer \
          ImgAlgos.NDArrAverage \
          ImgAlgos.Tahometer

[CSPadPixCoords.CSPadNDArrProducer]
source      = MecTargetChamber.0:Cspad.0
inkey       =
outkey      = cspad_ndarr
outtype     = int16
is_fullsize = yes
print_bits  = 3

[ImgAlgos.NDArrAverage]
source      = MecTargetChamber.0:Cspad.0
key         = cspad_ndarr
avefile     = cspad.0-ave
rmsfile     = cspad.0-rms
maskfile    = cspad.0-msk
hotpixfile  = cspad.0-hot
thr_rms_ADU = 10
thr_min_ADU = 4
thr_max_ADU = 10000
print_bits  = 29
#evts_stage1 = 100
#gate_width1 = 500.
#evts_stage2 = 200
#gate_width2 = 100.

[ImgAlgos.Tahometer]
dn          = 100
print_bits  = 7
```

- Module `ImgAlgos.Tahometer` is not required in this configuration file and is added for convenience to print timing statistics for this job.
- Evaluation of average intensity in 2 or 3 stages using gate-based algorithms excludes out-layers in intensity spectra and makes average more stable and reliable. However, the gate width is not an universal parameter. In order to set this parameter one has to look at spectrum of averaged intensities for particular device. The same is valid for evaluation of hot/bad pixel masks. One has to plot spectra of averaged intensity and rms values. Averaged intensity and rms spectra for `exp=meca1113:run=45` are shown on plots, respectively:

Then, one has to decide how to set parameters for `NDArrAverage` algorithms, for example, it is quite safe to use

- `thr_rms_ADU = 10` – to discard very noisy pixels,
- `thr_min_ADU = 4` – to discard presumably dead pixels with 0-intensity,
- `thr_max_ADU = 10000` – to discard pixels with intensity significantly exceeding average value. To be on safe side for `int16` data this parameter can be set to $2^{16}-4$, where 4 in both cases is just a small arbitrary number for spare safety gap.

Example for module `ImgAlgos::NDArrCalib`

- See [Module `ImgAlgos::NDArrCalib`](#), [Module `ImgAlgos::NDArrAverage`](#) The `NDArrCalib` module in combination with any device `NDArrProducer` (for example `PncddNDArrProducer`) can be used to apply intensity corrections to ndarray. Typical configuration files are shown below.

Example of `ImgAlgos::NDArrCalib` for `pnCCD`

```
# Run this script:
# psana -c psana-meca1113-r0045-cspad-cspad-dark-hotpix.cfg

[psana]
files = exp=meca1113:run=45
events = 400
#skip-events = 0

modules = CSPadPixCoords.CSPadNDArrProducer \
          ImgAlgos.NDArrAverage \
          ImgAlgos.Tahometer \
          ImgAlgos.NDArrCalib

[CSPadPixCoords.CSPadNDArrProducer]
source      = MecTargetChamber.0:Cspad.0
inkey       =
outkey      = cspad_ndarr
outtype     = int16
is_fullsize = yes
print_bits  = 3

[ImgAlgos.NDArrAverage]
source      = MecTargetChamber.0:Cspad.0
key         = cspad_ndarr
avefile     = cspad.0-ave
rmsfile     = cspad.0-rms
maskfile    = cspad.0-msk
hotpixfile  = cspad.0-hot
thr_rms_ADU = 10
thr_min_ADU = 4
thr_max_ADU = 10000
print_bits  = 29
#evts_stage1 = 100
#gate_width1 = 500.
#evts_stage2 = 200
#gate_width2 = 100.

[ImgAlgos.Tahometer]
dn          = 100
print_bits  = 7

[ImgAlgos.NDArrCalib]
source      = MecTargetChamber.0:Cspad.0
key         = cspad_ndarr
avefile     = cspad.0-ave
rmsfile     = cspad.0-rms
maskfile    = cspad.0-msk
hotpixfile  = cspad.0-hot
thr_rms_ADU = 10
thr_min_ADU = 4
thr_max_ADU = 10000
print_bits  = 29
#evts_stage1 = 100
#gate_width1 = 500.
#evts_stage2 = 200
#gate_width2 = 100.
```

```

[psana]
files = exp=amoa1214:run=7
#skip-events = 100
events      = 5

modules = ImgAlgos.Tahometer \
          ImgAlgos.PnccdNDArrProducer \
          ImgAlgos.NDArrCalib \
          ImgAlgos.PnccdImageProducer \
          ImgAlgos.ImgSaveInFile
#          ImgAlgos.NDArrAverage \

[ImgAlgos.PnccdNDArrProducer]
source      = DetInfo(Camp.0:pnCCD.0)
key_in      =
key_out     = pnccd-ndarr
outtype     = asdata
print_bits  = 13

[ImgAlgos.NDArrCalib]
source      = DetInfo(Camp.0:pnCCD.0)
key_in      = pnccd-ndarr
key_out     = calibrated
outtype     = float
do_peds     = yes
do_cmod     = yes
do_stat     = no
do_mask     = no
do_bkgd     = no
do_gain     = no
do_nrms     = no
do_thre     = no
fname_mask  =
fname_bkgd  =
masked_value      = 0
threshold_nrms    = 3
threshold          = 100
below_thre_value  = 0
bkgd_ind_min      = 0
bkgd_ind_max      = 1000
bkgd_ind_inc      = 10
print_bits        = 255

[ImgAlgos.PnccdImageProducer]
source      = DetInfo(Camp.0:pnCCD.0)
#inkey      = pnccd-ndarr
inkey       = calibrated
outimgkey   = pnccd-img
gap_size    = 20
gap_value   = 0
print_bits  = 1

[ImgAlgos.ImgSaveInFile]
source      = DetInfo(Camp.0:pnCCD.0)
key         = pnccd-img
fname       = pnccd-img-from-arr
ftype       = txt
saveAll     = true
print_bits  = 31

[ImgAlgos.Tahometer]
dn          = 100
print_bits  = 7

```


Example of ImgAlgos::NDArrCalib for CSPAD

```
[psana]
files = exp=cxi83714:run=136
events = 100

modules = ImgAlgos.Tahometer \
          CSPadPixCoords.CSPadNDArrProducer \
          ImgAlgos.NDArrCalib \
          ImgAlgos.NDArrAverage

[CSPadPixCoords.CSPadNDArrProducer]
source      = DetInfo(CxiDs1.0:Cspad.0)
inkey       =
outkey      = cspad_ndarr
outtype     = float
is_fullsize = yes
print_bits  = 3

[ImgAlgos.NDArrCalib]
source      = DetInfo(CxiDs1.0:Cspad.0)
key_in     = cspad_ndarr
key_out    = calibrated
outtype    = double
do_peds    = yes
do_cmod    = yes
do_stat    = yes
do_mask    = no
do_bkgd    = no
do_gain    = no
do_nrms    = no
do_thre    = no
fname_mask =
fname_bkgd =
masked_value = -10
threshold_nrms = 4
threshold   = 7
below_thre_value = 0
bkgd_ind_min = 0
bkgd_ind_max = 1000
bkgd_ind_inc = 10
print_bits  = 3

[ImgAlgos.NDArrAverage]
source      = DetInfo(CxiDs1.0:Cspad.0)
key         = calibrated
avefile     = cspad-ave
rmsfile     = cspad-rms
#maskfile   = cspad-msk
#hotpixfile = cspad-hot
thr_rms_ADU = 10
thr_min_ADU = 4
thr_max_ADU = 10000
print_bits  = 29

[ImgAlgos.Tahometer]
dn          = 10
print_bits  = 7
```

For test purpose we use exp=cxi83714:run=136 and loop over 100 exents.

Case 1: raw, non-calibrated images

do_peds = no; 50-60 ms/event

plots for average and rms value distribution for all pixels

Case 2: images with subtracted pedestals and applied status mask

do_peds = yes; 60-70 ms/event
do_stat = yes

Case 3: The same as 2 with common mode subtracted

do_cmod = yes; 107-117 ms/event

Andy's algorithm for CSPAD common mode correction with minor adaptive modifications is applied with parameters:

/reg/d/psdm/cxi/cxi83714/calib/CsPad::CalibV1/CxiDs1.0:Cspad.0/common_mode/135-136.data

1 10 10 100 - algorithm mode, allowed peak mean, allowed peak rms, threshold on number of pixels in ADU bin.

Comparison with Andy's module cspad_mod.CsPadCalib for common mode correction of int16_t data

```
[psana]
# Default calibration directory:
# calib-dir = /reg/d/psdm/mec/cxi83714/calib
# calib-dir = /reg/neh/home1/dubrovin/LCLS/CSPadAlignment-v01/calib-mec-2013-12-10/

files = exp=cxi83714:run=136
events = 100
#skip-events = 4000

modules = cspad_mod.CsPadCalib \
          CSPadPixCoords.CSPadNDArrProducer \
          ImgAlgos.NDArrAverage \
          ImgAlgos.Tahometer

[cspad_mod.CsPadCalib]
inputKey      =
outputKey     = calibrated_data
doPedestals   = yes
doPixelStatus = yes
doCommonMode  = yes

[CSPadPixCoords.CSPadNDArrProducer]
source        = DetInfo(CxiDs1.0:Cspad.0)
inkey         = calibrated_data
outkey        = cspad_ndarr
outtype       = float
is_fullsize   = yes
print_bits    = 3

[ImgAlgos.NDArrAverage]
source        = DetInfo(CxiDs1.0:Cspad.0)
key           = cspad_ndarr
avefile       = cspad-calib-ave
rmsfile       = cspad-calib-rms
#maskfile     = cspad-msk
#hotpixfile    = cspad-hot
thr_rms_ADU   = 10
thr_min_ADU   = 4
thr_max_ADU   = 10000
print_bits    = 29

[ImgAlgos.Tahometer]
dn            = 10
print_bits    = 7
```

doPedestals = yes

doCommonMode = no; 100-110ms

doPedestals = yes
doPixelStatus = yes
doCommonMode = yes; 150-160ms

Conclusion:

- Results of `ImgAlgos.NDArrCalib` well reproduce `cspad_mod.CsPadCalib`
- Common mode correction shrinks the width of averaged intensities from 0.60 to 0.54, and rms from 3.78 to 3.27.
- Algorithm is quite time expensive, it takes 30-40ms/event.

Example of `ImgAlgos::NDArrCalib` for `Fccd960`

- **Calibration**

To calibrate `fccd960` for dark runs use command

```
calibman
```

For test purpose a couple of dark runs were processed and constants were deployed:

```
Content of: /reg/d/psdm/XCS/xcsd7814/calib for detector: Fccd960
Camera::CalibV1
  XcsEndstation.0:Fccd960.0
    pixel_rms
      12-13.data      file is not used
      78-81.data      run range 0078 - 0081
    pedestals
      12-13.data      file is not used
      78-81.data      run range 0078 - 0081
    pixel_status
      12-13.data      file is not used
      78-81.data      run range 0078 - 0081
```

- **Access data in psana**

Example of the configuration file for psana [psana-xcsd7814-r0079-fccd960-aver.cfg](#) can be processed by the command

```
psana -c psana-xcsd7814-r0079-fccd960-aver.cfg
```

Two averaged images are saved in the text files for `exp=xcsd7814:run=79` for raw and calibrated (subtracted dark level) data.

Gain bit coded gain factor 8 is applied automatically.

? Unknown Attachment ? Unknown Attachment

Example for module `ImgAlgos::NDArrDropletFinder`

See description in [psana - Module Catalog](#) and examples in [Peak Finding Module](#)

Example for module `ImgAlgos::PixCoordsProducer`

See [Module `ImgAlgos::PixCoordsProducer`](#)

Psana configuration file which produces enough data to get CSPAD calibrated intensity and coordinate arrays:

```

[psana]
modules      = ImgAlgos.PixCoordsProducer \
               cspad_mod.CsPadCalib \
               CSPadPixCoords.CSPadNDarrProducer

[ImgAlgos.PixCoordsProducer]
source       = CxiDs1.0:Cspad.0
print_bits   = 0

[cspad_mod.CsPadCalib]
source       = CxiDs1.0:Cspad.0
inputKey     =
outputKey    = calibrated
doPedestals  = yes
doPixelStatus = no
doCommonMode = no

[CSPadPixCoords.CSPadNDarrProducer]
source       = CxiDs1.0:Cspad.0
inkey        = calibrated
outkey       = cspad_ndarr
outtype      = int16
is_fullsize  = yes
print_bits   = 3

```

Additional keywords need to be added to the list of module parameters in order to evaluate pixel area and coordinate indexes (for image) arrays. For example:

```

[ImgAlgos.PixCoordsProducer]
source       = CxiDs1.0:Cspad.0
key_out_area = pix_area
key_out_mask = pix_mask
key_out_ix   = pix_ix
key_out_iy   = pix_iy
x0_off_pix   = 1000
y0_off_pix   = 1000
mask_bits    = 15
print_bits   = 255

```

In python code these arrays can be obtained with `env.calibStore().get(...)` method:

```

env = ds.env()
cls = env.calibStore()
src = psana.Source('DetInfo(CxiDs1.0:Cspad.0)')
...
A = evt.get(psana.ndarray_int16_3, src, 'cspad_ndarr').flatten()

X = cls.get(psana.ndarray_float64_1, src, 'x-pix-coords')
Y = cls.get(psana.ndarray_float64_1, src, 'y-pix-coords')
Area = cls.get(psana.ndarray_float64_1, src, 'pix_area')
Mask = cls.get(psana.ndarray_int32_1, src, 'pix_mask')
iX = cls.get(psana.ndarray_uint32_1, src, 'pix_ix')
iY = cls.get(psana.ndarray_uint32_1, src, 'pix_iy')
fname= cls.get(str, src, 'geometry-fname')

// depricated method since ana-0.13.3:
path_nda = cls.get(psana.ndarray_uint8_1, src, 'geometry-calib')
path = ''.join(map(chr, path_nda)) if path_nda is not None else 'N/A'

```

Their shape=(32*185*388,) = (2296960,)

These arrays can be processed by users' code directly. In particular, it is easy to conver them in 2-d image numpy array and plot it:

```

import numpy as np
import matplotlib.pyplot as plt
from PSCalib.GeometryAccess import img_from_pixel_arrays

pix_size = 109.92
xmin, ymin = X.min()-pix_size/2, Y.min()-pix_size/2
iX, iY = np.array((X-xmin)/pix_size, dtype=np.uint), np.array((Y-ymin)/pix_size, dtype=np.uint)
img = img_from_pixel_arrays(iX,iY,W=A)

plt.imshow(img)
plt.draw()

```

Example for module ImgAlgos::ImgAverage

- See [Module ImgAlgos::ImgAverage](#)
- The ImgAverage module can be used for evaluation of averaged pedestals or background using dedicated runs. Typical configuration file may looks like this:

```

[psana]
files      = /reg/d/psdm/<instrument>/<experiment>/xtc/<file-name>.xtc
modules    = ImgAlgos.CameraImageProducer \
              ImgAlgos.ImgAverage
events     = 1000

[ImgAlgos.CameraImageProducer]
source      = DetInfo(:Opal1000)
key_in      =
key_out     = img
subtract_offset = true
print_bits  = 1

[ImgAlgos.ImgAverage]
source      = DetInfo(:Opal1000)
key         = img
avefile     = img-ave.dat
rmsfile     = img-rms.dat
print_bits  = 31
evts_stage1 = 100
evts_stage2 = 100
gate_width1 = 200
gate_width2 = 20

```

Example for module ImgAlgos::ImgMaskEvaluation

- See [Module ImgAlgos::ImgMaskEvaluation](#)
- Configuration parameters for psana:

```

[ImgAlgos.ImgMaskEvaluation]
source      = DetInfo(:Opal1000)
key         = img
file_mask_satu      = img-mask-satu.dat
file_mask_nois      = img-mask-nois.dat
file_mask_comb      = img-mask-comb.dat
file_frac_satu      = img-frac-satu.dat
file_frac_nois      = img-frac-nois.dat
thre_satu           = 400
frac_satu           = 0
dr                  = 1
thre_SoN            = 3
frac_nois           = 0.05
print_bits          = 63

```





In this example parameters were chosen in order to get a small number of "noisy" pixel just due to statistics.

Typical image:

Plots for fraction of noisy pixels:

Plots for fraction of saturated pixels:

Masks: **noisy**, **saturated**, and **combined**:

...

Example for module `ImgAlgos::ImgCalib`

See [Module `ImgAlgos::ImgCalib`](#)

```
[ImgAlgos.CameraImageProducer]
source      = DetInfo(:Opal1000)
key_in      = 
key_out     = img
subtract_offset = true
print_bits  = 1

[ImgAlgos.ImgCalib]
source      = DetInfo(:Opal1000)
key_in      = img
key_out     = calibrated
fname_peds  = <pedestal-file-name>
fname_mask  = <mask-file-name>
fname_bkgd  = 
fname_gain  = 
print_bits  = 31
```

Example of the mask file and resulting image:

-

Example for module `ImgAlgos::ImgRadialCorrection`

See [Module `ImgAlgos::ImgRadialCorrection`](#)

```

[psana]
files = /reg/d/psdm/<instrument>/<experiment>/xtc/<file-name-1>.xtc \
        /reg/d/psdm/<instrument>/<experiment>/xtc/<file-name-2>.xtc \
        ...
        /reg/d/psdm/<instrument>/<experiment>/xtc/<file-name-N>.xtc
skip-events = 500
events      = 10
modules = cspad_mod.CsPadCalib \
#          ImgAlgos.CSPadBkgdSubtract \
          CSPadPixCoords.CSPadImageProducer \
          ImgAlgos.ImgRadialCorrection \
          ImgAlgos.ImgSaveInFile:1

[cspad_mod.CsPadCalib]
inputKey      =
outputKey     = calibrated
doPedestals  = yes
doPixelStatus = no
doCommonMode = no

[ImgAlgos.CSPadBkgdSubtract]
source        = DetInfo(CxiDs1.0:Cspad.0)
inputKey      = calibrated
outputKey     = bkgd_subtracted_arr
bkgd_fname   = <the-file-name-with-background-array>
norm_sector   = 0
print_bits    = 0

[CSPadPixCoords.CSPadImageProducer]
calibDir      = /reg/d/psdm/<instrument>/<experiment>/calib
typeGroupName = CsPad::CalibV1
source        = CxiDs1.0:Cspad.0
key           = calibrated
imgkey        = current_img
#tiltIsApplied = true

[ImgAlgos.ImgRadialCorrection]
source        = DetInfo(CxiDs1.0:Cspad.0)
inkey         = current_img
outkey        = r_cor_img
xcenter       = 866
ycenter       = 857
rmin          = 100
rmax          = 810
n_phi_bins    = 60
event         = 0
print_bits    = 3

[ImgAlgos.ImgSaveInFile:1]
source        = CxiDs1.0:Cspad.0
key           = r_cor_img
fname         = <file-name-for-image-array>
#saveAll      = true
eventSave     = 8

```

Note: the option of the background subtraction (ImgAlgos.CSPadBkgdSubtract) is commented out in this configuration file . In order to evoke this option, the comment sign (#) should be removed from the list of modules and the key=bkgd_subtracted_arr should be used in CSPadPixCoords.CSPadImageProducer.

Calibrated image and spectrum:

-

Calibrated and radial-corrected image, spectrum, and subtracted r-phi65 distribution for n_phi_bins=65:

-

Calibrated and radial-corrected image, spectrum, and subtracted r-phi12 distribution for n_phi_bins=12:

-

Example for module ImgAlgos::ImgPeakFinder

See [Module ImgAlgos::ImgPeakFinder](#)

Configuration file example:

```
[psana]
files      = /reg/d/psdm/<instrument>/<experiment>/xrc/<file-name-1>.xrc \
              /reg/d/psdm/<instrument>/<experiment>/xrc/<file-name-2>.xrc \
              ...
              /reg/d/psdm/<instrument>/<experiment>/xrc/<file-name-N>.xrc

modules     = ImgAlgos.ImgPeakFinder PrintSeparator

#skip-events = 500
events       = 120

[ImgAlgos.ImgPeakFinder]
source       = DetInfo(:Opal1000)
key          =
peaksKey     = peaks
threshold_low = 20
threshold_high = 50
sigma        = 1.5
smear_radius = 2
peak_radius  = 3
xmin         = 200
xmax         = 800
ymin         = 100
ymax         = 900
testEvent    = 115
print_bits   = 0
finderIsOn   = true
```

This algorithm consumes ~15 ms/event on psana0101 for full Opal1000 (1024x1024) camera image.

Smearing algorithm use a "safety margin" which is currently set to 10 pixels (offset from each boarder of the full image size).

Image on different stages of this algorithm:

—

raw image,

image in the window with amplitudes above the threshold_low

- few peaks at the edges were discarded by the window limits,
- image still contains many 1-photon pixels, which need to be eliminated,

smearred image,

raw image with found peaks (marked by the red circles)

zoom of the previous plot.

Example for module ImgAlgos::ImgPeakFilter

```
[ImgAlgos.ImgPeakFilter]
source       = DetInfo(:Opal1000)
key          = peaks
threshold_peak = 5
threshold_total = 0
n_peaks_min  = 10
print_bits   = 11
fname        = img
selection_mode = SELECTION_ON
```

Example for module ImgAlgos::ImgPeakFinderAB

See [Module ImgAlgos::ImgPeakFinderAB](#)

Configuration file example:

```
[psana]
files   = /reg/d/psdm/<instrument>/<experiment>/xtc/<file-name>.xtc
modules = ImgAlgos.CameraImageProducer \
          ImgAlgos.ImgPeakFinderAB
events  = 10

[ImgAlgos.CameraImageProducer]
source      = DetInfo(:Opal1000)
key_in      =
key_out     = img
subtract_offset = true
print_bits  = 1

[ImgAlgos.ImgPeakFinderAB]
source      = DetInfo(:Opal1000)
key         = img
key_peaks_out = peaks
#key_signal_out = signal-arr
#hot_pix_mask_inp_file = ana-misc-exp/mask.dat
#hot_pix_mask_out_file = noise-mask-out.dat
#frac_noisy_evts_file = noise-frac-out.dat
evt_file_out = tmp/img-

rmin      = 10
dr        = 1
SoNThr_noise = 3
SoNThr_signal = 3
frac_noisy_imgs = 0.9
peak_npix_min = 3
peak_npix_max = 100
peak_amp_tot_thr = 0.
peak_SoN_thr = 4.
event_npeak_min = 5
event_npeak_max = 1000
event_amp_tot_thr = 0.
nevents_mask_update = 0
nevents_mask_accum = 50
selection_mode = SELECTION_ON
out_file_bits = 15
print_bits = 513
```

Results:

—

Example for module `ImgAlgos::ImgHitFinder`

See [Module ImgAlgos::ImgHitFinder](#)

`ImgHitFinder` in regular mode needs in file with pedestals (offset) to correct the image and file with threshold.

In amo74213 run 93 these files can be obtained directly from data, discarding signal hits as outliers using

`ImgAlgos.ImgAverage` module as follows with configuration file:

```
# File: psana-am074213-r0093-opal-img-average.cfg

[psana]
#files = /reg/d/psdm/AM0/am074213/xtc/e269-r0093-s05-c00.xtc
files = exp=am074213:run=93:xtc

modules = ImgAlgos.CameraImageProducer \
          ImgAlgos.ImgAverage
skip-events = 0
events      = 1000

[ImgAlgos.CameraImageProducer]
source      = DetInfo(:Opal1000)
key_in      =
key_out     = img
subtract_offset = true
print_bits  = 1

[ImgAlgos.ImgAverage]
source      = DetInfo(:Opal1000)
key         = img
avefile     = img-ave-for-peds
rmsfile     = img-rms-for-thre
evts_stage1 = 100
evts_stage2 = 100
gate_width1 = 50
gate_width2 = 10
print_bits  = 31
```

Run psana using command:

```
psana -c psana-am074213-r0093-opal-img-average.cfg
```

At the end of this procedure two files will be created:

img-ave-for-peds-r0093.dat - may be used for pedestal subtraction:

-

img-rms-for-thre-r0093.dat - may be used multiplied by number of rms as a threshold:

-

The file with accumulated pixel hits can be obtained using configuration file:

```
# File: psana-am074213-r0093-opal-img-hit-finder.cfg

[psana]
#files = /reg/d/psdm/AM0/amo74213/xtc/e269-r0093-s05-c00.xtc
files = exp=amo74213:run=93:xtc

modules = ImgAlgos.CameraImageProducer \
          ImgAlgos.ImgHitFinder \
          ImgAlgos.ImgAverage
skip-events = 0
events      = 1000

[ImgAlgos.CameraImageProducer]
source      = DetInfo(:Opal1000)
key_in      =
key_out     = img
subtract_offset = true
print_bits  = 1

[ImgAlgos.ImgHitFinder]
source      = DetInfo(:Opal1000)
key_in      = img
key_out     = img_hits
fname_peds  = img-ave-for-peds-r0093.dat
fname_mask  =
fname_gain  =
fname_thre  = img-rms-for-thre-r0093.dat
masked_value = 0
thre_mode   = 3
#thre_mode  = 2
thre_param  = 5
thre_below_value = 0
thre_above_value = 1
win_row_min = 10
win_row_max = 1000
win_col_min = 10
win_col_max = 1000
print_bits  = 39

[ImgAlgos.ImgAverage]
source      = DetInfo(:Opal1000)
key         = img_hits
sumfile     = img-sum-result
print_bits  = 25
```

and run it by the command:

```
psana -c psana-am074213-r0093-opal-img-hit-finder.cfg
```

Which creates the file: img-sum-result-r0093.dat for thre_mode = 2 and thre_mode = 3, respectively:

..

Example for module `ImgAlgos::ImgSpectra`

See [Module `ImgAlgos::ImgSpectra`](#)

Configuration file for psana:

```

[psana]
files = /reg/d/psdm/<instrument>/<experiment>/xtc/<file-name>.xtc

modules      =  ImgAlgos.CameraImageProducer \
                 ImgAlgos.ImgSpectra \
                 ImgAlgos.ImgSaveInFile:2
#             ImgAlgos.ImgSaveInFile
#             psana_examples.DumpOpal1k \

events       = 100

[ImgAlgos.CameraImageProducer]
source       = DetInfo(:Opal1000)
key_in       =
key_out      = img
subtract_offset = true
print_bits   = 1

[ImgAlgos.ImgSpectra]
source       = DetInfo(:Opal1000)
key_in       = img
key_out      = spectra
sig_band_rowc = 512.
ref_band_rowc = 552.
sig_band_tilt = 0.
ref_band_tilt = 0.
sig_band_width = 10
ref_band_width = 10
print_bits   = 3

[ImgAlgos.ImgSpectraProc]
source       = DetInfo(:Opal1000)
key_in       = spectra
print_bits   = 15

[ImgAlgos.ImgSaveInFile:2]
source       = DetInfo(:Opal1000)
key          = spectra
fname        = spec-xppi0412
saveAll      = true

[ImgAlgos.ImgSaveInFile]
source       = DetInfo(:Opal1000)
key          = img
fname        = img-xppi0412
saveAll      = true

```

This script can be used in order to produce text files with image and spectral array:

—

or graphics for several images:

—

Example for module `ImgAlgos::ImgSpectraProc`

See [Module `ImgAlgos::ImgSpectraProc`](#)

Configuration file for psana:

```

[ImgAlgos.ImgSpectraProc]
source       = DetInfo(:Opal1000)
key_in       = spectra
print_bits   = 15

```

For each event it prints something similar to:

```
[info:ImgAlgos.ImgSpectraProc] Spectral array shape =3, 1024
[info:ImgAlgos.ImgSpectraProc] Image spectra for run=0060 Evt=000100
Column:      0      100      200      300      400      500      600      700      800      900      1000
Signal:    1211    4062    11150    17070    16406    12949    7991    5168    3968    3542    3811
Refer.:     933    3485    10425    17128    17791    13522    8315    5000    3390    2967    3193
Diff. :    0.259    0.153    0.067   -0.003   -0.081   -0.043   -0.040    0.033    0.157    0.177    0.176
[info:ImgAlgos.ImgSpectraProc] Run=0060 Evt=000100 Time=20120507-125420.982421325 done...
```

Example for module ImgAlgos::ImgSaveInFile

```
modules = ... ImgAlgos.ImgSaveInFile:1 ...

[ImgAlgos.ImgSaveInFile:1]
source      = DetInfo(:0pal1000)  # or CxiDs1.0:Cspad.0
key         = img
fname       = my-img
#ftype      = txt
#ftype      = bin
#ftype      = png
ftype       = tiff
#eventSave  = 5
saveAll     = true
```

See [Module ImgAlgos::ImgSaveInFile](#)

Example for module ImgPeakFinder and ImgPeakFilter for CSPad

Module `ImgAlgos::ImgPeakFinder` works on image. In order to apply this algorithm to CSPad the image should be produced. In next example the image is produced using consecutive modules `cspad_mod.CsPadCalib`, `ImgAlgos.CSPadMaskApply`, and `CSPadPixCoords`.
`CSPadImageProducer`:

```
[psana]
files = /reg/d/psdm/<instrument>/<experiment>/xtc/<file-name-1>.xtc \
...
      /reg/d/psdm/<instrument>/<experiment>/xtc/<file-name-N>.xtc

events      = 1000
experiment  = cxii0212
calib-dir   = ana-cxii0212/calib

modules = cspad_mod.CsPadCalib \
          ImgAlgos.CSPadMaskApply \
          CSPadPixCoords.CSPadImageProducer \
          ImgAlgos.ImgPeakFinder \
          ImgAlgos.ImgPeakFilter \
          ImgAlgos.ImgSaveInFile:1

[cspad_mod.CsPadCalib]
inputKey    =
outputKey   = calibrated
doPedestals = yes
doPixelStatus = no
doCommonMode = yes

[ImgAlgos.CSPadMaskApply]
source      = DetInfo(CxiDs1.0:Cspad.0)
inkey       = calibrated
outkey      = masked_arr
mask_fname  = <your-local-directory>/<mask-file-name>.dat
masked_amp  = 0
print_bits  = 5
mask_control_bits = 15

[CSPadPixCoords.CSPadImageProducer]
calibDir    = /reg/d/psdm/<instrument>/<experiment>/calib
```

```

typeGroupName = Cspad::CalibV1
source        = CxiDs1.0:Cspad.0
key           = masked_arr
imgkey        = img
print_bits    = 0
#tiltIsApplied = true

[ImgAlgos.ImgPeakFinder]
source        = DetInfo(CxiDs1.0:Cspad.0)
key           = img
peaksKey      = peaks
threshold_low = 2
threshold_high = 5
sigma         = 1.5
smear_radius  = 5
peak_radius   = 7
xmin          = 20
xmax          = 1700
ymin          = 20
ymax          = 1700
#testEvent    = 5
print_bits    = 3
#finderIsOn   = true

[ImgAlgos.ImgPeakFilter]
source        = DetInfo(CxiDs1.0:Cspad.0)
key           = peaks
threshold_peak = 5
threshold_total = 0
n_peaks_min   = 10
print_bits    = 11
fname         = cspad-img
selection_mode = SELECTION_ON

[ImgAlgos.ImgSaveInFile:1]
source        = CxiDs1.0:Cspad.0
key           = img
fname         = cspad-img
#eventSave    = 1
saveAll       = true

```

Example for module ImgAlgos::CSPadArrAverage

See [Module ImgAlgos::CSPadArrAverage](#)

Configuration file example for evaluation of pedestals:

```

[psana]
modules = ImgAlgos.CSPadArrAverage
files   = <path-to-the-dark-run-file>.xtc

[ImgAlgos.CSPadArrAverage]
source = DetInfo(CxiDs1.0:Cspad.0)
key    =
avefile = cspad-pedestals-ave.dat
rmsfile = cspad-pedestals-rms.dat
print_bits = 15
evts_stage1 = 100
evts_stage2 = 100
gate_width1 = 100
gate_width2 = 10

```

Configuration file example for evaluation of background:

```
[psana]
files      = <path-to-the-background-run-file>.xtc
modules    = cspad_mod.CsPadCalib ImgAlgos.CSPadArrAverage
skip-events = 500
events     = 1000000

[cspad_mod.CsPadCalib]
inputKey    =
outputKey   = calibrated
doPedestals = yes
doPixelStatus = no
doCommonMode = no

[ImgAlgos.CSPadArrAverage]
source      = DetInfo(CxiDs1.0:Cspad.0)
key         = calibrated
avefile     = cspad-background-ave.dat
rmsfile     = cspad-background-rms.dat
print_bits  = 15
```

Images of the CSPad arrays for averaged and rms values, respectively, in one of the CXI runs:

-

Example for module `ImgAlgos::CSPadBkgdSubtract`

See [Module `ImgAlgos::CSPadBkgdSubtract`](#)

```
[psana]
files = /reg/d/psdm/<instrument>/<experiment>/xtc/<file-name-1>.xtc \
        /reg/d/psdm/<instrument>/<experiment>/xtc/<file-name-2>.xtc \
        ...
        /reg/d/psdm/<instrument>/<experiment>/xtc/<file-name-N>.xtc
skip-events = 500
events      = 10
modules     = cspad_mod.CsPadCalib ImgAlgos.CSPadBkgdSubtract

[cspad_mod.CsPadCalib]
inputKey    =
outputKey   = calibrated
doPedestals = yes
doPixelStatus = no
doCommonMode = no

[ImgAlgos.CSPadBkgdSubtract]
source      = DetInfo(CxiDs1.0:Cspad.0)
inputKey    = calibrated
outputKey   = bkgd_subtracted
bkgd_fname  = <the-file-name-with-background-array>
norm_sector = 0
print_bits  = 3
```

The file with the background array, `bkgd_fname`, was obtained by averaging 1000 events using module [CSPadArrAverage](#). Subtraction is done with normalization for `norm_sector=0`.

Event image and pixel amplitude spectrum before and after the background subtraction are shown in plots:

-

-

Other event with better subtracted background:

-

Example for Module `ImgAlgos::CSPadMaskApply`

See [Module `ImgAlgos::CSPadMaskApply`](#)

The array for mask contains zeros and ones for masked and passed pixels, respectively, and has a shape of full-size CSPad array `4*8*185388`.

For example, it can be generated by the command

```
./MakePixelMask.py <input-background-cspad-arr-file-name> <threshold> <output-file-name>
```

for the averaged background amplitude array <input-background-cspad-arr-file-name> obtained as a result of `ImgAlgos::CSPadArrAverage` module.

Plots show the averaged background, and the mask arrays generated from this background for three thresholds 10, 20, and 30 EDU:

—
The best results in filtering can be achieved in combination of modules:

```
modules = cspad_mod.CsPadCalib \
          ImgAlgos.CSPadBkgdSubtract \
          ImgAlgos.CSPadMaskApply \
          ...
[cspad_mod.CsPadCalib]
inputKey      =
outputKey     = calibrated_arr
doPedestals   = yes
doPixelStatus = no
doCommonMode  = no

[ImgAlgos.CSPadBkgdSubtract]
source        = DetInfo(CxiDs1.0:Cspad.0)
inputKey      = calibrated_arr
outputKey     = bkgd_subtracted_arr
bkgd_fname    = ana-cxi49012/cspad-cxi49012-r0025-background-ave.dat
norm_sector   = 0
print_bits    = 0

[ImgAlgos.CSPadMaskApply]
source        = DetInfo(CxiDs1.0:Cspad.0)
inkey         = bkgd_subtracted_arr
outkey        = masked_arr
mask_fname    = ana-cxi49012/cspad-cxi49012-r0025-mask-40.dat
masked_amp    = 0
print_bits    = 3
mask_control_bits = 1
```

where

`cspad_mod.CsPadCalib` - subtracts the pedestals from raw CSPad data,

`ImgAlgos.CSPadBkgdSubtract` - subtracts the background,

`ImgAlgos.CSPadMaskApply` - apply the mask.

In the test with images for background represented by the water and solvent rings this filter provides the background suppression factor about 100.

The background images that still pass this filter have significantly larger intensity with respect to averaged background:

—
Input parameter `mask_control_bits` allows to control masking regions of 2x1. For example, if all edges need to be masked, then use `mask_control_bits = 15`, which gives image array like:

•
where red regions/lines of pixels of amplitude=8 are masked.

Example for module `ImgAlgos::CSPadArrNoise`

See [Module `ImgAlgos::CSPadArrNoise`](#)


```
[psana]
modules = cspad_mod.CsPadCalib ImgAlgos.CSPadArrNoise

files = /reg/d/psdm/cxi/cxi49012/xtc/e158-r0020-s00-c00.xtc \
        /reg/d/psdm/cxi/cxi49012/xtc/e158-r0020-s01-c00.xtc \
        /reg/d/psdm/cxi/cxi49012/xtc/e158-r0020-s02-c00.xtc \
        /reg/d/psdm/cxi/cxi49012/xtc/e158-r0020-s03-c00.xtc \
        /reg/d/psdm/cxi/cxi49012/xtc/e158-r0020-s04-c00.xtc \
        /reg/d/psdm/cxi/cxi49012/xtc/e158-r0020-s05-c00.xtc

#skip-events = 1000
events = 10

[cspad_mod.CsPadCalib]
inputKey =
outputKey = calibrated
doPedestals = yes
doPixelStatus = no
doCommonMode = no

[ImgAlgos.CSPadArrNoise]
source = DetInfo(CxiDs1.0:Cspad.0)
key = calibrated
statusfile = ana-cxi49012/cspad-cxi49012-r0020-noise-status.dat
maskfile = ana-cxi49012/cspad-cxi49012-r0200-noise-mask.dat
print_bits = 255
rmin = 3
dr = 1
SoNThr = 3
frac_noisy_imgs = 0.15
```

Index map in median algorithm for rmin=3, dr=1:

```
CSPadArrNoise::printMatrixOfIndexesForMedian():
 0 0 0 0 1 0 0 0 0
 0 0 1 1 1 1 1 0 0
 0 1 0 0 0 0 0 1 0
 0 1 0 0 0 0 0 1 0
 1 1 0 0 + 0 0 1 1
 0 1 0 0 0 0 0 1 0
 0 1 0 0 0 0 0 1 0
 0 0 1 1 1 1 1 0 0
 0 0 0 0 1 0 0 0 0
```

Pixel status (fraction of events where $S/N > \text{SoNThr}$):

-

For cspad-cxi49012-r0020 with parameters from configuration file (frac_noisy_imgs=0.15) we get, depending on number of events:

Nnoisy, Ntotal, Nnoisy/Ntotal pixels =94585 2296960 0.041 for 10 events

Nnoisy, Ntotal, Nnoisy/Ntotal pixels =2112 2296960 0.00092 for 100 events

Pixel mask for noisy pixels with $|S/N| > \text{SoNThr}$:

-

Example for Module `ImgAlgos::CSPadArrPeakFinder`

See [Module `ImgAlgos::CSPadArrPeakFinder`](#)

```

[psana]
files = \
    /reg/d/psdm/cxi/cxi49012/xtc/e158-r0150-s00-c00.xtc \
    /reg/d/psdm/cxi/cxi49012/xtc/e158-r0150-s01-c00.xtc \
    /reg/d/psdm/cxi/cxi49012/xtc/e158-r0150-s02-c00.xtc \
    /reg/d/psdm/cxi/cxi49012/xtc/e158-r0150-s03-c00.xtc \
#    /reg/d/psdm/cxi/cxi49012/xtc/e158-r0150-s04-c00.xtc \ suddenly it became unavailable...
    /reg/d/psdm/cxi/cxi49012/xtc/e158-r0150-s05-c00.xtc

#skip-events = 1000
#events      = 200

modules = cspad_mod.CsPadCalib \
          ImgAlgos.CSPadMaskApply \
          ImgAlgos.CSPadArrPeakFinder

[cspad_mod.CsPadCalib]
inputKey      =
outputKey     = calibrated
doPedestals  = yes
doPixelStatus = no
doCommonMode = no

[ImgAlgos.CSPadMaskApply]
source        = DetInfo(CxiDs1.0:Cspad.0)
inkey         = calibrated
outkey        = masked_arr
mask_fname    = ana-cxi49012/cspad-cxi49012-r0150-mask-badregs.dat
#mask_fname   = ana-cxi49012/cspad-cxi49012-r0150-mask-bkgd.dat
#mask_fname   = ana-cxi49012/cspad-cxi49012-r0150-mask-rects.dat
masked_amp    = 8
print_bits    = 1
mask_control_bits = 15

[ImgAlgos.CSPadArrPeakFinder]
source        = DetInfo(CxiDs1.0:Cspad.0)
key           = masked_arr
key_peaks_out = peaks

hot_pix_mask_inp_file = ana-cxi49012/cspad-cxi49012-r0150-noise-mask.dat
hot_pix_mask_out_file = ana-cxi49012/cspad-cxi49012-r0150-noise-mask-out.dat
frac_noisy_evts_file  = ana-cxi49012/cspad-cxi49012-r0150-noise-frac.dat

evt_file_out = tmp/cspad-ev-

rmin          = 3
dr            = 1
SoNThr        = 3
frac_noisy_imgs = 0.1

peak_npix_min = 4
peak_npix_max = 25
peak_amp_tot_thr = 100.

event_npeak_min = 10
event_amp_tot_thr = 1000.

nevents_mask_update = 100
nevents_mask_accum  = 50

selection_mode = SELECTION_ON
out_file_bits  = 15
print_bits     = 512

```

Results:

```
[info:TimeInterval::startTime] Start time: 2012-06-12 15:32:02
[info:ImgAlgos.CSPadArrPeakFinder] N processed events = 1000 N selected = 55 Fraction of selected = 0.055
[info:ImgAlgos.CSPadArrPeakFinder] N processed events = 2000 N selected = 62 Fraction of selected = 0.031
[info:ImgAlgos.CSPadArrPeakFinder] N processed events = 3000 N selected = 81 Fraction of selected = 0.027
[info:ImgAlgos.CSPadArrPeakFinder] N processed events = 4000 N selected = 95 Fraction of selected = 0.02375
[info:ImgAlgos.CSPadArrPeakFinder] N processed events = 5000 N selected = 150 Fraction of selected = 0.03
[info:ImgAlgos.CSPadArrPeakFinder] N processed events = 6000 N selected = 265 Fraction of selected =
0.0441667
[info:ImgAlgos.CSPadArrPeakFinder] N processed events = 7000 N selected = 404 Fraction of selected =
0.0577143
[info:ImgAlgos.CSPadArrPeakFinder] ===== JOB SUMMARY =====
[info:TimeInterval::stopTime] Time to process 7945 events is 3747.48 sec, or 0.471678 sec/event
```

Selected events

ev-007713:

-

ev-008944:

-

Next example shows how to use `ImgAlgos::CSPadArrPeakFinder` in combination with python module `pyimgalgos.ex_peaks_nda`:

```

[psana]
#skip-events = 50
events      = 10
files = exp=cxi83714:run=136

modules = cspad_mod.CsPadCalib \
          ImgAlgos.CSPadArrPeakFinder \
          pyimgalgos.ex_peaks_nda
#          EventKeys

[cspad_mod.CsPadCalib]
inputKey      =
outputKey     = calibrated
doPedestals  = yes
doPixelStatus = yes
doCommonMode = yes

[ImgAlgos.CSPadArrPeakFinder]
source      = DetInfo(CxiDs1.0:Cspad.0)
key         = calibrated
key_peaks_out =
key_peaks_nda = peaks_nda
hot_pix_mask_inp_file = ana-cxi83714/cspad-cxi83714-r0136-noise-mask-ini.dat
hot_pix_mask_out_file = ana-cxi83714/cspad-cxi83714-r0136-noise-mask-out.dat
frac_noisy_evts_file = ana-cxi83714/cspad-cxi83714-r0136-noise-frac.dat
evt_file_out      = cxi83714/cspad-ev-
rmin              = 8
dr                = 1
SoNThr_noise      = 3
SoNThr_signal     = 4
frac_noisy_imgs   = 0.9
peak_npix_min     = 3
peak_npix_max     = 500
peak_amp_tot_thr  = 0.
peak_SoN_thr      = 5.
event_npeak_min   = 5
event_npeak_max   = 1000
event_amp_tot_thr = 0.
nevents_mask_update = 0
nevents_mask_accum = 50
selection_mode    = SELECTION_ON
out_file_bits     = 15
print_bits        = 3681

[pyimgalgos.ex_peaks_nda]
source      = DetInfo(CxiDs1.0:Cspad.0)
key_in      = peaks_nda
print_bits  = 255

```

Example for module `ImgAlgos::CSPadArrPeakAnalysis`

See [Module `ImgAlgos::CSPadArrPeakAnalysis`](#)

Example of the psana configuration file:

```

modules = cspad_mod.CsPadCalib \
          ImgAlgos.CSPadMaskApply \
          ImgAlgos.CSPadArrPeakFinder \
          ImgAlgos.CSPadArrPeakAnalysis

# ...configuration parameters of other modules...

[ImgAlgos.CSPadArrPeakAnalysis]
source      = DetInfo(CxiDs1.0:Cspad.0)
key         = peaks
print_bits  = 7
fname_root  = file.root

```

After execution in psana the file `.root` containing histogram(s) and ntuple(s) will be produced. Then, auxiliary script in root, running by the commend

```
root -q -f proc.C
```

produces the plots with histograms:

-

Example for TimeStampFilter and XtcOutputModule

This example demonstrates how to run psana with the **time stamp filter** and **event writer** in xtc file.

Both modules are available in psana library and they need only to be described in the configuration file. For example, the configuration file `tstamp-filter-and-event-writer.cfg` may looks like:

```

[psana]
files      = /reg/d/psdm/<INSTRUMENT>/<experiment>/xtc/e158-r0021-s00-c00.xtc \
             /reg/d/psdm/<INSTRUMENT>/<experiment>/xtc/e158-r0021-s01-c00.xtc \
             /reg/d/psdm/<INSTRUMENT>/<experiment>/xtc/e158-r0021-s02-c00.xtc \
             /reg/d/psdm/<INSTRUMENT>/<experiment>/xtc/e158-r0021-s03-c00.xtc \
             /reg/d/psdm/<INSTRUMENT>/<experiment>/xtc/e158-r0021-s04-c00.xtc \
             /reg/d/psdm/<INSTRUMENT>/<experiment>/xtc/e158-r0021-s05-c00.xtc

skip-events = 10
events      = 100
modules     = ImgAlgos.TimeStampFilter  PSXtcOutput.XtcOutputModule

[PSXtcOutput.XtcOutputModule]
dirName     = ./test_out

[ImgAlgos.TimeStampFilter]
tsinterval  = 2012-02-02 18:17:00.409143728 / 2012-02-02 18:17:00.525853474
filterIsOn  = true
print_bits  = 11

```

Command to run psana is:

```
psana -c ./tstamp-filter-and-event-writer.cfg
```

For this configuration file psana will skip 10 events and loop over the next 100 events from the `/reg/d/psdm/<INSTRUMENT>/<experiment>/xtc/e158-r0021-s0*-c00.xtc` files and run consecutively modules = `ImgAlgos.TimeStampFilter PSXtcOutput.XtcOutputModule`. Parameters of these modules are described in the bottom part of the configuration file. Module `TimeStampFilter` passes events from the specified time interval and prints some useful information. Module `XtcOutputModule` will write passed events in the file with auto-generated name `./test_out/e158-r0021.xtcf`.

See also: [Module ImgAlgos::TimeStampFilter](#) and [PSXtcOutput::PSXtcOutput](#)

Example for module ImgAlgos::UsdUsbEncoderFilter

Example of configuration file:

```
[psana]
files = exp=amo75113:run=156
events = 10

modules = ImgAlgos.UsdUsbEncoderFilter

[ImgAlgos.UsdUsbEncoderFilter]
source = DetInfo(:USDUSB)
mode = 1
ifname = tstamp-code-in.txt
ofname = tstamp-code-out.txt
print_bits = 31
```

See also: [Module ImgAlgos::UsdUsbEncoderFilter](#)

Examples for Package pyimgalgos

Package **pyimgalgos** contains python modules which work with both frameworks **pyana** and **psana**.

Configuration file for **pyana** and/or **psana** should have relevant sections with parameters for [pyana] and/or [psana]. Alean section is ignored in each framework at run time. This is the only difference between two frameworks in the configuration file. All module descriptions are the same for two frameworks, as shown in examples below.

See description of modules in [Package pyimgalgos](#).

Example of configuration file for CSPAD

File py-xcs72913-r0049-cspad.cfg:

```

# Run this script:
# psana -c py-xcs72913-r0049-cspad.cfg
# pyana -c py-xcs72913-r0049-cspad.cfg
#
# Useful commands:
# psana -m EventKeys -n 5 /reg/d/psdm/xcs/xcs72913/xtc/e265-r0049-*
# psana -m psana_examples.dump_cspad -n 5 exp=xcs72913:run=49
# pyana -m pyana_examples.dump_cspad -n 5 /reg/d/psdm/xcs/xcs72913/xtc/e265-r0049-*

[pyana]
files = /reg/d/psdm/xcs/xcs72913/xtc/e265-r0049-s00-c00.xtc /reg/d/psdm/xcs/xcs72913/xtc/e265-r0049-s04-c00.
xtc /reg/d/psdm/xcs/xcs72913/xtc/e265-r0049-s05-c00.xtc
num-events = 5
#skip-events = 0
#num-cpu = 1
verbose = 1 ; logging output: 0-nothing?, 1+INFO, 2+DEBUG, ...
modules = pyimgalgos.tahometer pyimgalgos.cspad_arr_producer pyimgalgos.cspad_image_producer pyimgalgos.
image_save_in_file

[psana]
files = exp=xcs72913:run=49
events = 5
#skip-events = 0
modules = pyimgalgos.tahometer pyimgalgos.cspad_arr_producer pyimgalgos.cspad_image_producer pyimgalgos.
image_save_in_file
verbose = 1

[pyimgalgos.tahometer]
dn = 10
print_bits = 255

[pyimgalgos.cspad_arr_producer]
#source = -|Cspad-*
source = XcsEndstation-0|Cspad-0
data_type = double
#data_type = float
#data_type = unsigned
#data_type = uint16
val_miss = 0
key_out = cspad_array
print_bits = 1

[pyimgalgos.cspad_image_producer]
calib_dir = /reg/d/psdm/xcs/xcs72913/calib/CsPad::CalibV1/XcsEndstation.0:Cspad.0/
key_in = cspad_array
key_out = cspad_image
print_bits = 1

# Supported output file formats tiff, gif, png, eps, jpg, jpeg, txt, npy(default), npz
[pyimgalgos.image_save_in_file]
key_in = cspad_image
ofname = img-for-cspad.txt
print_bits = 255

```

To run this script use command

```
psana -c py-xcs72913-r0049-cspad.cfg
```

or

```
pyana -c py-xcs72913-r0049-cspad.cfg
```

Example of configuration file for CSPAD2x2

File py-meca6113-r0028-cspad2x2.cfg

```
# Run this script:
# psana -c py-meca6113-r0028-cspad2x2.cfg
# pyana -c py-meca6113-r0028-cspad2x2.cfg
#
# Useful commands:
# psana -m EventKeys -n 5 /reg/d/psdm/mec/meca6113/xtc/e332-r0028-s03-c00.xtc
# psana -m psana_examples.dump_cspad -n 5 exp=meca6113:run=28
# pyana -m pyana_examples.dump_cspad -n 5 /reg/d/psdm/mec/meca6113/xtc/e332-r0028-s03-c00.xtc

[pyana]
files = /reg/d/psdm/mec/meca6113/xtc/e332-r0028-s03-c00.xtc
num-events = 5
#skip-events = 0
#num-cpu = 1
verbose = 0 ; logging output: 0-nothing?, 1+INFO, 2+DEBUG, ...
modules = pyimgalgos.tahometer pyimgalgos.cspad_arr_producer pyimgalgos.cspad_image_producer pyimgalgos.
image_save_in_file

[psana]
files = exp=meca6113:run=28
events = 5
#skip-events = 0
modules = pyimgalgos.tahometer pyimgalgos.cspad_arr_producer pyimgalgos.cspad_image_producer pyimgalgos.
image_save_in_file

[pyimgalgos.tahometer]
dn = 10
print_bits = 255

[pyimgalgos.cspad_arr_producer]
#source = -|Cspad-*
source = MecTargetChamber-0|Cspad2x2-3
#data_type = double
data_type = float
#data_type = unsigned
#data_type = uint16
#data_type = uint32
val_miss = 0
key_out = cspad2x2_array
print_bits = 255

[pyimgalgos.cspad_image_producer]
calib_dir = /reg/d/psdm/mec/meca6113/calib/CsPad2x2::CalibV1/MecTargetChamber.0:Cspad2x2.3/
key_in = cspad2x2_array
key_out = cspad2x2_image
print_bits = 1

# Supported output file formats tiff, gif, png, eps, jpg, jpeg, txt, npy(default), npz
[pyimgalgos.image_save_in_file]
key_in = cspad2x2_image
ofname = img-for-cspad2x2.tiff
print_bits = 255
```

To run this script use command

```
psana -c py-meca6113-r0028-cspad2x2.cfg
```

or

```
pyana -c py-meca6113-r0028-cspad2x2.cfg
```


References

[psana - Module Catalog](#)