

Administrative rules for LCLS detectors

Content

- [Content](#)
- [Detector name labels](#)
- [Optical metrology](#)
 - [epix10ka2m](#)
 - [Use the same measurement pattern for all epix10ka2m detectors](#)
 - [Camera orientation in optical metrology](#)
 - [Measurement order](#)
 - [MEC Quad Orientation](#)
- [Label on connectors and cables](#)
- [Detector name in DAQ](#)
- [Summary](#)

In this note epix10ka2m is the specific example, but the conventions adopted would apply to any future multi-panel detector.

Detector name labels

Add **label with unique detector name** on each camera body

- preliminary agreed on unique detector names: 2M.0, 2M.1, 2M.2, ...
- if expected number of detectors more than 10 it is better to use 2M.00, 2M.01, 2M.02, ... (better for files sorting in catalogs).

Why it is important? We need to distinguish detectors in optical metrology and geometry constant files. It does not matter for intensity correction, where detector will be identified by the panels Id from data.

Optical metrology

epix10ka2m



Use the same measurement pattern for all epix10ka2m detectors

- orientation of all cameras in optical metrology should always be the same
- order of measured points should always be the same

Why it is important? We have to unambiguously map panels in optical metrology and in data from DAQ. Different detector orientation or different order of measurement need in special treatment in the processing script, that is error-prone and always needs in verification on real data.

Camera orientation in optical metrology

Preferred optical metrology pattern from 2018-11-18:



- orientation of camera should be consistent with its beam-view in hutch, quad **Q0 in the top left corner**.
 - orientation of quad labels should be human readable (from left to right) in hutch position; on current photo labels need to be rotated upside-down. Right photo shows how embarrassing labels look in the hutch.
 - add labels for vector of gravity or x-y for optical metrology in order to exclude wrong orientation
 - microscope origin should stay the same during the measurements.
Origin location does not matter, but it can optionally be labeled somewhere for uniformity of metrology data.

Measurement order

- measurement order of points per quad should be done in style of Chris Kenney, as shown on photo, **clock-wise order for panels and corners, beginning from bottom-left**.
- metrology data is saved in spreadsheets xlsx or text format.
- units - signed micrometers, without decimal dot.
- example of data:

Example of optical metrology data in text format

```
quad 0
point  X      Y      Z
1       0       0       0
2      403     39683  -147
3     37088    39308  -472
4     36680   -378    -326
5       322    40433  -141
6       366    80129  -302
7     37056    80079  -657
8     37007    40391  -452
9     42053    40578  -470
10     42046    80274  -669
11     78743    80269  -961
12     78733    40576  -725
13     42083     219   -297
14     42078    39472  -413
15     78774    39473  -711
16     78767     214   -640

quad 1
point  X      Y      Z
1     81539   -4620  -731
2     81352   32062  -814
3    121034   32267  -1079
4    121225  -4420   -1015
5     81390   37579  -882
6     81497   74238  -993
7    121177   74148  -1219
8    121073   37468  -1127
9    122083   37596  -1104
10   122171   74281  -1329
11   161853   74189  -1620
12   161768   37499  -1433
13   122081  -4398   -916
14   122183   32283  -988
15   161877   32161  -1333
16   161765  -4520  -1286

Quad 2
...
Quad 3
...
```

MEC Quad Orientation

< Insert documentation on MEC Quad Orientation here >

Label on connectors and cables

- All **labels on detector connectors and cables** should show **THE ONLY WAY** to connect camera in hutch.

Detector name in DAQ

- Latest incident revealed by Philip on 2020-07-16 - "We've been mostly using the 2M.1 camera but the daq calls it 2M.0."
- Detector name like DetLab.0:Epix10ka2M.0 is used to access data in xtc file.
- it should be consistent with calibration constants, like
/reg/d/psdm/det/detdaq18/calib/Epix10ka2M::CalibV1/DetLab.0:Epix10ka2M.0/pedestals/<run-range>.data
- To prevent mess with names it would be nice to use the same detector name on camera body (2M.0) and in DAQ (DetLab.0:Epix10ka2M.0).

Summary

- Engrave **detector name on camera body**.
- Add **labels for unambiguous detector orientation**:
 - vector of gravity,
 - x-y (and optional origin) for optical metrology,
 - all labels should be left-to-right readable in the hutch (now for 2M.0 they are upside-down).
- Use the same detector orientation in optical metrology, **Q0 is in upper left corner** in the hutch and in microscope x-y Cartesian frame.
 - Origin location does not matter, but it is better to label for uniformity of constants.
- Use the **same order of measuring points per quad**, as numerated on photo.
- Plug-in **properly labeled cables to the properly labeled camera connectors**.
- Use **the same detector name on camera body** (2M.0) **and in DAQ** (DetLab.0:Epix10ka2M.0).