

EPIX10KA Data Access

- Calibration constants
- Run, Step, and Event Loops
- Methods of det.raw
 - Methods for AMI interface
 - raw
 - calib
 - image
 - segments
 - Other useful methods hidden for AMI
 - _mask and _mask_comb
 - Separate masks
 - Pixel coordinate arrays
- References

Calibration constants

- [Jungfrau and Epix10ka Calibration](#)

Run, Step, and Event Loops

```
class Arguments:  
    expt = 'ueddaq02'  
    run = 66  
    evtmax = 5  
    detname = 'epixquad'  
  
args = Arguments()  
  
from psana.pyalgos.generic.NDArrUtils import info_ndarr  
from psana import DataSource  
ds = DataSource(exp=args.expt, run=args.run, dir=f'/cds/data/psdm/{args.expt[:3]}/{args.expt}/xtc')  
  
for irun,run in enumerate(ds.runs()):  
    print('\n==== %02d run: %d exp: %s detnames: %s' % (irun, run.runnum, run.expt, ','.join(run.detnames)))  
  
    print('make %s detector object' % args.detname)  
    det = run.Detector(args.detname)  
  
    for istep,step in enumerate(run.steps()):  
        print('\nStep %ld' % istep)  
  
        for ievt,evt in enumerate(step.events()):  
            if ievt>args.evtmax: exit('exit by number of events limit %d' % args.evtmax)  
  
            print('%s\nEvent %04d' % (80*' ',ievt))  
            segs = det.raw.segments(evt)  
            raw = det.raw.raw(evt)  
  
            print(info_ndarr(segs, 'segments'))  
            print(info_ndarr(raw, 'raw '))
```

output of the above script

```
==== 00 run: 66 exp: ueddaq02 detnames: timing,epixquad
make epixquad detector object

Step 0

Event 0000
segments shape:(4,) size:4 dtype:uint16 [0 1 2 3]
raw shape:(4, 352, 384) size:540672 dtype:uint16 [3304 3385 3401 3407 3291...]

Event 0001
segments shape:(4,) size:4 dtype:uint16 [0 1 2 3]
raw shape:(4, 352, 384) size:540672 dtype:uint16 [3278 3381 3385 3393 3298...]

Event 0002
segments shape:(4,) size:4 dtype:uint16 [0 1 2 3]
raw shape:(4, 352, 384) size:540672 dtype:uint16 [3317 3431 3431 5256 5790...]

Event 0003
segments shape:(4,) size:4 dtype:uint16 [0 1 2 3]
raw shape:(4, 352, 384) size:540672 dtype:uint16 [3325 3439 3426 3462 3361...]

Event 0004
segments shape:(4,) size:4 dtype:uint16 [0 1 2 3]
raw shape:(4, 352, 384) size:540672 dtype:uint16 [3328 3436 3434 3456 3348...]


Event 0005
segments shape:(4,) size:4 dtype:uint16 [0 1 2 3]
raw shape:(4, 352, 384) size:540672 dtype:uint16 [3335 3428 3426 3454 3347...]
exit by number of events limit 5
```

Methods of det.raw

Methods for AMI interface


raw

```
raw = det.raw.raw(evt)
```



calib


```
clb = det.raw.calib(evt)
```



```
user_mask = np.ones_like(det.raw.raw(evt), dtype=DTYPE_MASK) #np.uint8
user_mask[0,100:150,200:250] = 0
```

```
clb = det.raw.calib(evt, cmpars=(7,2,100,10), mbits=0o7, mask=user_mask, edge_rows=10, edge_cols=10, center_rows=5, center_cols=5, **kwargs)
```

common mode correction, user defined mask, status mask and mask of edges are applied:



image

```
img = det.raw.image(evt, **kwargs)
```

```
img = det.raw.image(evt)
```

```
img = det.raw.image(evt, nda=arr, pix_scale_size_um=args.pscsize, mapmode=args.mapmode)
```

segments

```
segs = det.raw.segments(evt)
```


Other useful methods hidden for AMI

_mask and _mask_comb

```
msk = det.raw._mask(calib=False, status=False, edges=False, **kwargs)
```

```
msk = det.raw._mask_comb(mbits=0o7, mask=None, edge_rows=10, edge_cols=10, center_rows=5, center_cols=5)
```

where mbits &1 - calib, &2 - status, &4 -edges , &8 - neighbors



Separate masks

```

o = det.raw

m = o._mask_default(dtype=DTYPE_MASK)
m = o._mask_calib_or_default(dtype=DTYPE_MASK)
m = o._mask_from_status(**kwa)
m = o._mask_edges(edge_rows=10, edge_cols=10, center_rows=5, center_cols=5, **kwa)

```

Pixel coordinate arrays

```

a = o._pixel_coord_indexes(**kwa) # 'pix_scale_size_um',None; 'xy0_off_pix',None; # do_tilt',True; 'cframe',0
a = o._pixel_coords(**kwa) # do_tilt',True; 'cframe',0

```

References

- EPIX10KA2M References
- Jungfrau and Epix10ka Calibration
- Common mode correction algorithms
- <https://github.com/slac-lcls/lcls2/blob/master/psana/detector/areadetector.py>