

Jungfrau dark processing

Content

- [Content](#)
- [Intro](#)
- [Repository organization](#)
 - [Calibration file naming conventions](#)
 - [Log file naming conventions](#)
- [Dark processing](#)
- [Merging and deployment of constants](#)
- [Example for parallel processing with time comparison](#)
- [Optional parameters on 2021-12-10](#)
- [References](#)

Intro

In April 2021 we introduce a new approach to dark runs processing for Jungfrau detectors. This approach is pretty much similar to epix10ka dark calibration based on **per-panel** processing of constants. Entire dark calibration is split for two stages. At first stage specialized dark runs for three gain modes or regular runs with "drop-shots" are processed and results per-panel are saved in the repository. At second stage, per-panel constants are merged for particular detector or set of panels (for DRP) and deployed in the regular calibration directory. This approach allows to manage each panel and its gain ranges separately.

Repository organization

- `/reg/g/psdm/detector/gains/jungfrau/panels/` # repository for jungfrau per-panel constants
 - `190408-181206-50c246df50010d/` # e.g. - for one of the panel ids
 - `rms/` # subdirectory for calibration constant's type
 - `status/`
 - `dark_max/`
 - `dark_min/`
 - `pedestals/`
 - `jungfrau_0001_20201201073354_cxilu9218_r0230_pedestals_gm0-Normal.dat` # e.g. - panel constants for gain mode 0: Normal
 - `jungfrau_0001_20201201073425_cxilu9218_r0231_pedestals_gm1-ForcedGain1.dat`
 - `jungfrau_0001_20201201073453_cxilu9218_r0232_pedestals_gm2-ForcedGain2.dat`
 - `jungfrau_0001_20201201085333_cxilu9218_r0238_pedestals_gm0-Normal.dat`
 - `jungfrau_0001_20201201085333_cxilu9218_r0238_pedestals.txt` # e.g. - per-panel. constants
 - `merged for three gain modes valid for runs >=238`
 - `merge_tmp/`
 - `CxiDs1.0:Jungfrau.0-20201201085333-cxilu9218-r0238-pedestals.txt` # merged constants for detector in run 238
 - `logs/`
 - `2021_log_jungfrau_dark_proc.txt` # log files with single record per job for each script
 - `2021_log_jungfrau_deploy_constants.txt`
 - `2021/` # annual directory for job
 - `submission commands`
 - `2021-04-23T142354_log_jungfrau_dark_proc_dubrovin.txt` # log file for each job
 - `2021-04-23T154102_log_jungfrau_deploy_constants_dubrovin.txt`

Calibration file naming conventions

File name example: `jungfrau_0001_20201201073354_cxilu9218_r0230_pedestals_gm0-Normal.dat`

Fields meaning: `<panel-type>_<panel-unique-numeric-alias>_<run-time-stamp>_<experiment>_r<run-number>_<constants-type>_gm<gain-mode-index-and-name>.dat`

Log file naming conventions

File name example: `2021-04-23T142354_log_jungfrau_dark_proc_dubrovin.txt`

Fields meaning: `<job-submission-time-stamp>_<script-name>_<login-name>.txt`

Dark processing

`jungfrau_dark_processing`

Processes jungfrau dark data

`jungfrau_dark_proc -d <dataset> -s <source> ...`

Ex.1: `jungfrau_dark_proc -d exp=cxilu9218:run=242,243,244:smd -s CxiDs1.0:Jungfrau.0 # regular dark processing`

Ex.2: `jungfrau_dark_proc -d exp=cxilu9218:run=242,243,244:smd -s CxiDs1.0:Jungfrau.0 -l1 # speed-up for single panel`

Ex.3: `jungfrau_dark_proc -d exp=cxilu9218:run=238:smd -s CxiDs1.0:Jungfrau.0 --evcode 162 # use drop-shots as dark events`

Help: `jungfrau_dark_proc -h`

Merging and deployment of constants

`jungfrau_merge_constants`

Deployment jungfrau calibration parameters

`jungfrau_deploy_constants -e <experiment> -d <detector> -r <run-number> [-D] [-L <logging-mode>] [...]`

TEST COMMAND:

Ex.1: `jungfrau_deploy_constants -e cxilu9218 -d CxiDs1.0:Jungfrau.0 -r238`

Ex.2: `jungfrau_deploy_constants -e cxilu9218 -d CxiDs1.0:Jungfrau.0 -r238 -D -c ./calib`

REGULAR COMMAND:

Ex.3: `jungfrau_deploy_constants -e cxilu9218 -d CxiDs1.0:Jungfrau.0 -r238 -D`

Help: `jungfrau_deploy_constants -h`

Example for parallel processing with time comparison

Dataset `exp=xpplw4319:run=1` contains three calib-cycles with dark data for 2-panel `XppEndstation.0:Jungfrau.0`

process all steps (all gain ranges) and all panels takes 306 sec on psanagpu103 with command

`jungfrau_dark_proc -d exp=xpplw4319:run=1:smd:stream=0-79 -s XppEndstation.0:Jungfrau.0`

process all steps for a single panel (`--segind=1` is selected from 0,1) takes 103 sec

`jungfrau_dark_proc -d exp=xpplw4319:run=1:smd:stream=0-79 -s XppEndstation.0:Jungfrau.0 --segind=1`

process a single step (`--stepnum=2` is selected from 0,1,2) for all panels takes 83 sec with command

`jungfrau_dark_proc -d exp=xpplw4319:run=1:smd:stream=0-79 -s XppEndstation.0:Jungfrau.0 --stepnum=2`

process a single step (`--stepnum=2`) for a single panel (`--segind=1`) takes 44 sec with command

`jungfrau_dark_proc -d exp=xpplw4319:run=1:smd:stream=0-79 -s XppEndstation.0:Jungfrau.0 --stepnum=2 --segind=1`

Then results have to be deployed in the calib directory using command

`jungfrau_deploy_constants -e xpplw4319 -d XppEndstation.0:Jungfrau.0 -r1 -D`

Optional parameters on 2021-12-10

jungfrau_dark_proc -h

Processes jungfrau dark data

Options:

- h, --help show this help message and exit
- d DSNNAME, --dsname=DSNAME
dataset name, default = None
- s SOURCE, --source=SOURCE
input ndarray source name, default = None
- n EVENTS, --events=EVENTS
maximal number of events total (in runs, steps),
default = 100000
- m EVSKIP, --evskip=EVSKIP
number of events to skip in the beginning of each
step, default = 0
- e EVSTEP, --evstep=EVSTEP
maximal number of events to process in each step,
default = 2000
- b INT_LO, --int_lo=INT_LO
intensity low limit, default = 1
- t INT_HI, --int_hi=INT_HI
intensity high limit, default = 16000
- B RMS_LO, --rms_lo=RMS_LO
rms low limit, default = 0.001
- T RMS_HI, --rms_hi=RMS_HI
rms high limit, default = 16000
- F FRACLM, --fraclm=FRACLM
allowed fraction limit, default = 0.1
- p PLOTIM, --plotim=PLOTIM
control bit-word to plot images, default = 0
- D INTNLO, --intnlo=INTNLO
number of sigma from mean for low limit on INTENSITY,
default = 6.0
- U INTNHI, --intnhi=INTNHI
number of sigma from mean for high limit on INTENSITY,
default = 6.0
- L RMSNLO, --rmsnlo=RMSNLO
number of sigma from mean for low limit on RMS,
default = 6.0
- H RMSNHI, --rmsnhi=RMSNHI
number of sigma from mean for high limit on RMS,
default = 6.0
- c EVCODE, --evcode=EVCODE
comma separated event codes for selection as OR
combination, any negative code inverts selection,
default = None
- u, --upload upload files in calib directory, default = False
- N STEPNUM, --stepnum=STEPNUM
step/calibcycle number [0,2] or all by default,
default = None
- M STEPMAX, --stepmax=STEPMAX
maximal number of steps or all by default, default =
None
- I SEGIND, --segind=SEGIND
segment index to process, default = None
- nrecs=NRECS number of records to collect data, default = 1000
- nrecs1=NRECS1 number of records for 1st stage processing, default =
50
- fraclo=FRACLO fraction of statistics [0,1] below low limit of the
gate, default = 0.050000
- frachi=FRACHI fraction of statistics [0,1] below high limit of the
gate, default = 0.950000
- o DIRREPO, --dirrepo=DIRREPO
repository for calibration results, default =
/reg/g/psdm/detector/gains/jungfrau/panels
- logmode=LOGMODE logging mode, one of INFO, CRITICAL, WARN, WARNING,
ERROR, DEBUG, NOTSET, default = INFO

jungfrau_deploy_constants -h

Deployment jungfrau calibration parameters

Options:

- h, --help show this help message and exit
- e EXP, --exp=EXP experiment name, default = None
- d DET, --det=DET detector name, default = None
- r RUN, --run=RUN run number for beginning of the validity range, default = None
- t TSTAMP, --tstamp=TSTAMP non-default time stamp (<YYYYmddHHMMSS) or run number(<10000) for constants selection in repo. By default run time is used, default = None
- x DIRXTC, --dirxtc=DIRXTC non-default xtc directory which is used to access run start time, default = None
- o DIRREPO, --dirrepo=DIRREPO non-default repository of calibration results, default = /reg/g/psdm/detector/gains/jungfrau/panels
- c DIRCALIB, --dircalib=DIRCALIB deployment calib directory if different from standard one, default = None
- I PANINDS, --paninds=PANINDS comma-separated (str) panel indexds to generate constants for subset of panels (ex. 2 panel of 8 in the detector), default = None
- D, --deploy deploy constants to the calib dir, default = False
- L LOGMODE, --logmode=LOGMODE logging mode, one of INFO, CRITICAL, WARN, WARNING, ERROR, DEBUG, NOTSET, default = INFO

References

- [Jungfrau References](#)
- [Jungfrau and Epix10ka Calibration](#)
- [Jungfrau naming and calibration files](#)
- [EPIX10KA2M References](#)