ImgMan User Guide

ImgMan User Guide

- 1. #Overview
- 2. #Starting ImgMan from Matlab Desktop
- 3. #Viewing Live Images
 - #Image Processing Panel
- 4. #Collecting Images
 - #Saving and Loading Datasets
- 5. #Browsing Image Datasets
- 6. #Analyzing Single Image
- 7. Online Help
 - #Image Acquisition
 - #Image Analysis
 - #Image Browser

Overview

ImgMan (short for Image Management) is a suite of Matlab applications for collecting and analyzing beam images from laser cameras and profile monitors at LCLS. ImgMan consists of the Image Acquisition, Image Browser, and Image Analysis applications.

Starting ImgMan from Matlab Desktop

To start ImgMan from Matlab, type

imgAcq_main

in the Command Window.

Viewing Live Images

Mage Acquisition (DEVELOPMENT) Version 1.0.2 (04-13-20	(70		- ×
YAG01 Measure Retrieve images	Cancel Save Load Bro	owse Help	Exit Matlab
3 Background images 6 Beam images	Dataset #1 (OTR1)		0 datasets
Screen is IN In Out		Background image f	rom 04-11-2007 20:02:30.816
			Show live image

In the Image Acquisition window, select a camera and check the "Show Live Image" box. After a short delay, the window expands, and you can see the live image as well as the #Image Processing Panel, which can be found on all three main windows of ImgMan.



Image Processing Panel

Slice plane 🛛 🗙 📃	Y # Slices 1	Slice # 1 💌	Subtract background	Saved	Calculated	Gaussian 👻
Use filter(s)	Floor	Median	Crop	Auto	Custom	XMEAN 3550 VMEAN 2462
Chairing an an an aid (a)	Current	Save	Color map 🛛 jet 🛛 👻	256		XRMS 933.44 VRMS 1064.27
	Laser beam	Beam size units	um 👻		CORR 19865.89 SUM 143842019	

Image Processing Panel contains widgets for setting image processing parameters, e.g. slicing, background subtraction, filters, cropping, color mapps, annotations (centroids), unit scales, the image processing algorithm. In general, ImgMan handles changes to these parameters immediately. However, if an "Apply" button is present, you must press it before some costly image processing is started (alternatively, you can check the "Instant" Apply box).

Collecting Images

✓ Image Acquisition (DEVELOPMENT))	/ersion 1.0.2 (04-13-200)7)					- ×
YAG01 Measure	Retrieve images	Cancel	Save	Load	Browse	Help	Exit Matlab
3 Background images 6 Bea	m images		Dataset ;	#1 (OTR1)			0 datasets
Screen is IN In Ou	t				Back <u>c</u>	ground image fro	m 04-11-2007 20:02:30.816
							Show live image

In the Image Acquisition window, select a camera, specify the desired number of background and beam images, and press the "Measure" button. If you don't want to save images locally, uncheck the "Retrieve images" box. You can monitor the progress of the measurement either in the Image Acquisition window or in *cmlog*. Each measurement results in one dataset of images.

Saving and Loading Datasets

Mage Acquisition (DEVELOPMENT) Version 1.0.2 (04-13-2007)			WWWWWWWW					
YAG01 -	Measure Retrieve images	Cancel	Save	Load	Browse	Help	Exit Matlab	
Background images	6 Beam images		Dataset ;	#1 (OTR1)			0 datasets	
Screen is IN In	Out				Backg	ground image fro	m 04-11-2007 20:02:30.816	
							Show live image	

Press "Save" button to open a "Save As..." dialog and save your image datasets. Press "Load" button to open a "Load File..." dialog and load image datasets.

Browsing Image Datasets

Mage Acquisition (DEVELOPMENT) Version 1.0.2 (04-13-20	07)		- ×
YAG01 Measure Retrieve images	Cancel Save Load	Browse) Help	Exit Matlab
3 Background images 6 Beam images	Dataset #1 (OTR1)		0 datasets
Screen is IN In Out		Background i	mage from 04-11-2007 20:02:30.816
			Show live image

In the Image Acquisition window, press the "Browse..." button to open the #Image Browser. You can select a dataset, browse images, and specify image processing parameters on an #Image Processing Panel.



Analyzing Single Image



In the Image Browser window, press the "Analyze..." button to open an #Image Analysis application for the corresponding image. Image Analysis also contains an #Image Processing Panel.

Online Help For ImgMan

ImgMan (short for Image Management) is a suite of Matlab applications for collecting and analyzing beam images from laser cameras and profile monitors at LCLS. ImgMan consists of the Image Acquisition, Image Browser, and Image Analysis applications. For more details, check the User's Guide.

Image Acquisition

Image Acquisition is used to acquire background and beam images from cameras in the LCLS beam line. Some cameras are associated with a screen, which you must put IN first. After specifying the number of the desired background and beam images, press the "Measure" button to start the acquisition. You can monitor the status of your measurement on the main panel.

When the acquisition is complete, you can browse through your image datasets.

♥ Image Acquisition (DEVELO	PMENT) Version 1.0.2 (04-13-20	07)					- ×
YAG01 -	leasure Retrieve images	Cancel	Save	Load	Browse	Help	Exit Matlab
3 Background images	6 Beam images		Dataset	#1 (OTR1)			0 datasets
Screen is IN In	Out				Backg	pround image from	n 04-11-2007 20:02:30.816
							Show live image

You also have the ability to process live images and extract beam data from them.

Image Analysis

Image Analysis is used to analyze single images from your dataset. Images are processed, and beam data is extracted according to the parameters that you specify. You can trim the width or the height of the cropped area- depending on which slice plane you have selected.

Image Browser

Image Browser is used to browse through image datasets. Images are processed according to the parameters that you specify. You can set a master crop area for all images in a dataset.

