

LCLS2 Data Diagnostic Tools

- Set LCLS2 environment
- Data diagnostic commands
 - [datinfo](#)
 - [epix10ka_raw_calib_image](#)
 - [epix10ka_calib_components](#)
- [References](#)

Set LCLS2 environment

- [Set production and development environment](#)

Data diagnostic commands

Commands below can be run in lcls2 environment on psana nodes having access to experimental data under `/cds/data/psdm/<INSTRUMENT>/<EXPERIMENT>/xtc/*`.



Command examples are valid for code release ps-4.5.5, default version on 2021-10-25. They may evolve in other releases. Use "-h" option to find latest option of each command.

datinfo

```
USAGE:  
  datinfo -d <detector> -e <experiment> -r <run-number(s)> [kwargs]  
COMMAND EXAMPLES:  
  datinfo -d epixquad -e ueddaq02 -r 27 -td -L DEBUG  
  datinfo -d epixquad -e ueddaq02 -r 30-82 <--- DOES NOT WORK - missconfigured  
  datinfo -d epixquad -e ueddaq02 -r 83 <--- dark  
  datinfo -d epixquad -e ueddaq02 -r 84 <--- PARTLY WORKS charge injection  
  datinfo -d epixquad -f /cds/data/psdm/ued/ueddaq02/xtc/ueddaq02-r0065-s001-c000.xtc2  
  datinfo -d epixquad -f /cds/data/psdm/ued/ueddaq02/xtc/ueddaq02-r0086-s001-c000.xtc2  
  datinfo -d tmoopal -e tmoc00118 -r 123 -td  
HELP: datinfo -h
```

epix10ka_raw_calib_image

```

USAGE:
epix10ka_raw_calib_image [test-name] [optional-arguments]
where test-name:
  raw - test_raw W/O GRAPHICS
  calib - test_calib W/O GRAPHICS
  [image] - test_image WITH GRAPHICS - optional default
  mask, peds, rms, status, gain, z/y/xcoords - test_single_image WITH GRAPHICS
COMMAND EXAMPLES:
  epix10ka_raw_calib_image raw -e ueddaq02 -d epixquad -r66 # raw
  epix10ka_raw_calib_image calib -e ueddaq02 -d epixquad -r66 # calib
  epix10ka_raw_calib_image mask -e ueddaq02 -d epixquad -r66 # mask
  epix10ka_raw_calib_image -e ueddaq02 -d epixquad -r66 -N1000 # image
  epix10ka_raw_calib_image -e ueddaq02 -d epixquad -r108 -N1 -S grind
  epix10ka_raw_calib_image -e ueddaq02 -d epixquad -r140 -N100 -M2 -S calibcm8
  epix10ka_raw_calib_image -e ueddaq02 -d epixquad -r140 -N100 -M2 -S calibcm8 -o img-ueddaq02-epixquad-r140-
ev0002-cm8-7-100-10.png -N3
  epix10ka_raw_calib_image -e ueddaq02 -d epixquad -r211 -N1 -M0 -Speds -g0 # - plot pedestals for gain group
0/FH
  epix10ka_raw_calib_image -e ueddaq02 -d epixquad -r211 -N100 -Sraw-peds -M2 -g2 # - plot calib[step=2] -
pedestals[gain group 2]
  epix10ka_raw_calib_image -e rixx45619 -d epixhr -r118 --gramin 1 --gramax 32000 -Sraw
  epix10ka_raw_calib_image -e rixx45619 -d epixhr -r118 --gramin 1 --gramax 32000 -Speds -g0
  epix10ka_raw_calib_image -e rixx45619 -d epixhr -r118 --gramin -100 --gramax 100 -Sraw-peds -g0
  epix10ka_raw_calib_image -e rixx45619 -d epixhr -r119 -Scalib
  epix10ka_raw_calib_image -e rixx45619 -d epixhr -r119 -Sones
  epix10ka_raw_calib_image -e rixx45619 -d epixhr -N10000 -J200 --gramin 0 --gramax 10 -Sgrind
  epix10ka_raw_calib_image mask -e rixx45619 -d epixhr -r119
  epix10ka_raw_calib_image peds -e rixx45619 -d epixhr -r119 -gl
  epix10ka_raw_calib_image gains -e rixx45619 -d epixhr -r119 -gl
  epix10ka_raw_calib_image xcoords -e rixx45619 -d epixhr -r119 --cframe=1
HELP: epix10ka_raw_calib_image -h

```

epix10ka_calib_components

```

USAGE:
epix10ka_calib_components -r554 -t1
epix10ka_calib_components -e ueddaq02 -d epixquad -r554 -t1
-t, --tname - test number:
  1 - segment numeration
  2 - gain range index
  3 - gain, ADU/keV
  4 - pedestals
  5 - rms
  6 - raw
  7 - raw-peds
  8 - (raw-peds)/gain, keV
  9 - calib, keV
  10 - status
  11 - gain factor = 1/gain, keV/ADU
  ---
  21 - run 401 two-threshold selection issue
  22 - (raw-peds)/gain, keV hot - specific isuue test
  23 - (raw-peds)/gain, keV cold - specific isuue test

epix10ka_calib_components - test per-event components of the det.raw.calib method

HELP: epix10ka_calib_components -h

```

References

- [EPIX10KA2M References](#)
- [Set production and development environment](#)
- [github: psana/app](#)