

# 3.1 Post-processing functions

- Common
- Azimuthal integration
  - pyFAI
- Fourier transform
- Autocorrelation

## Common

The additional analysis function are passed to the cube workflow by adding a `det_proc` keys to the detector definition dictionary. The content of this key must be another dictionary with the following structure:

```
{  
    name: '<function_name>' # must be the exact DetObjectFunc name,  
    'func_kwarg1': <...>,  
    'func_kwarg2': <...>,  
    ...  
}
```

The function `<function_name>` will then be applied to the summed binned detector data and its output will be saved under `<detector_name>_variable` in the binned hdf5 file, in addition to the full detector data.

## Azimuthal integration

### pyFAI

Detector definition section with azimuthal integration of the Rayonix:

```
# ##### DETECTORS #####  
# List detectors to be cubed. Area detector have additional options such as threshold  
# Full images will always be saved.  
# More area detector can be defined following the same syntax, and adding them to varList  
detDict = {'source':'Rayonix',  
           'full':1,  
           'image':1,  
           'thresADU':-1e5,  
           'common_mode':0}  
  
pix_size = 176e-6  
func_kwargs = {  
    'name': 'azav_pyfai', # must be the name of the smalldata_tools DetObjectFunc  
    'ai_kwargs': {'dist':1, 'poni1':960*pix_size, 'poni2':960*pix_size},  
    'npts': 512,  
    'int_units' : '2th_deg',  
    'return2d' : False  
}  
det_proc = [func_kwargs]  
detDict['det_proc'] = det_proc  
  
# make list of all variables to be added to the cube  
varList = ['ipm_dg1/sum','ipm_dg2/sum', detDict]
```

The h5 resulting from this detector definition then contains:

- Rayonix\_azav (Array)
- Rayonix\_data (Array)
- Rayonix\_nEntries (Array)
- Rayonix\_q (Array)

The available keyword arguments to pass to the `azav_pyfai` function are

## azav\_pyfai documentation

```
"""
Parameters
-----
name: str
    Function name

mask: array, optional
    User defined mask. 1 for valid pixels.

return2d: bool, optional
    Return a 2d (q,phi). Default: False

poni_file: str, Path object, optional
    Path to a pyFAI calibration file

ai_kwargs: dict, optional
    Arguments to pyFAI.AzimuthalIntegrator. Either this parameter or a calib file is necessary
    For arguments see: https://pyfai.readthedocs.io/en/master/api/pyFAI.html#module-pyFAI.
azimuthalIntegrator

pol_factor: float, optional
    Polarization factor. Default 1. Passed to integrate1d or integrate2d.

npts_radial: int, optional
    Number of points for the radial binning. Default 256.

npts_az: int, optional
    Number of points for the azimuthal binning. Default 360. Only used for the 2d integration.

azav_kwargs: dict, optional
    Additional arguments to pass to integrate1d or integrate2d.
    See https://pyfai.readthedocs.io/en/master/api/pyFAI.html#pyFAI.azimuthalIntegrator.AzimuthalIntegrator.
integrate1d
"""


```

## Fourier transform

TBD

## Autocorrelation

TBD