New big compressor alignment procedure




Alignment procedure originated following 2022 reconfiguration project

This procedure is for full alignment of the SPL big compressor. A camera for viewing the IR beam reflection on the iris is recommend.


Step-by-step guide

- 1. send the beam to grating 1 (G1): align incident beam through iris 1 &2
- 2. normalize G1: turn G1 to retro beam back to iris1. Record the reading on G1-R (5.5°)
- 3. send the beam to grating 2 (G2): turn G1 by 17° so that the beam is reflected onto the center of G2. Record the reading on G1-R (-11.5° = 5.5° -
- 4. normalize grating 2 (G2): turn G2 to retro back to the iris1. Record the reading on G2-R(-33.3º)
- 5. zero G1 & G2: turn both G1 and G2 to zero angles (5.5° & -33.3°+34°)
- 6. working angles: turn both G1 and G2 to 1st order angles $(5.5^{\circ} 56.1^{\circ} = -50.6^{\circ} \& 0.7^{\circ} 56.1^{\circ} = 55.4^{\circ})$.
- 7. rotate the roof mirror so that the "plumb bob" shadows overlap to each other
- 8. adjust the linear stage (G2 T)to get the shortest pulse on SSA

	G1 rotation (deg)	G2 rotation(deg)
Step 2: Retro G1	5.5	
Step 3: center @ G2	5.5-17	
Step 4: Retro G2	5.5-17	-33.3
Step 5: Zero G1 & G2	5.5	-33.3+34=0.7
Step 6: 1st order	5.5-56.1=-50.6	0.7-56.1=-55.4



clockwise (+); counterclockwise (-)



Related articles

- Fully power down the SPLCold start the SPL
- New big compressor alignment procedure