Running CCD in the ASC (in the test chamber)

RIXS-CCD Controls and Interlock Scheme (from 11/30 meeting with controls/DAQ)

Meeting with Patrick/Kaz/Bhavna (2/22):

- Cryotel IOC (Patrick's update)
 - · IOC screen has PID controls, but since we only use PID in archon, that section will be wiped before we give this to users • For cryotel screen:
 - Hit only POWER mode: constant cooling, 60W
 - ON mode aka temperature mode: uses PID to get there
 - $^{\circ}$ 3/8" tubing for chiller, raise flow rate
 - Software wise, everything is in place (trigger from the archon, timing, etc)
 - RIX:STA:CCD:01 as the base name for the PVs
- · Waiting to hear from Dan about archon IOC and testing sensor in RIXS daq
- Next steps: meet with Dan, test basic functions, do a full cool-down/data collection/warmup at RIXS

This is a how-to for a single sensor. Ultimately, there will be two of these sensors at qRIXS.

Devices involved:

- Archon controller
- CryoTel controller
- Water chiller (small)
- Sensor with cryotel coldhead
- Laptop (PCxxxxx)

8 cable connections:

- Water lines to/from chiller x2
- Glenair signal cable from flange to archon
- 1-to-3 breakout controls cable between CryoTel coldhead/flange and cryotel controller
 - On the controller, one connection.
 - 0 On the coldhead/flange side, three connections: Two black connectors which go to coldhead, and one db9 connector that goes to the flange.
- · Cable that goes from CryoTel controller to coldhead
- Cat5 from Archon controller to laptop
- RS232 to USB cable from cryotel controller (RS232) to laptop (USB)
- Trigger cable to bottom BNC port on Archon

Navigating ArchonGUI:

- 1. Run "./archongui -small" (the -small option allows for a re-sizeable window)
- 2. File Open *whatever_name_of_config_file*.acf
- Hit "Connect", "Auto Fetch", "Apply all"
 Turn "Power On" on lower right, this is bias voltage
- 5. In "Slot 3: HeaterX"
 - a. Sensor C should be RTD1000
 - b. Heater A parameters:
 - i. Target: -100C
 - ii. Sensor: Sensor C
 - iii. Force: unchecked
 - iv. Limit(V): 25.0
 - v. P: 50
 - vi. I: 1
 - vii. D: 10
 - viii. IL: 20
 - ix. Ramp enable: checked
 - x. Ramp rate: 20mK/tick
 - c. Hit "Apply" (*note: hitting "Apply" will reset any values the PID had accumulated up to that point, like a refresh)

Cooldown takes ~2hrs.

CryoTel Software:

*note: The cryotel has a PID, but we don't use it. We use the other PID (via the Archon) that regulates the heater near/on the sensor. Keeping pwout to a minimum (i.e 60) keeps vibrations to a minimum.

- 1. Open new terminal: sudo putty /dev/ttyUSB0 -serial -sercfg 9600,8,n,1,N
- 2. Commands
 - a. View status: "status'

- b. Check cooler power: "pwout" (this should be 60. If not, set with "pwout=60")
- c. Turn cryotel on: "cooler=power"
 d. Turn cryotel off: "cooler=OFF" (might need to enter this command more than once to actually turn power off)

-----Only applicable for bench testing------

Parameters on waveform generator:

- Waveform: Pulsed mode
- Frequency: 200Hz (or anything under 250Hz)
- Pulse Height: 3.3V
- Pulse Offset: 1.65V
- Pulse Width: 1us
- Burst Mode: ON (lately it seems you must have it on to see anything on detector)

Turning things off (warming up CCD):

- 1. Turn cryotel off (cooler=off)
- 2. Turn bias voltage off (power off on gui)
- 3. Set target temp to +10C
- 4. That's it

Troubleshooting issues and lessons learned:

- If frames aren't updating, it's a trigger issue.
 How to save data: NO TILDES, e.g. "/home/bnayak/test" and data gets saved as "test_4800x300_xx.raw" where xx is the frame count
- Detector constantly running, "Save Sequence" just determines how many get data files get written
 Make sure temperature is actually at the target temp (DO NOT CLICK ramp enable)