

# DAQ User Guide (&troubleshooting)

- [Common features](#)
  - [DAQ epicsArch - Archived Data](#)
    - [Adding a single variable to the file](#)
    - [Add a whole file with a bunch of variables](#)
- [LCLS2 DAQ](#)
- [LCLS1 DAQ](#)
  - [DAQ Troubleshooting](#)
  - [Operating the DAQ](#)
  - [AMI Online Monitoring](#)
  - [Python Scripting](#)

## Common features

### DAQ epicsArch - Archived Data

You can add/remove any data you want in the epicsArchive data group readout.

The path is /cds/group/pcds/dist/pds/tmo/misc/

The tmo file we use is epicsArch\_tmo.txt. Let's look at two features of this

#### Adding a single variable to the file

Add a comment with a #

# Power meter at IM2K4

Add a pointer name of your choosing with \*

\*IM2K4\_XrayPower

For the **LCLS2 DAQ**, you should also specify which protocol is used to access the PV:

Add the pv name specifying whether it is read using channel access (ca), or pv access (pva) at the end

IM2K4:PPM:SPM:VOLT\_RBV ca

For the **LCLS1 DAQ**, you only list the PV

IM2K4:PPM:SPM:VOLT\_RBV

#### Add a whole file with a bunch of variables

You can make file with a subset of PVs, give it a meaningful name like epicsArch\_vmi.txt for all the vmi variables to record. You can then reference that file with '<' in the main file as below.

# MPOD/VMI information

< epicsArch\_vmi.txt

## LCLS2 DAQ

[LCLS-II Data Acquisition and Analysis](#)

## LCLS1 DAQ

[Online Analysis Tutorial \(AMI\)](#) - presented at 2014 SSR/LCLS Users Meeting: LCLS Data Analysis Workshop

[DAQ Overview and Online Analysis](#) - presented at 2010 SSRL/LCLS Users Meeting : LCLS Data Analysis Workshop

[Online Analysis Design](#) - Matt's presentation for Jan 2012 DAQ Meeting

## DAQ Troubleshooting

Trouble Shooting guide for scientists or first level responders.

## Operating the DAQ

- Editing a configuration
  - Example: Opal camera configuration
  - Changing a configuration while running
- Selecting detectors for readout
- Starting/stopping a run
- Watching progress of a run
- Running scans
  - Launching a Scan from DAQ Control GUI (deprecated)
  - Launching Scans Remotely (scripted)
- Running the sequencer
- Configuring the EVR
- Detector timing settings
- Recovering from errors: restarting the DAQ

## AMI Online Monitoring

- Online Monitoring GUI
  - Using the Online Monitoring GUI
  - Writing a plug-in to the Online Monitoring GUI
  - Writing a user application (reads from shared memory) (deprecated old C++ psana)
  - Online Monitoring and Simulation Using Files (deprecated old C++ psana)
  - Writing a user application, offline analysis style (reads from a file) (deprecated old C++ psana)
- XTC playback (a.k.a Offline AMI)

## Python Scripting

- Editing configurations - the 'pycdb' module
- Controlling the DAQ - the 'pydaq' module
- Monitoring the data - the 'pyami' module