


# Differential spectrum from Opal camera image

- Motivation
- Psana configuration file to get image:
- Project of analysis-specific modules
  - Module producing spectra
    - Functionality
    - Module parameters
  - Module example for analysis of spectra
    - Functionality
    - Module parameters

## Motivation

- 2012-09-13 14:00 - Mtg with Ryan Coffee, Josef Frisch, Nick Hartmann
- Would like to have psana application to process xppi0412 data



Psana configuration file to get image:

```

[psana]
files      = \
/reg/d/psdm/xpp/<experiment>/xtc/<file-1>.xtc \
/reg/d/psdm/xpp/<experiment>/xtc/<file-2>.xtc \
...
/reg/d/psdm/xpp/<experiment>/xtc/<file-N>.xtc

modules    = \
            ImgAlgos.CameraImageProducer \
            ImgAlgos.SaveImageInFile
#          psana_examples.DumpOpallk \
events     = 3

[ImgAlgos.CameraImageProducer]
source     = DetInfo(:Opall000)
key_in     =
key_out    = img
subtract_offset = true
print_bits = 15

[ImgAlgos.SaveImageInFile]
source     = DetInfo(:Opall000)
key        = img
fname      = img-spectra
saveAll    = true
#eventSave = 1

```

In order to run this script, the xtc file path and names should be substituted, then use commands:

```

sit_setup
psana -c <the-name-of-the-configuration-file>

```

## Project of analysis-specific modules

### Module producing spectra

#### Functionality

- Define and cut the spectral regions for signal and reference band.
- Integrate in columns and evaluate two 1024-bin spectra, fs and fr, and their relative difference  $r = 2(fs - fr)/(fs + fr)$ .
- Put 3 spectra in the event as `ndarray<double, 2>` object.
- Save spectra in the text file, 1024 x 3 lines of float numbers, if necessary.

#### Module parameters

- source
- key\_in
- key\_out
- Signal band Y-center
- Signal band width
- Signal band tilt angle
- Reference band Y-center
- Reference band width
- Reference band tilt angle
- Output file name (text format, 1024 x 3 lines of float numbers). If specified, then the file is saved for all selected events.

### Module example for analysis of spectra

#### Functionality

- get the `ndarray<double, 2>` object from the event.
- do iteration over vectors and print elements for example.



Then users extend this module functionality depending on their goal...

## Module parameters

- source
- key\_in
- other parameters for analysis or accumulation of statistics etc.