

# Myana Reference Manual

- [main](#)
  - [Namespaces](#)
  - [Defines](#)
  - [Enumerations](#)
  - [Typedefs](#)
  - [Classes](#)
  - [Functions](#)
- [myana.cc and myana.hh](#)
- [myana\\_morefeatures.cc and myana.hh](#)

Here is a list of all files with brief descriptions:

main.hh - declares enums and functions

main.cc - implements getters and other utility functions and classes.

myana.hh - declaration of user functions (beginrun, event etc)

myana.cc - implementation of user functions (beginrun, event etc)

myana\_morefeatures.cc - same as myana.cc, but has more examples.

## main

### Namespaces

These are the namespaces used.

```
namespace      Pds
namespace      Pds : CsPad
```

### Defines

Define some aliases.

```
#define      getOpallkConfig(x)      getFrameConfig(x)
#define      getTm6740Config(x)      getFrameConfig(x)
#define      getOpallkValue(w, x, y, z)      getFrameValue(w, x, y, z)
#define      getTm6740Value(w, x, y, z)      getFrameValue(w, x, y, z)

#define      FNAME_LEN      256
#define      OpallkDetectorIndex(x, y)      FrameDetectorIndex(x, DetInfo::Opall000, y)
#define      Tm6740DetectorIndex(x, y)      FrameDetectorIndex(x, DetInfo::TM6740 , y)
#define      FccdDetectorIndex(x, y)      FrameDetectorIndex(x, DetInfo::Fccd, y)
```

### Enumerations

These are "named integers" that are used to control e.g. which detector to fetch information from.

Most of these are copied from pdsdata/xtc/DetInfo.hh.

```

enum          AcqDetector {
    AmoIms = 0, AmoGasdet = 1, AmoETof = 2, AmoITof = 3,
    AmoMbes = 4, AmoVmiAcq = 5, AmoBpsAcq = 6, Camp = 7,
    SxrBeamlineAcq1 = 0, SxrBeamlineAcq2 = 1, SxrEndstationAcq1 = 2, SxrEndstationAcq2 = 3,
    NumAcqDetector = 8
}
enum          FrameDetector {
    AmoVmi = 0, AmoBps1 = 1, AmoBps2 = 2, SxrBeamlineOpall = 0,
    SxrBeamlineOpal2 = 1, SxrEndstationOpall = 2, SxrEndstationOpal2 = 3, SxrFccd = 4,
    XppSb1PimCvd = 0, XppMonPimCvd = 1, XppSb3PimCvd = 2, XppSb4PimCvd = 3,
    NumFrameDetector = 5
}
enum          PrincetonDetector {
    SxrBeamlinePrinceton1 = 0, SxrBeamlinePrinceton2 = 1,
    SxrEndstationPrinceton1 = 2, SxrEndstationPrinceton2 = 3,
    AmoPrinceton1 = 0, AmoPrinceton2 = 1,
    AmoPrinceton3 = 2, AmoPrinceton4 = 3,
    NumPrincetonDetector = 4
}
enum          IpimbDetector {
    SxrBeamlineIpimb1 = 0, SxrBeamlineIpimb2 = 1, SxrBeamlineIpimb3 = 2,
    SxrBeamlineIpimb4 = 3, SxrBeamlineIpimb5 = 4, SxrBeamlineIpimb6 = 5,
    SxrBeamlineIpimb7 = 6, SxrBeamlineIpimb8 = 7, SxrBeamlineIpimb9 = 8,
    SxrBeamlineIpimb10 = 9, SxrBeamlineIpimb11 = 10, SxrBeamlineIpimb12 = 11,
    SxrBeamlineIpimb13 = 12, SxrBeamlineIpimb14 = 13, SxrBeamlineIpimb15 = 14,
    SxrBeamlineIpimb16 = 15,
    SxrEndstationIpimb1 = 16, SxrEndstationIpimb2 = 17, SxrEndstationIpimb3 = 18,
    SxrEndstationIpimb4 = 19, SxrEndstationIpimb5 = 20, SxrEndstationIpimb6 = 21,
    SxrEndstationIpimb7 = 22, SxrEndstationIpimb8 = 23, SxrEndstationIpimb9 = 24,
    SxrEndstationIpimb10 = 25, SxrEndstationIpimb11 = 26, SxrEndstationIpimb12 = 27,
    SxrEndstationIpimb13 = 28, SxrEndstationIpimb14 = 29, SxrEndstationIpimb15 = 30,
    SxrEndstationIpimb16 = 31,
    XppSb1Ipm = 0, XppSb1Pim = 1, XppMonPim = 2, XppSb2Ipm = 3,
    XppSb3Ipm = 4, XppSb3Pim = 5, XppSb4Pim = 6, NumIpimbDetector = 32
}
enum          PnCcdDetector { PnCcd0 = 0, PnCcd1 = 1, NumPnCcdDetector = 2 }
enum          EncoderDetector { SxrBeamlineEncoder1 = 0, NumEncoderDetector = 1 }
enum          Result { OK, End, Error }

```

## Typedefs

These are defined in main.cc.

```

typedef vector< AcqWaveForm >  AcqChannelList
typedef int(*)  TEpicValueFuncPtr (const EpicsPvHeader &epicsPv1, const void **ppValue,
                                   int *piDbrype, struct tm *pTmTimeStamp, int *piNanoSec)

```

## Classes

These are all defined in main.cc.

- **struct AcqWaveForm**  
Public Attributes

```

double          vfTrig
vector< double >  vfTime
vector< double >  vfVoltage

```

- **class myLevelIter**  
Public Member Functions

```

// Constructor
myLevelIter (Xtc *xtc, unsigned depth, int iDebugLevel)

void          process( const DetInfo &di, const ... &cfg) // lots of these
// one process function for each device configuration.

int           process (Xtc *xtc)

```

- **class XtcSlice**  
Public Member Functions

```

XtcSlice (std::string fname)
~XtcSlice ()
bool      add_file (std::string fname)
void      init ()
Result    next (Pds::Dgram *dg)
Result    skip ()
const Pds::Dgram &      hdr () const

```

- **class XtcRun**  
Public Member Functions

```

XtcRun ()
~XtcRun ()
void      reset (std::string fname)
bool      add_file (std::string fname)
const char *      base () const
void      init ()
Result      next (Pds::Dgram *dg)

```

## Functions

These are all declared in `main.hh` and implemented in `main.cc`.

Time functions and file name

```

int      getTime (int &seconds, int &nanoSeconds)
int      getLocalTime (const char *&time)
int      getFiducials (unsigned &fiducials)
const char *      getInputName ()
const char *      getFileName ()
void      saveName (const char *file, char *dest)

```

Configuration retrieval functions:

```

int      getAcqConfig (AcqDetector det, int &numChannels, int &numSamples, double &sampleInterval)
int      getFrameConfig (FrameDetector det)
int      getPrincetonConfig (Pds::DetInfo::Detector det, int iDevId,
                             int &width, int &height,
                             int &orgX, int &orgY,
                             int &binX, int &binY)
int      getIpimbConfig (Pds::DetInfo::Detector det, int iDevId, uint64_t &serialID,
                             int &chargeAmpRange0, int &chargeAmpRange1,
                             int &chargeAmpRange2, int &chargeAmpRange3)
int      getEncoderConfig (Pds::DetInfo::Detector det, int iDevId)
int      getFccdConfig (FrameDetector det, uint16_t &outputMode,
                             bool &ccdEnable, bool &focusMode, uint32_t &exposureTime,
                             float &dacVoltage1, float &dacVoltage2, float &dacVoltage3,
                             float &dacVoltage4, float &dacVoltage5, float &dacVoltage6,
                             float &dacVoltage7, float &dacVoltage8, float &dacVoltage9,
                             float &dacVoltage10, float &dacVoltage11, float &dacVoltage12,
                             float &dacVoltage13, float &dacVoltage14, float &dacVoltage15,
                             float &dacVoltage16, float &dacVoltage17,
                             uint16_t &waveform0, uint16_t &waveform1, uint16_t &waveform2,
                             uint16_t &waveform3, uint16_t &waveform4, uint16_t &waveform5,
                             uint16_t &waveform6, uint16_t &waveform7, uint16_t &waveform8,
                             uint16_t &waveform9, uint16_t &waveform10, uint16_t &waveform11,
                             uint16_t &waveform12, uint16_t &waveform13, uint16_t &waveform14)
int      getDiodeFexConfig (Pds::DetInfo::Detector det, int iDevId, float *base, float *scale)
int      getIpmFexConfig (Pds::DetInfo::Detector det, int iDevId,
                             float *base0, float *scale0, float *base1, float *scale1,
                             float *base2, float *scale2, float *base3, float *scale3,
                             float &xscale, float &yscale)
int      getCspadConfig (Pds::DetInfo::Detector det, unsigned &quadMask, unsigned &asicMask)
int      getCspadConfig (Pds::DetInfo::Detector det, Pds::CsPad::ConfigV1 &cfg)

```

#### L1Accept (event/shot) Data retrieval functions

```

int      getAcqValue (AcqDetector det, int channel, double *time, double *voltage)
int      getAcqValue (AcqDetector det, int channel, double *time, double *voltage, double &trigtime)
int      getFrameValue (FrameDetector det, int &frameWidth, int &frameHeight, unsigned short *image)
int      getPrincetonValue (Pds::DetInfo::Detector det, int iDevId, unsigned short *image)
int      getPrincetonTemperature (Pds::DetInfo::Detector det, int iDevId, float &temperature)
int      getIpimbVolts (Pds::DetInfo::Detector det, int iDevId,
                             float &channel0, float &channel1, float &channel2, float &channel3)
int      getFeeGasDet (double *shotEnergy)
int      getEBeam (double &charge, double &energy, // charge and energy
                             double &posx, double &posy, // position (x,y)
                             double &angx, double &angy) // angle (x,y)
int      getEBeam (double &charge, double &energy, // charge and energy
                             double &posx, double &posy, // position (x,y)
                             double &angx, double &angy, // angle (x,y)
                             double &pkcurr) // current
int      getPhaseCavity (double &fitTime1, double &fitTime2, double &charge1, double &charge2)
int      getPvInt (const char *pvName, int &value)
int      getPvFloat (const char *pvName, float &value)
int      getPvString (const char *pvName, char *value)
int      getPnCcdValue (int deviceId, unsigned char *image, int &width, int &height)
int      getEvrDataNumber ()
int      getEvrData (int id, unsigned int &eventCode, unsigned int &fiducial, unsigned int &timeStamp)
int      getEncoderCount (Pds::DetInfo::Detector det, int iDevId, int &encoderCount)
int      getDiodeFexValue (Pds::DetInfo::Detector det, int iDevId, float &value)
int      getIpmFexValue (Pds::DetInfo::Detector det, int iDevId, float *channels,
                             float &sum, float &xpos, float &ypos)
int      getCspadQuad (Pds::DetInfo::Detector det, unsigned quad, const uint16_t *pixels)
int      getCspadQuad (Pds::DetInfo::Detector det, unsigned quad, const Pds::CsPad::ElementV1 *data)

```

#### Control data retrieval functions

```

int      getControlPvNumber ()
int      getControlPvName (int pvId, const char *&pvName, int &arrayIndex)
int      getControlValue (const char *pvName, int arrayIndex, double &value)
int      getMonitorPvNumber ()
int      getMonitorPvName (int pvId, const char *&pvName, int &arrayIndex)
int      getMonitorValue (const char *pvName, int arrayIndex, double &hilimit, double &lolimit)

```

#### Advanced interface

```

int      getEpicsPvNumber ()
int      getEpicsPvConfig (int iPvId, const char *&sPvName, int &iType, int &iNumElements)
int      getEpicsPvValue (int pvId, const void *&value,
                          int &dbtype, struct tm &tmTimeStamp, int &nanoSec)
void      fillConstFrac (double *t, double *v, unsigned numSamples,
                        float baseline, float thresh, TH1 *hist)
void      fillConstFrac (double *t, double *v, unsigned numSamples,
                        float baseline, float thresh, double *edge, int &n, int maxhits)

```

#### Program control functions

```

void      usage (char *programe)
void      makeoutfilename (char *filename, char *outfilename)
void      anarun (XtcRun &run, unsigned &maxevt, unsigned &skip, int iDebugLevel)
int       main (int argc, char *argv[])

```

## myana.cc and myana.hh

#### Functions

```

void      beginjob ()
void      beginrun ()
void      begincalib ()
void      event ()
void      endcalib ()
void      endrun ()
void      endjob ()

```

## myana\_morefeatures.cc and myana.hh

#### Functions

```

void      beginjob ()
void      beginrun ()
void      begincalib ()
void      event ()
void      endcalib ()
void      endrun ()
void      endjob ()

```