

# Optical Metrology Quality Check

This note is about application (script) which allows to check quality of the optical metrology input during measurement.

## Content

- [Content](#)
- [Introduction](#)
- [Get application up and running](#)
  - [Direct usage](#)
  - [Use from psana release](#)
- [Example](#)
  - [Output](#)
  - [Parameters to check](#)
  - [Plots](#)
- [References](#)

## Introduction

Optical metrology measurements contains enough information for self check. Historically this quality check is done offline when entire detector or its quad is measured and entire file is available. Sometimes this file contains typo, which can be fixed using over-defined rectangular geometry of the detector segments. However, it could be much better to fix this problem during optical measurements with microscope. Simple application which can check partial optical metrology input is now available and described in this note.

## Get application up and running

### Direct usage

Optical metrology quality check (QC) application is implemented as a self consistent script which uses standard python features with numpy and matplotlib libraries. It can be downloaded from

[https://github.com/lcls-psana/CalibManager/blob/master/app/optical\\_metrology\\_check](https://github.com/lcls-psana/CalibManager/blob/master/app/optical_metrology_check)

and used as is by the command

```
python optical_metrology_check <opt-met-file.txt>
```

### Use from psana release

For development purpose it could be better to use it through psana-conda release system. Code is available since ana-1.2.8. Example:

#### Login to psana and set psana environment

```
# login to one of psana nodes
ssh -Y pslogin.slac.stanford.edu
ssh -Y psana.slac.stanford.edu

# set psana environment if it is missing in your .bash_profile or .bash_login
source /reg/g/psdm/etc/ana_env.sh
```

#### Setup conda environment, directory etc.

```
# setup conda environment
source conda_setup

# go to any directory
mkdir my-test-dir
cd my-test-dir
```

### Run script on optical metrology file

```
optical_metrology_check -h  
optical_metrology_check <opt-met-file.txt>
```

## Example

A few complete optical metrology files can be found under directory /reg/g/psdm/detector/alignment/cspad/, e.g. [calib-xpp-2016-08-22-camera3/2016-08-22-CSPAD-XPP-Camera3-Metrology.txt](#)

However, script can run on partial metrology file, e.d.: [metrology\\_test.txt](#), which can be downloaded and used for test purpose:

```
optical_metrology_check metrology_test.txt
```

## Output

Running script prints a few section of results, e.g.:

### Command line and application input/default parameters

```
(ana-1.2.8) [dubrovin@psanagpu104:~/my-test-dir]$ optical_metrology_check metrology_test.txt  
Command:  
  /reg/g/psdm/sw/conda/inst/miniconda2-prod-rhel7/envs/ana-1.2.8/bin/optical_metrology_check metrology_test.txt  
Argument list: ['metrology_test.txt']  
Optional parameters:  
  <key>      <value>      <default>  
  psz       109.92      109.92  
  vrb       15          15  
  ifn       ./optical_metrology.txt  ./optical_metrology.txt  
  plt       True        True  
  ofp       ./om        ./om  
-----
```

### Input file parsing line-by-line

```
In read_optical_metrology_file: metrology_test.txt
EMPTY LINE IS IGNORED
IS RECOGNIZED QUAD: 1
EMPTY LINE IS IGNORED
COMMENT IS IGNORED : "# this is comment"
EMPTY LINE IS IGNORED
TITLE LINE:      Point      X      Y      Z
EMPTY LINE IS IGNORED
ACCEPT RECORD:   1          0          0      697
ACCEPT RECORD:   2         -15      20903      604
ACCEPT RECORD:   3      43529      20930      172
ACCEPT RECORD:   4      43542          23      297
EMPTY LINE IS IGNORED
ACCEPT RECORD:   9          159      46961      508
ACCEPT RECORD:  10           78      90506      321
ACCEPT RECORD:  11      20986      90547      156
ACCEPT RECORD:  12      21070      47003      343
EMPTY LINE IS IGNORED
ACCEPT RECORD:  13      23423      46973      341
ACCEPT RECORD:  14      23540      90519      116
ACCEPT RECORD:  15      44447      90457      -67
ACCEPT RECORD:  16      44342      46917       93
EMPTY LINE IS IGNORED
ACCEPT RECORD:  21      47105      69413        27
ACCEPT RECORD:  22      47103      90323       -40
ACCEPT RECORD:  23      90653      90322     -435
ACCEPT RECORD:  24      90651      69415     -365
EMPTY LINE IS IGNORED
ACCEPT RECORD:  25      44480          259      284
ACCEPT RECORD:  26      44381      43805      129
ACCEPT RECORD:  27      65305      43846      -50
ACCEPT RECORD:  28      65389          303      109
EMPTY LINE IS IGNORED
ACCEPT RECORD:  29      67741          149       98
ACCEPT RECORD:  30      67705      43731     -48
EMPTY LINE IS IGNORED
-----
```

### Point numeration check (consecutive 4 per segment)

```
check_points_numeration
WARNING - numeration of points in table is NOT sequential or started from non-x4 number
Measured segment 0 point numbers: ( 1  2  3  4) OK - points in segment are sequential
Measured segment 1 point numbers: ( 9 10 11 12) OK - points in segment are sequential
Measured segment 2 point numbers: (13 14 15 16) OK - points in segment are sequential
Measured segment 3 point numbers: (21 22 23 24) OK - points in segment are sequential
Measured segment 4 point numbers: (25 26 27 28) OK - points in segment are sequential
```

### Per segment in-plane (X-Y) quality check

X-Y quality check for optical metrology measurements

```
-----
segm:      S1      S2      dS1      dS2      L1      L2      dL1      dL2      angle(deg)  D1      D2      dD      d
(dS)    d(dL)
-----
segm: 0    20903    20907      23      27      43542    43544      -15      -13      0.03290    48302    48304
-1      -4      -2
segm: 1    20911    20908     -81     -84      43545    43544       42       41      0.10855    48307    48302
4        3        1
segm: 2    20919    20907     117     105      43546    43540     -56     -62     -0.14606    48301    48310
-8      12        6
segm: 3    20910    20907       2      -1      43546    43550      -2       2      0.00066    48309    48308
1        3      -4
segm: 4    20909    20924     -99     -84      43546    43543      44       41      0.12040    48307    48309
-1     -15        3
-----
```

### Quality check in Z measurement

Z quality check for optical metrology measurements

```
-----
segm:      SA      LA  XSize  YSize  dZS1 dZS2 dZL1 dZL2  dZSA dZLA ddZS ddZL  dZX  dZY
angXZ(deg) angYZ(deg) dz3(um)
-----
segm: 0    20905    43543    43543    20905     -93  -125  -400  -432    -109  -416   32   32   -416  -109
-0.54737  -0.29874   -31.962
segm: 1    20909    43544    20909    43544    -165  -165  -187  -187    -165  -187    0    0   -165  -187
-0.45212  -0.24605    0.028
segm: 2    20913    43543    20913    43543    -248  -183  -225  -160    -215  -192   -65  -65   -215  -192
-0.59039  -0.25330   -64.821
segm: 3    20908    43548    43548    20908     -67   -70  -392  -395     -68  -393    3    3   -393   -68
-0.51771  -0.18771   -2.973
segm: 4    20916    43544    20916    43544    -175  -179  -155  -159    -177  -157    4    4   -177  -157
-0.48484  -0.20658    3.885
-----
```

## Parameters to check

See for detail [CSPAD Geometry and Alignment](#)

dD (um) - difference between diagonals

d(dS) (um) - difference between two short sides

d(dL) (um) - difference between two long sides

dz3 (um) - deviation one of the measured points from the plane of three other points




All these parameters should be about zero within fraction of the pixel size. If not... re-measurement is desirable.





## Plots




Application generates and saves in current directory a few useful plots;

map of measured sensors for [metrology\\_test.txt](#) and [2016-08-22-CSPAD-XPP-Camera3-Metrology.txt](#) metrology files:



frequency distributions of a few important parameters, as shown for [2016-08-22-CSPAD-XPP-Camera3-Metrology.txt](#)





## References

- [CSPAD Geometry and Alignment](#)
- [CSPAD Alignment](#)